# Immunglobulin G ELISA

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
# Immunoglobulin G ELISA

*For the in vitro determination of immunoglobulin G  
in serum, plasma and urine*

Valid from 2022-01-05

**REF** K 6510

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>17</b>
<b>2. INTRODUCTION</b>	<b>17</b>
<b>3. MATERIAL SUPPLIED</b>	<b>17</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>18</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>19</b>
<i>EDTA plasma and serum samples</i>	19
<i>Urine</i>	20
<b>7. ASSAY PROCEDURE</b>	<b>20</b>
<i>Principle of the test</i>	20
<i>Test procedure</i>	20
<b>8. RESULTS</b>	<b>22</b>
<b>9. LIMITATIONS</b>	<b>22</b>
<b>10. QUALITY CONTROL</b>	<b>23</b>
<i>Reference range</i>	23
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>23</b>
<i>Accuracy – Precision</i>	23
<i>Accuracy – Trueness</i>	24
<i>Linearity</i>	24
<i>Analytical sensitivity</i>	25
<i>Analytical specificity</i>	25
<b>12. PRECAUTIONS</b>	<b>25</b>
<b>13. TECHNICAL HINTS</b>	<b>26</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>26</b>

## 1. INTENDED USE

This Immundiagnostik AG assay is an enzyme immunoassay intended for the quantitative determination of Immunoglobulin G (IgG) in plasma, serum and urine. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Immunoglobulin G (IgG) is the most common type of antibody found in human blood (~75% of serum antibodies). IgG molecules protect the body from infection by immobilisation and agglutination of antigens, opsonisation, activation of the classical complement pathway and other means. IgG molecules are also involved in the regulation of allergic reactions, and its small size allows it to easily perfuse tissues.

IgG molecules are secreted by plasma B cells. IgG is a monomer consisting of four peptide chains with a total molecular weight of about 150 kDa. It contains two identical class  $\gamma$  heavy chains of about 50 kDa and two identical light chains of about 25 kDa. The two heavy chains are linked to each other and to a light chain each by disulfide bonds. The resulting tetramer has two identical halves, which together form the Y-like shape. Each end of the fork contains an identical antigen binding site.

### Indications

- autoimmune hepatitis
- monoclonal gammopathy type IgG kappa or IgG lambda
- amyloidosis
- chronic infections

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6510	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 6510	NACL	0.9 % NaCl solution	1 x 30 ml
K 6510	CONJBUF	Conjugate dilution buffer, ready-to-use	1 x 22ml
K 6510	CONJ	Conjugate concentrate, (rabbit-anti-IgG, peroxidase-labelled)	1 x 50 $\mu$ l
K 6510	STD	Standards, lyophilised (0; 0.32; 0.8; 2; 5 mg/l)	5 x 1 vial

Cat. No.	Label	Kit components	Quantity
K 6510	CTRL	Control, lyophilised (see specification for range)	1 x 1 vial
K 6510	SAMPLEBUF	Sample dilution buffer, ready-to-use	2x 100 ml
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	2x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is

stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

- The **lyophilised standards (STD)** and **control (CTRL)** are stable at **2–8 °C** until the expiry date stated on the label. Before use, the STD and CTRL have to be reconstituted with **250 µl of ultrapure water** and mixed by gentle inversion to ensure complete reconstitution. Allow the vial content to dissolve for 10 minutes and then mix thoroughly. **Standards and control** (reconstituted STD and CTRL) **can be stored at -20 °C until the expiry date given on the label. Avoid repeated thawing and freezing.**
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:1000** in **conjugate dilution buffer** (10 µl CONJ + 10 ml CONJBUF). The **CONJ** is stable at **2–8 °C** until the expiry date stated on the label. **Conjugate** (1:1000 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *EDTA plasma and serum samples*

#### **Sample storage**

Freshly collected EDTA plasma or serum can be stored for up to 2 weeks at 2–8 °C or for longer storage at -20 °C.

#### **Sample dilution**

EDTA plasma or serum samples must be diluted **1:10 000** in sample dilution buffer (SAMPLEBUF) before performing the assay, e.g.

- **10 µl** sample + **990 µl** SAMPLEBUF, mix well = 1:100 (dilution I)
- **10 µl** dilution I + **990 µl** SAMPLEBUF, mix well = 1:100 (dilution II). This results in a final dilution of 1:10 000.

**10 µl** of dilution II are used in the test per well.

## Urine

### Sample storage

Adjust the urine to a pH of 6 to 8 with 1 N NaOH. Samples can be stored at 2–8°C up to 2 weeks or at –20°C for longer storage.

### Sample dilution

Samples with IgG-concentration higher than 5 mg/l must be diluted **1:10** with sample dilution buffer (SAMPLEBUF), e. g.

**100 µl** sample + **900 µl** SAMPLEBUF, mix well.

**10 µl** of the dilution are used in the test per well.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of immunoglobulin G.

In a first incubation step, the Immunoglobulin G in the samples is bound to polyclonal rabbit antibodies (in excess) immobilised to the surface of the microtitre wells. After removal of all unbound substances, a peroxidase-labelled anti Immunoglobulin G antibody is added. The second washing step is followed by incubation with the substrate, tetramethyl–benzidine (TMB). The reaction is terminated by an acidic stop solution converting the colour from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of Immunoglobulin G in the sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Immunoglobulin G, present in the patient samples, is determined directly from this curve.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Mark the positions of standards/control/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8°C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	<b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Add <b>200 µl</b> 0.9 % NaCl (NACL) into each well.
3.	Add each <b>10 µl standards/control/prepared or diluted samples</b> into the respective wells.
4.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal shaker*</b> .
5.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
6.	Add <b>200 µl conjugate</b> (diluted CONJ) into each well.
7.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal shaker*</b> .
8.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
9.	Add <b>200 µl substrate</b> (SUB) into each well.
10.	Incubate for <b>10–20 min**</b> at room temperature (15–30 °C) in the <b>dark</b> .
11.	Add <b>50 µl stop solution</b> (STOP) into each well and mix well.
12.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### EDTA-plasma and serum

The obtained results have to be multiplied by the **dilution factor of 10 000** to get the actual concentrations.

Urine (1:10 dilution)

The obtained results have to be multiplied by the **dilution factor of 10** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## 9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

Samples with concentrations lower than the measurement range cannot be clearly quantified.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

Plasma or serum: 8–18g/l

Urine: 0.5–3.2mg/24 hours

We recommend each laboratory to establish its own reference range. Above mentioned values are only for orientation and may vary from other published data.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay)**

The repeatability was assessed with 3 serum samples under **constant** parameters (same operator, instrument, day and kit lot). The following values have been obtained without considering the sample dilution factors:

Sample	Mean value [mg/l]	CV [%]
1 (n = 26)	0,479	5,7
2 (n = 26)	0,396	6,3
3 (n = 28)	0,278	6,3

#### **Reproducibility (Inter-Assay); n = 33**

The reproducibility was assessed with 2 serum samples under **varying** parameters (different operators, instruments, days and kit lots).

Sample	Mean value [mg/l]	CV [%]
1	0,244	13,2
2	1,312	12,4

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, one sample with a known concentration was tested 4 times:

Sample [g/l]	Obtained [g/l]	Recovery [%]
9.17	9.35	101.92
	9.22	100.55
	9.22	100.56
	8.72	95.12

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A with a serial dilution of 5 different serum samples.

For IgG in serum, plasma and urine, the method has been demonstrated to be linear from 0.079 to 2.000 mg/l showing a non-linear behaviour of less than  $\pm 20\%$  in this interval.

Sample	Dilution	Expected [mg/l]	Obtained [mg/l]	Recovery [%]
A	1:5 000	1.329	1.329	100.00
	1:12 500	0.532	0.493	92.80
	1:31 250	0.213	0.195	91.93
	1:78 125	0.085	0.073	85.91
B	1:5 000	1.234	1.234	100.00
	1:12 500	0.493	0.464	93.99
	1:31 250	0.197	0.175	88.84
	1:78 125	0.079	0.065	81.71
C	1:10 000	0.837	0.837	100.00
	1:25 000	0.335	0.347	103.53
	1:62 500	0.134	0.125	93.52

Sample	Dilution	Expected [mg/l]	Obtained [mg/l]	Recovery [%]
D	1:1 600	2.000	2.000	100.00
	1:4 000	0.800	0.763	95.43
	1:10 000	0.320	0.313	97.84
	1:25 000	0.128	0.132	102.76
E	1:1 600	1.900	1.900	100.00
	1:4 000	0.760	0.728	95.84
	1:10 000	0.304	0.304	99.92
	1:25 000	0.122	0.128	105.02

### *Analytical sensitivity*

The following value has been estimated based on the concentrations of the standards without considering possibly used sample dilution factors.

Limit of blank, LoB

0.062 mg/l

### *Analytical specificity*

There was no cross-reactivity observed with mouse or rat IgG.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still should be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant