# L-Citrullin

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
4. Yu, Jin-Ju, and Suk-Heung Oh. 2010. "Isolation and Characterization of Lactic Acid Bacteria Strains with Ornithine Producing Capacity from Natural Sea Salt." *Journal of Microbiology (Seoul, Korea)* **48** (4) (August): 467–72. doi:10.1007/s12275-010-0204-9.

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
# L-Citrulline

*For the in vitro determination of L-citrulline in urine, serum  
and plasma*

Valid from 2021-05-11

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com) [www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>13</b>
<b>2. INTRODUCTION</b>	<b>13</b>
<b>3. MATERIAL SUPPLIED</b>	<b>13</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>14</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>14</b>
<b>6. SPECIMEN COLLECTION AND PREPARATION</b>	<b>15</b>
<b>7. ASSAY PROCEDURE</b>	<b>16</b>
<i>Principle of the test</i>	16
<i>Test procedure</i>	16
<b>8. RESULTS</b>	<b>17</b>
<b>9. LIMITATIONS</b>	<b>17</b>
<b>10. QUALITY CONTROL</b>	<b>17</b>
<i>Reference range</i>	18
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>18</b>
<i>Precision and reproducibility</i>	18
<i>Dilution recovery</i>	18
<i>Analytical sensitivity</i>	18
<b>12. PRECAUTIONS</b>	<b>19</b>
<b>13. TECHNICAL HINTS</b>	<b>19</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>19</b>
<b>15. REFERENCES</b>	<b>20</b>

## 1. INTENDED USE

This Immundiagnostik AG assay is intended for the quantitative determination of L-citrulline in urine, serum and plasma. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

**Nitric oxide (NO)** is an intra- and intercellular signaling molecule. It reacts with free radicals, metalloproteins and specific amino acid residues of proteins. NO plays an important role in the regulation of vascular tone. **Endothelial NO (eNO)** is produced by the vascular endothelium. It diffuses to neighbouring vascular smooth muscle cells (VSMC), where NO activates soluble guanylate cyclase (sGC), which subsequently increases the intracellular cGMP production from GTP, and which in turn causes relaxation of smooth muscle and vasodilatation.

Thus, functional changes of the endothelium in coronary artery disease may be an important factor in the development of vasospasm, ischaemia and thrombosis.

### L-citrulline as surrogate marker for NO

NO is synthesised in the **citrulline-NO-cycle** when L-arginine is oxidised to citrulline by NO synthase (NOS). In the second part of the urea cycle, arginine is re-synthesised from citrulline. The NOS catalysed formation of L-citrulline and NO proceeds in two steps, whereby the product stoichiometry of L-citrulline and NO is 1:1. Thus, the conversion of L-arginine to **L-citrulline can be used as a surrogate marker for the NO synthesis.**

Pathologic high levels of citrulline serve as an indicator of nitrosative stress.

### Indications

- Estimation of NOS activity (NO production)
- Detection of nitrosative stress due to an enhanced synthesis of inducible nitric oxide (iNO)

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6600	PREC	Precipitation reagent	1 x 20 ml
K 6600	STD	Standard concentrate (40 mM/l L-citrulline)	1 x 50 µl
K 6600	STDBUF	Standard dilution buffer	1 x 20 ml
K 6600	SOL A	Solution A	1 x 10 ml

Cat. No.	Label	Kit components	Quantity
K 6600	SOL B	Solution B	1 x 40 ml
K 6600	PLATE	Microtiter plate	2 x
K 6600	SAMDIL	Sample dilution buffer, lyophilised	4 x 2 ml
K 6600	CTRL1	Control 1, ready-to-use	4 x 200 µl
K 6600	CTRL2	Control 2, ready-to-use	4 x 200 µl

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1 000 µl single-use tips
- Foil to cover the microtiter plate
- Multi-channel pipets or repeater pipets
- Centrifuge, 3 000 g
- Heated incubator at 90 °C
- Water bath at 37 °C
- Metal frame for the microtiter plate strips
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 540 nm

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Reconstitute lyophilised **sample dilution buffer (SAMDIL)** with **2 ml of ultrapure water**, mix gently, allow the vial content to dissolve for 5 minutes at room temperature and mix again. **Sample dilution buffer** (reconstituted SAMDIL) **is not stable and cannot be stored.**

- **Standard curve preparation**

Prepare from the L-citrulline **stock solution (STD)** a standard curve according to the following scheme:

- **Standard 1** (400 µM/l): dilute L-citrulline stock solution (STD, 40 mM/l) **1:100** with standard dilution buffer (STDBUF), e.g. **10 µl** stock solution (40 mM/l) + **990 µl** STDBUF, mix well.
  - **Standard 2** (200 µM/l): Standard 1 1:2 diluted with STDBUF, e.g. **250 µl** standard 1 (40 mM/l) + **250 µl** STDBUF, mix well.
  - **Standard 3** (100 µM/l): Standard 2 1:2 diluted with STDBUF
  - **Standard 4** (50 µM/l): Standard 3 1:2 diluted with STDBUF
  - **Standard 5** (25 µM/l): Standard 4 1:2 diluted with STDBUF
  - **Standard 6** (12.5 µM/l): Standard 5 1:2 diluted with STDBUF
  - **Standard 7** (6.25 µM/l): Standard 6 1:2 diluted with STDBUF
  - For **standard 8** standard dilution buffer is used.
- **Preparation of the colour solution**

Mix **1 part solution A (SOL A)** with **3 parts solution B (SOL B)**. Prepare fresh colour solution for each assay, because it is stable only for around 30 minutes. Store SOL A and SOL B at 2–8 °C and bring it to room temperature before use.

- **L-Citrulline stock solution (STD), controls (CTRL) and sample dilution buffer (SAMDIL)** should be stored at **-20 °C** before use. They can be freeze-thawed up to four times.
- All other test reagents are ready-to-use. Test reagents can be used until the expiry date (see label) when stored at **2–8 °C**.

## 6. SPECIMEN COLLECTION AND PREPARATION

1.	Pipet <b>500 µl of sample</b> in 1.5 ml reaction vial
2.	Add <b>100 µl of reconstituted sample dilution buffer (SAMDIL)</b> to the sample
3.	Mix well
4.	Incubate for <b>1 h at 37 °C</b>
5.	Add <b>150 µl of cold (2–8 °C) precipitation reagent (PREC)</b>

6.	Mix well
7.	Incubate for <b>30 min at 2–8 °C</b>
8.	<b>Centrifuge</b> at 3 000 g for <b>10 min</b>
9.	Use the <b>supernatant</b> in the test

## 7. ASSAY PROCEDURE

### *Principle of the test*

After a sample pre-treating to eliminate the interference of other substances, a development solution composed of two components is added to the sample. The colour changes to intensive red due to the reaction of L-citrulline with DAMO. The interference of reaction byproducts is reduced by TSC-treating.

The colour intensity is proportional to the analyte concentration. The absorbance is measured at 540 nm. The concentration of the samples is estimated using a standard curve. **In order to eliminate the effect of the sample matrix on the absorption, an individual sample blank should be run. The obtained blank value is subtracted from the sample result.**

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of standards/controls/samples/blank for sample on a protocol sheet.

We recommend to carry out the tests in duplicate.

It is recommended to switch on the oven at 90 °C and place the metal frame for the microtiter plate modules in it before the start of the test procedure.

### **Note**

Standards, controls and samples should be pipetted without air bubbles.

1.	Loosen the strips of the microtiter plate (PLATE), so that they can be easily taken out. Pipet <b>2 x 60 µl standards/controls</b> into the into the respective wells (2 wells per standard/control; 60 µl into each).
2.	Pipette <b>4 x 60 µl of the prepared sample</b> into the respective wells (4 wells per sample; 60 µl into each).

3.	Add <b>200 µl of colour solution</b> into <b>each well for standards/controls</b> and into <b>2 of the sample wells</b> .
4.	Add <b>200 µl of SOL B</b> into the <b>2 remaining wells for the individual sample blank</b> (these without colour solution).
5.	<b>Cover</b> microtiter plate strips, take them out and bring them in a <b>metal holder pre-heated to 90 °C</b> .
6.	Incubate at <b>90 °C for 15 minutes</b> .
7.	Take the microtiter plate strips <b>out of the heater</b> and place them in the original holder.
8.	Let the strips <b>cool down</b> to room temperature <b>for 10 minutes</b> (the samples are stable for ~ 30 minutes).
9.	Read <b>absorption</b> with an microtiter plate reader at <b>540 nm</b> .

## 8. RESULTS

The final L-citrulline concentration in µmol/l is calculated as a difference between the sample concentration with the colour solution and the concentration of the sample blank (sample with SOL B) multiplied by 1.5.

$$\text{L-citrulline } [\mu\text{mol/l}] = ([\text{measured content}_{\text{sample}}] - [\text{measured content}_{\text{blank}}]) * 1.5$$

## 9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### Reference range

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### Precision and reproducibility

#### Intra-Assay (n = 12)

Sample	Citrullin [ $\mu\text{mol/l}$ ]	Standard deviation (SD) [%]
1	32.9	1.3

#### Inter-Assay (n = 7)

Sample	Citrullin [ $\mu\text{mol/l}$ ]	Standard deviation (SD) [%]
1	54.7	3.0
2	48.9	3.7

### Dilution recovery

Linearity of the test was determined by diluting a patient sample.

Dilution	Measured [ $\mu\text{mol/l}$ ]	Expected [ $\mu\text{mol/l}$ ]	Recovery [%]
original	100	100	100.0
1:1.2	87.4	83.3	104.9
1:1.5	71.2	66.6	106.9
1:3.0	37.2	33.3	111.7

### Analytical sensitivity

The Zero-standard was measured 20 times. The detection limit was set as  $B_0 + 2 \text{ SD}$ .

Sample	Citrullin mean value [OD]	2 standard deviations (2 x SD)[%]	Detection limit [ $\mu\text{mol/l}$ ]
Zero-standard	1.1	0.2	1.5

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Solution B (SOL B) contains a strong acid and must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immunodiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Knipp, M, and M Vasák. 2000. "A Colorimetric 96-Well Microtiter Plate Assay for the Determination of Enzymatically Formed Citrulline." *Analytical Biochemistry* **286** (2) (November 15): 257–64. doi:10.1006/abio.2000.4805.
2. Kuklinski, B. 2005. "Zur Praxisrelevanz von Nitrosativem Stress." *Umwelt-Medizin-Gesellschaft* **18** (2): 89–172.
3. Wanchu, A, M Khullar, K Sud, V Sakhuja, K Thennarasu, A Sud, and P Bambery. 2001. "Serum and Urine Nitrite and Citrulline Levels among Patients with Systemic Lupus Erythematosus: A Possible Addition to Activity Parameters?" *Journal of Clinical Rheumatology : Practical Reports on Rheumatic & Musculoskeletal Diseases* **7** (1) (February): 10–5; discussion 15.
4. Yu, Jin-Ju, and Suk-Heung Oh. 2010. "Isolation and Characterization of Lactic Acid Bacteria Strains with Ornithine Producing Capacity from Natural Sea Salt." *Journal of Microbiology (Seoul, Korea)* **48** (4) (August): 467–72. doi:10.1007/s12275-010-0204-9.

### Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)