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Leerwert (*limit of blank, LoB*) 6,4 ng/ml

Nachweisgrenze (*limit of detection, LoD*) 8,1 ng/ml

Bestimmungsgrenze (*limit of quantitation, LoQ*) 9,0 ng/ml

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## 15. LITERATUR

1. Gershon MD. 5-Hydroxytryptamine (serotonin) in the gastrointestinal tract. *Curr Opin Endocrinol Diabetes Obes.* 2013 Feb; 20(1):14-21. doi: 10.1097/MED.0b013e32835bc703
2. Farzaei MH, Bahramsoltani R, Abdollahi M, Rahimi R. The Role of Visceral Hypersensitivity in Irritable Bowel Syndrome: Pharmacological Targets and Novel Treatments. *J Neurogastroenterol Motil.* 2016, 22, 558–74. doi: 10.5056/jnm16001

3. Coates MD, et al. Molecular defects in mucosal serotonin content and decreased serotonin reuptake transporter in ulcerative colitis and irritable bowel syndrome. *Gastroenterology*. 2004, 126, 1657–1664
4. Yazar A, Buyukafpar K, Polat G, et al. The urinary 5-hydroxyindole acetic acid and plasmanitric oxide levels in irritable bowel syndrome: a preliminary study. *Scott Med J*. 2005, 50, 27–29
5. Dunlop SP, Coleman NS, Blackshaw E, et al. Abnormalities of 5-hydroxytryptamine metabolism in irritable bowel syndrome. *Clin Gastroenterol Hepatol*. 2005, 3(4), 349-357. doi: 10.1016/s1542-3565(04)00726-8
6. Jin DC, et al. Regulation of the serotonin transporter in the pathogenesis of irritable bowel syndrome. *World J Gastroenterol*. 2016, 22, 8137–8148
7. Yu FY, Huang SG, Zhang HY, Ye H, Chi HG, Zou Y, et al. Comparison of 5-hydroxytryptophan signaling pathway characteristics in diarrhea-predominant irritable bowel syndrome and ulcerative colitis. *World J Gastroenterol*. 2016, 22, 3451–9
8. Fu R, Chen M, Chen Y, Mao G, Liu S. Expression and clinical significance of 5-HT and 5-HT3R in the intestinal mucosa of patient with diarrhea-type irritable bowel syndrome. *Exp Ther Med*. 2019, 13, 3077-3082. doi: 10.3892/etm.2019.7297
9. Cremon C, Carini G, Wang B, Vasina V, Cogliandro RF, De Giorgio R, Stanghellini V, Grundy D, Tonini M, De Ponti F, et al. Intestinal serotonin release, sensory neuron activation, and abdominal pain in irritable bowel syndrome. *Am. J. Gastroenterol*. 2011, 106, 1290–1298.
10. Luo M, Zhuang X, Tian Z, Xiong L. Alterations in short-chain fatty acids and serotonin in irritable bowel syndrome: a systematic review and meta-analysis. *BMC Gastroenterol*. 2021 Jan 6;21(1):14. doi: 10.1186/s12876-020-01577-5

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
*For professional use only*

# **IDK® Serotonin ELISA**

***For the in vitro determination of serotonin (5-HT) in stool***

Valid from 2025-02-04

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel: + 49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

**Safety information**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit.

Important safety information for this product can be found in the chapter WARNINGS AND PRECAUTIONS.

# Table of Contents

<b>1. INTENDED USE</b>	<b>20</b>
<b>2. INTRODUCTION</b>	<b>20</b>
<b>3. MATERIAL SUPPLIED</b>	<b>20</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>21</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>21</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>22</b>
<i>Stability and storage of samples</i>	22
<i>Extraction of the stool samples</i>	22
<b>7. ASSAY PROCEDURE</b>	<b>23</b>
<i>Principle of the test</i>	23
<i>Derivatisation procedure</i>	24
<i>Test procedure</i>	25
<b>8. RESULTS</b>	<b>26</b>
<b>9. LIMITATIONS</b>	<b>26</b>
<i>Biotin interference</i>	27
<b>10. QUALITY CONTROL</b>	<b>27</b>
<i>Reference Range</i>	27
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>27</b>
<i>Accuracy – Precision</i>	27
<i>Accuracy – Trueness</i>	28
<i>Linearity</i>	28
<i>Analytical sensitivity</i>	29
<i>Analytical specificity</i>	29
<b>12. WARNINGS AND PRECAUTIONS</b>	<b>30</b>
<b>13. TECHNICAL HINTS</b>	<b>30</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>31</b>
<b>15. REFERENCES</b>	<b>31</b>
<b>16. SYMBOLS</b>	<b>32</b>

## 1. INTENDED USE

The Serotonin assay K 6881 is an enzyme-linked immunosorbent assay (ELISA) intended for professional laboratory users for the quantitative detection of the free Serotonin concentration in stool from patients of any age and gender.

The assay is an *in vitro* diagnostic medical device, which can be used manually or by an automated platform.

The detection of a decreased or elevated gastrointestinal Serotonin concentration in the stool of a patient is used as a diagnostic aid to inform the treating therapist about the necessity to prescribe appropriate dietetic interventions such as the intake of probiotics or the recommendation of a tryptophan-rich diet.

## 2. INTRODUCTION

Serotonin (5-hydroxytryptamine, 5-HT) is an important neurotransmitter in the two-way communication in the brain-gut axis. In the human body, this molecule is mainly produced (up to 95%) by the enterochromaffin cells of the intestinal mucosa. 5-HT is closely related to IBS, as it affects directly the gut motility, which has been found in numerous studies<sup>[1]</sup>. It is known, that IBS is characterised by aberrant gut motility and visceral hypersensitivity that collectively contribute to seizures and increased pain perception during digestion<sup>[2]</sup>. In this context, enteric 5-HT appears to be one of the key molecules in IBS pathology. It is known, that decreased plasma 5-HT levels are reported in patients with low motility and chronic functional constipation<sup>[3;4;5]</sup>. Vice-versa, elevated 5-HT plasma or serum levels and high spontaneous 5-HT release in mucosa are observed in patients with high gut motility or diarrhea<sup>[6;7;8;9;10]</sup>. In summary, two different kinds of situations can be identified. Constipation patients are characterized by low 5-HT levels. Diarrhea patients are characterized by high 5-HT levels.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6881	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0015	COPLATE	Plate for derivatisation	12 x 8 wells
K 6881	STD	Standards, ready-to-use (0, 10, 30, 100, 250, 750 ng/ml)	6 x 500 µl
K 6881	CTRL 1	Control, ready-to-use (see specification for range)	1 x 500 µl

K 6881	CTRL 2	Control, ready-to-use (see specification for range)	1 x 500 µl
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10 x	2 x 100 ml
K 7999.100	IDK® Amino Extract	Extraction buffer <i>IDK® Amino Extract</i> , ready-to-use	1 x 100 ml
K 6881	AB	Serotonin antibody, peroxidase- labelled, ready-to-use	1 x 6 ml
K 6881	REABUF	Reaction buffer, ready-to-use	1 x 70 ml
K 6881	DER	Derivatisation reagent, lyophilised	1 x 100 mg
K 0008.07	DMSO	Dimethylsulfoxide (DMSO)	1 x 7 ml
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipets and 10-1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml

ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** can be used until the expiry date stated on the label when stored at **2-8 °C**. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.

- Store **standards and controls (STD/CTRL)** frozen at **-20 °C**. They are stable at -20 °C until the expiry date stated on the label. Thaw before use in the test and mix well. Re-freeze standards and controls after use.
- **DMSO** crystallises at 2-8 °C. Before use, bring to room temperature to dissolve the crystals.
- The **lyophilised derivatisation reagent (DER)** can be used until the expiry date stated on the label when stored at **2-8 °C**. Bring to room temperature before opening and reconstitute the DER (100 mg) with **6 ml DMSO**. Allow to dissolve for 10 minutes and mix thoroughly with a vortex-mixer. The **derivatisation reagent** (reconstituted DER) can be stored at **2-8 °C for 2 months**. Bring to room temperature before reuse. Please note: DMSO attacks all plastics but not polypropylene products and laboratory glass.
- All other test reagents are ready-to-use. Test reagents can be used until the expiry date stated on the label when stored at **2-8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *Stability and storage of samples*

**Stool extract** is stable for 7 days at room temperature or at 2-8 °C. For longer storage keep frozen at -20 °C.

### *Extraction of the stool samples*

We recommend using the **stool sample preparation system filled with IDK® Amino extract** (extraction buffer), **Cat No. K 7999**.

The extraction buffer IDK® Amino Extract is ready-to-use. We recommend the following sample preparation:

#### ***Stool sample tube – Instructions for use***

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer:

**Stool sample tube with 0.75 ml buffer:**

Applied amount of stool:	15 mg
Buffer volume (IDK® Amino Extract):	0.75 ml
Dilution factor:	1:50

Please follow the instructions for the preparation of stool samples using the stool sample tubes as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenization using an applicator, inoculation loop or similar device.
- b) Unscrew the tube (yellow part of cap) to open. Insert the yellow dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place the dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped of, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- c) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with buffer for ~10 minutes improves the result.
- d) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- e) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

**Dilution: 1:50**

To **25 µl** of this suspension a derivatisation reagent is added (see derivatisation procedure).

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of serotonin in stool. The assay is based on the method of competitive enzyme linked immunoassays.

The sample preparation includes the addition of a derivatisation reagent for serotonin derivatisation. Afterwards, the treated samples and a peroxidase-

conjugated polyclonal serotonin antibody are incubated in wells of a microtiter plate coated with serotonin derivative (tracer). During the incubation period, the target serotonin in the sample competes with the tracer, immobilised on the wall of the microtiter wells, for the binding of the polyclonal antibodies.

After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The colour changes from blue to yellow, and the absorbance is measured in the photometer at 450 nm. The intensity of the yellow colour is inverse proportional to the serotonin concentration in the sample; this means, high serotonin concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Serotonin, present in the patient samples, is determined directly from this curve.

### *Derivatisation procedure*

Bring **all reagents and samples to room temperature** (15-30 °C) and mix well.

Derivatisation of standards, controls and extracted stool samples is carried out in reaction vials (e.g. 1.5 ml polypropylene vials). Alternatively, the derivatisation can be carried out in the wells of the COPLATE.

We recommend preparing one derivatisation per standard, control and sample and transferring it in duplicate determinations into the wells of the microtiter plate.

1.	Add <b>25 µl standard</b> (STD)/ <b>control</b> (CTRL)/ <b>extracted sample</b> into the respective vials, or into the wells of the COPLATE.
2.	Add <b>250 µl reaction buffer</b> (REABUF) into each vial (STD, CTRL, sample), or into each well of the COPLATE.
3.	Add <b>25 µl derivatisation reagent</b> into each vial (STD, CTRL, sample) and <b>mix thoroughly</b> by repeated inversion or several seconds on a vortex mixer. Incubate for <b>30 min</b> at room temperature (15-30 °C) on a <b>horizontal shaker</b> .  <i>Alternatively:</i> Add <b>25 µl derivatisation reagent</b> into each well (STD, CTRL, sample) of the COPLATE and incubate immediately on a <b>horizontal shaker</b> for <b>30 min</b> at room temperature (15-30 °C).

2 x 50 µl of the derivatised standards, controls and samples are used in the ELISA as duplicates.

## Test procedure

Mark the positions of standards/controls/samples in duplicate on a protocol sheet. Take as many microtiter strips (PLATE) as needed from the kit. Store unused strips covered with foil at 2-8 °C. Strips are stable until expiry date stated on the label.

4.	For the analysis in duplicate take <b>2 x 50 µl</b> of the <b>derivatised standards/controls/samples</b> out of the vials, or the COPLATE, and add into the respective wells of the microtiter plate.
5.	Add <b>50 µl serotonin antibody</b> (AB) into each well of the microtiter plate.
6.	Cover the plate tightly with foil and incubate for <b>1 hour at room temperature</b> (15-30 °C) on a <b>horizontal shaker</b> .
7.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8.	Add <b>100 µl substrate</b> (SUB) into each well.
9.	Incubate for <b>12-18 min*</b> at room temperature (15-30 °C) in the <b>dark</b> .
10.	Add <b>100 µl stop solution</b> (STOP) into each well and mix well.
11.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm (690 nm) as a reference.

\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. Attention should be paid to the following points:

Derivatisation of the samples is possible in 30 minutes on a horizontal shaker or, after thorough mixing, in 90 minutes without shaking in an automated processor. In automated processing, the volumes of samples and diluents may be scaled while maintaining the respective dilutions. Device-specific minimum and maximum volumes must be taken into account. The homogeneity of the resulting dilution must be ensured.

For further details please contact your supplier or Immundiagnostik AG.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the program used, the duplicate values should be evaluated manually.

### Stool samples

The obtained results have to be multiplied by the **dilution factor of 50** to get the actual concentrations.

In case another dilution factor has been used, multiply the obtained result by the dilution factor used.

## 9. LIMITATIONS

Samples with concentrations above the measurement range can be diluted with extraction buffer (IDK® Amino Extract) and re-assayed. Please consider this higher dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

The lower limit of the measurement range can be calculated as:

$$LoB \times \text{sample dilution factor to be used}$$

LoB see chapter "Performance Characteristics".

### *Biotin interference*

Samples containing a biotin concentration of  $\leq 1200$  ng/ml show a change of the results of  $< 25$  %. Higher concentrations of biotin can lead to false results. Patients taking  $> 5$  mg biotin per day should wait at least 24 hours after taking biotin to have their samples collected. Results of patients taking biotin supplements or receiving a high-dose biotin therapy should generally be interpreted along with the total clinical picture.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

### *Reference Range*

#### **1 g stool is equivalent to 1 ml.**

Based on in-house studies with samples of apparently healthy persons, a median of 1830 ng/g was determined ( $n = 40$ ). The 10<sup>th</sup> and 90<sup>th</sup> percentile were 734 ng/g and 2434 ng/g, respectively, resulting in the following normal range:

Serotonin concentration in stool: 734 - 2434ng/g

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n = 16**

The repeatability was assessed with 4 samples under constant parameters (same operator, measurement system, day and kit lot).

sample	mean value [ng/ml]	CV [%]
1	4354	4,5
2	1229	6,1
3	1059	9,8
4	829	8,9

### Reproducibility (Inter-Assay); n = 12

The reproducibility was assessed with 3 samples under varying parameters (different operators, measurement systems, days and kit lots).

sample	mean value [ng/ml]	CV [%]
1	5180	8.3
2	1238	8.9
3	676	12.7

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, three high level samples (spikes) were added to 3 low level samples. The expected values result from the mean of sample and spike.

sample [ng/ml]	spike [ng/ml]	expected [ng/ml]	obtained [ng/ml]	recovery [%]
1273	5553	3413	3564	104.4
	3644	2459	2442	99.3
	2627	1950	1872	96.0
774	5553	3164	3313	104.7
	3644	2209	2348	106.3
	2627	1700	1634	96.1
807	5553	3180	3560	111.9
	3644	2226	2038	91.6
	2627	1717	1744	101.6

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed with a serial dilution of 3 extracted stool samples.

For Serotonin in stool, the method has been demonstrated to be linear from 7.0 to 89.9 ng/ml based on the standard curve without considering possibly used sample

dilution factors (352 to 4495 ng/ml when considering the dilution factor of 50), showing a recovery rate of 88.7 – 114.8 % in this interval.

sample [ng/ml]	dilution	expected [ng/ml]	obtained [ng/ml]	recovery [%]
A	1:50		89.9	
	1:75	59.9	54.5	91.0
	1:100	44.9	44.6	99.1
	1:150	30.0	31.8	106.0
	1:200	22.5	19.9	88.7
	1:300	15.0	14.4	96.2
	1:400	11.2	12.5	111.6
B	1:50		42.2	
	1:75	28.1	29.5	105.0
	1:100	21.1	22.2	105.0
	1:150	14.1	14.3	102.0
	1:200	10.5	11.3	107.0
	1:300	7.0	8.1	114.8
C	1:50		30.1	
	1:75	20.1	20.1	99.9
	1:100	15.1	16.9	112.1
	1:150	10.0	11.2	111.8
	1:200	7.5	8.4	111.3

### *Analytical sensitivity*

The following values have been estimated based on the concentrations of the standard curve without considering possibly used sample dilution factors.

Limit of blank, LoB 6.4 ng/ml

Limit of detection, LoD 8.1 ng/ml

Limit of quantitation, LoQ 9.0 ng/ml

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20 % CV.

### *Analytical specificity*

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to serotonin. The specificity is calculated in percent in relation to the serotonin-binding activity:

3-Indoleacrylic acid < 0.04%

Indole-3-pyruvic acid	< 0.03%
3-Indoleacetic acid	< 0.05 %
5-Methoxytryptophol	< 0.15 %
L-5-OH-tryptophan	< 0.21 %

## 12. WARNINGS AND PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are harmful to health and the environment. Substrates for enzymatic colour reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact. **Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.

- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE


- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDK® is a trade mark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Gershon MD. 5-Hydroxytryptamine (serotonin) in the gastrointestinal tract. *Curr Opin Endocrinol Diabetes Obes.* 2013 Feb; 20(1):14-21. doi: 10.1097/MED.0b013e32835bc703
2. Farzaei MH, Bahramsoltani R, Abdollahi M, Rahimi R. The Role of Visceral Hypersensitivity in Irritable Bowel Syndrome: Pharmacological Targets and Novel Treatments. *J Neurogastroenterol Motil.* 2016, 22, 558–74. doi: 10.5056/jnm16001
3. Coates MD, et al. Molecular defects in mucosal serotonin content and decreased serotonin reuptake transporter in ulcerative colitis and irritable bowel syndrome. *Gastroenterology.* 2004, 126, 1657–1664
4. Yazar A, Buyukafpar K, Polat G, et al. The urinary 5-hydroxyindole acetic acid and plasmanitric oxide levels in irritable bowel syndrome: a preliminary study. *Scott Med J.* 2005, 50, 27–29
5. Dunlop SP, Coleman NS, Blackshaw E, et al. Abnormalities of 5-hydroxytryptamine metabolism in irritable bowel syndrome. *Clin Gastroenterol Hepatol.* 2005, 3(4), 349-357. doi: 10.1016/s1542-3565(04)00726-8
6. Jin DC, et al. Regulation of the serotonin transporter in the pathogenesis of irritable bowel syndrome. *World J Gastroenterol.* 2016, 22, 8137–8148

7. Yu FY, Huang SG, Zhang HY, Ye H, Chi HG, Zou Y, et al. Comparison of 5-hydroxytryptophan signaling pathway characteristics in diarrhea-predominant irritable bowel syndrome and ulcerative colitis. *World J Gastroenterol.* 2016, 22, 3451–9
8. Fu R, Chen M, Chen Y, Mao G, Liu S. Expression and clinical significance of 5-HT and 5-HT3R in the intestinal mucosa of patient with diarrhea-type irritable bowel syndrome. *Exp Ther Med.* 2019, 13, 3077-3082. doi: 10.3892/etm.2019.7297
9. Cremon C, Carini G, Wang B, Vasina V, Cogliandro RF, De Giorgio R, Stanghellini V, Grundy D, Tonini M, De Ponti F, et al. Intestinal serotonin release, sensory neuron activation, and abdominal pain in irritable bowel syndrome. *Am. J. Gastroenterol.* 2011, 106, 1290–1298.
10. Luo M, Zhuang X, Tian Z, Xiong L. Alterations in short-chain fatty acids and serotonin in irritable bowel syndrome: a systematic review and meta-analysis. *BMC Gastroenterol.* 2021 Jan 6;21(1):14. doi: 10.1186/s12876-020-01577-5

## 16. SYMBOLS

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin