IDK[®] S100A12 ELISA

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15. LITERATUR

1. Hofmann MA, Drury S; Fu C, Qu W, Taguchi A, Lu Y, Avila C, Kambham N, Bierhaus A, Nawroth P, Neurath MF, Slattery T, Beach, D, McClary J, Nagashima M, Morser J, Stern D, Schmidt AM. RAGE mediates a novel proinflammatory axis: a central cell surface receptor for S100/calgranulin polypeptides. *Cell*. **97**: 889-901, 1999
2. Karl J, Wild N, Tacke M, Andres H, Garczarek U, Rollinger W, Zolg W. Improved diagnosis of colorectal cancer using a combination of fecal occult blood and novel fecal protein markers. *Clin Gastroenterol Hepatol*. 2008 Oct;**6**(10):1122-8.
3. Foell D, Wittkowski H, Ren Z, Turton J, Pang G, Daebritz J, Ehrchen J, Heidemann J, Borody T, Roth J, Clancy R. Phagocyte-specific S100 proteins are released from affected mucosa and promote immune responses during inflammatory bowel disease. *J Pathol*. 2008 Oct;**216**(2):183-92.

4. Sidler MA, Leach ST, Day AS. Fecal S100A12 and fecal calprotectin as noninvasive markers for inflammatory bowel disease in children. *Inflamm Bowel Dis*. 2008 Mar;**14**(3):359-66.
5. Kaiser T, Langhorst J, Wittkowski H, Becker K, Friedrich AW, Rueffer A, Dobos GJ, Roth J, Foell D. Faecal S100A12 as a non-invasive marker distinguishing inflammatory bowel disease from irritable bowel syndrome. *Gut* 2007;**56**:1706-1713
6. de Jong Naomi SH, Leach ST, Day AS. Fecal S100A12: A Novel Noninvasive Marker in Children with Crohn`s Disease. *Inflamm Bowel Dis*. July 2006;**12**(7):566–572

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
IDK[®] S100A12 ELISA

For the in vitro determination of S100A12 in stool

Valid from 2022-01-04

REF K 6938

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. STORAGE AND PREPARATION OF REAGENTS	19
6. STORAGE AND PREPARATION OF SAMPLES	20
<i>Storage of stool samples</i>	20
<i>Extraction of the stool samples</i>	20
<i>Dilution of stool extract</i>	21
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Test procedure</i>	21
8. RESULTS	23
9. LIMITATIONS	23
10. QUALITY CONTROL	24
11. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Spiking Recovery</i>	24
<i>Dilution recovery</i>	25
<i>Analytical Sensitivity</i>	25
12. PRECAUTIONS	25
13. TECHNICAL HINTS	26
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
15. REFERENCES	27

1. INTENDED USE

This IDK® assay is an enzyme immunoassay intended for the quantitative determination of S100A12 in stool. For *in vitro* diagnostic use only.

2. INTRODUCTION

Alternative names: Calgranulin C, EN-RAGE

Members of the S100 protein family are low molecular mass acidic proteins characterised by cell-type-specific expression and the presence of 2 EF-hand calcium-binding domains. The S100 proteins or calgranulins are expressed in activated granulocytes under conditions of chronic inflammation. Hofmann et al. (1999) reported that RAGE is a central cell surface receptor for S100A12, which they referred to as EN-RAGE (Extracellular Newly identified RAGE-binding protein), and related members of the S100/calgranulin super family. Interaction of S100A12 with cellular RAGE on endothelium, mononuclear phagocytes, and lymphocytes triggered cellular activation, with generation of key pro-inflammatory mediators. In murine models, blockade of S100A12 quenched delayed-type hypersensitivity and inflammatory colitis by arresting activation of central signaling pathways and expression of inflammatory mediators.

Kaiser et al. (2007) determined the S100A12 concentration in stool samples of 171 patients suffering from infectious gastroenteritis, Crohn's disease, ulcerative colitis or irritable bowel syndrome and 24 healthy controls. Faecal S100A12 correlated with the degree of intestinal inflammation. Patients with inflammatory bowel disease had significantly higher S100A12 levels compared to these with irritable bowel syndrome or healthy controls. The S100A12 content in stool correlated with the degree of intestinal inflammation.

Faecal S100A12 is a novel non-invasive marker distinguishing inflammatory bowel disease (IBD) from irritable bowel syndrome (IBS) or healthy individuals with a high sensitivity and specificity. S100A12 stool concentrations can be used as an indicator of inflammation activity in Morbus Crohn's disease and Colitis ulcerosa. Since S100A12 is secreted by activated neutrophils and acts as a pro-inflammatory mediator, it can be used as marker for diagnosis and therapeutic control of inflammatory bowel diseases.

Indications

- Marker for acute inflammation
- Estimation of gastrointestinal inflammation degree
- Parameter for monitoring Morbus Crohn's disease and Colitis ulcerosa

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6938	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 6938	EXBUF	Extraction buffer concentrate, 2.5x	1 x 100 ml
K 6938	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 100 ml
K 6938	STD	S100A12 standards, lyophilised (0; 0.66; 2.0; 6.0; 18.0; 54 ng/ml)	4 x 6 vials
K 6938	CTRL1	Control, lyophilised (see specification for range)	4 x 1 vial
K 6938	CTRL2	Control, lyophilised (see specification for range)	4 x 1 vial
K 6938	CONJ	Conjugate, peroxidase-labelled, concentrate, 100x	1 x 200 µl
K 0002.15	SUB	Substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Stool sample application system such as cat. no.: K 6998SAS
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37°C. The **WASHBUF** is stable at **2–8°C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8°C for 1 month**.
- **Preparation of the extraction buffer:** The **extraction buffer concentrate (EXBUF)** has to be diluted with ultrapure water **1:2.5** before use (100 ml EXBUF + 150 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals must be redissolved at 37°C in a water bath. The **EXBUF** is stable at **2–8°C** until the expiry date stated on the label. Extraction buffer (1:2.5 diluted EXBUF) can be stored in a closed flask **for 3 months at 2–8°C**.
- The **lyophilised standards (STD)** and **controls (CTRL)** are stable at **2–8°C** until the expiry date stated on the label. **Reconstitution** details are given in the **specification data sheet**. **Standards and controls** (reconstituted STD and CTRL) **are not stable and cannot be stored**.
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in **wash buffer** (100 µl CONJ + 10 ml wash buffer). The **CONJ** is stable at **2–8°C** until the expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored**.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8°C**.

6. STORAGE AND PREPARATION OF SAMPLES

Storage of stool samples

Raw stool can be stored up to 7 days at 4 °C^[6] or for longer storage at -20 °C.

Extraction of the stool samples

Extraction buffer (1:2.5 diluted *EXBUF*) is used as a **sample extraction buffer**. We recommend the following sample preparation:

Stool Sample Application System (SAS) (Cat. No.: K 6998SAS)

Stool sample tube – Instructions for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

SAS with 0.75 ml sample extraction buffer:

Applied amount of stool: 15 mg

Buffer Volume: 0.75 ml

Dilution Factor: 1:50

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Fill the **empty stool sample tube** with 0.75 ml **sample extraction buffer** (1:2.5 diluted *EXBUF*) before using it with the sample. **Important:** Allow the sample extraction buffer to reach room temperature.
- c) Unscrew the tube (yellow part of cap) to open. Insert the yellow dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- d) Shake the tube well until no stool sample remains in the notches. Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with sample extraction buffer for ~ 10 minutes improves the result.

- e) Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- f) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

Dilution I: 1:50

Dilution of stool extract

The supernatant of the extraction (dilution I) is diluted **1:20** with **sample dilution buffer** (SAMPLEBUF). For example:

50 µl dilution I + **950 µl** SAMPLEBUF = 1:50 (dilution II)

This results in a final dilution of **1:1000**.

For analysis, pipette **100 µl** of the supernatant of **dilution II** per well.

We recommend to store the extracted samples (dilution I) for no longer than **24 h at 2–8 °C**. The **extracted samples are not stable and cannot be frozen** and stored for later use.

7. ASSAY PROCEDURE

Principle of the test

This ELISA is designed for the quantitative determination of S100A12. A monoclonal antibody specific for S100A12 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and the immobilised antibody binds any S100A12 present. After washing away any unbound substances, an HRP-conjugated polyclonal antibody specific for S100A12 is added to the wells. Following a wash to remove any unbound antibody HRP-conjugate, the remaining conjugate is allowed to react with the HRP substrate tetramethylbenzidine. The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured at 450 nm. The absorbance is proportional to the concentration of S100A12. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. S100A12, present in the patient samples, is determined directly from this curve.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Before use wash each well 5 x with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
2.	Add each 100 µl standards/controls/samples into the respective wells.
3.	Cover the strips and incubate for 1 hour at room temperature (15–30 °C) on a horizontal mixer* .
4.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5.	Add 100 µl conjugate (CONJ) in each well.
6.	Cover the strips and incubate for 1 hour at room temperature (15–30 °C) on a horizontal mixer* .
7.	Discard the contents of each well and wash 5 times with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add 100 µl TMB substrate (SUB) in each well.
9.	Incubate for 10–20 min** at room temperature (15–30 °C) in the dark .
10.	Add 100 µl ELISA stop solution (STOP) and mix well.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

** The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the “4 parameter algorithm”.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

Stool samples

The obtained S100A12 levels of stool samples have to be multiplied by the dilution factor of **1000**.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be further diluted and re-assayed. Please consider this higher dilution when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

analytical sensitivity × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 16)

Sample	S100A12 [ng/ml]	CV [%]
1	3718.4	7.1
2	918.3	3.5

Inter-Assay (n = 22)

Sample	S100A12 [ng/ml]	CV [%]
1	3003.64	15.4
2	1363.49	20.2

Spiking Recovery

Two samples were spiked with different S100A12 concentrations and measured using this assay (n = 2).

Sample	Spike [ng/ml]	S100A12 expected [ng/ml]	S100A12 measured [ng/ml]
A	0.56	0.892	0.803
	0.28	0.602	0.56
B	1.125	1.642	1.701
	0.56	1.077	0.98

Dilution recovery

Two patient samples were diluted and analysed. The results are shown below (n = 2):

Sample	Dilution	S100A12 expected [ng/ml]	S100A12 measured [ng/ml]
A	1:2500	8797.3	8797.3
	1:5000	4399	4239.8
	1:10 000	2199.5	1986
B	1:2500	2016.1	2016.1
	1:5000	1008	932.2
	1:10 000	504	362

Analytical Sensitivity

The Zero-standard was measured 10 times. The detection limit was set as $B_0 + 2 SD$ and estimated to be 0.016 ng/ml.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

Warning: Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still should be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Hofmann MA, Drury S, Fu C, Qu W, Taguchi A, Lu Y, Avila C, Kambham N, Bierhaus A, Nawroth P, Neurath MF, Slattery T, Beach, D, McClary J, Nagashima M, Morser J, Stern D, Schmidt AM. RAGE mediates a novel proinflammatory axis: a central cell surface receptor for S100/calgranulin polypeptides. *Cell*. **97**: 889-901, 1999
2. Karl J, Wild N, Tacke M, Andres H, Garczarek U, Rollinger W, Zolg W. Improved diagnosis of colorectal cancer using a combination of fecal occult blood and novel fecal protein markers. *Clin Gastroenterol Hepatol*. 2008 Oct;**6**(10):1122-8.
3. Foell D, Wittkowski H, Ren Z, Turton J, Pang G, Daebritz J, Ehrchen J, Heidemann J, Borody T, Roth J, Clancy R. Phagocyte-specific S100 proteins are released from affected mucosa and promote immune responses during inflammatory bowel disease. *J Pathol*. 2008 Oct;**216**(2):183-92.
4. Sidler MA, Leach ST, Day AS. Fecal S100A12 and fecal calprotectin as noninvasive markers for inflammatory bowel disease in children. *Inflamm Bowel Dis*. 2008 Mar;**14**(3):359-66.
5. Kaiser T, Langhorst J, Wittkowski H, Becker K, Friedrich AW, Rueffer A, Dobos GJ, Roth J, Foell D. Faecal S100A12 as a non-invasive marker distinguishing inflammatory bowel disease from irritable bowel syndrome. *Gut* 2007;**56**:1706-1713
6. de Jong Naomi SH, Leach ST, Day AS. Fecal S100A12: A Novel Noninvasive Marker in Children with Crohn`s Disease. *Inflamm Bowel Dis*. July 2006;**12**(7):566-572.

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant