# IDK<sup>®</sup> S100A12 ELISA

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
## 15. LITERATUR

1. Hofmann MA, Drury S; Fu C, Qu W, Taguchi A, Lu Y, Avila C, Kambham N, Bierhaus A, Nawroth P, Neurath MF, Slattery T, Beach, D, McClary J, Nagashima M, Morser J, Stern D, Schmidt AM. RAGE mediates a novel proinflammatory axis: a central cell surface receptor for S100/calgranulin polypeptides. *Cell*. **97**: 889-901, 1999
2. Meijer, B., Gearry, R.B. & Day, A.S., 2012. The role of S100A12 as a systemic marker of inflammation. *International journal of inflammation*, **2012**, p.907078.

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
# **IDK<sup>®</sup> S100A12 ELISA**

***For the in vitro determination of S100A12 in serum***

Valid from 2022-01-04

**REF** **K 6939**

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>16</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>17</b>
<i>Sample storage</i>	17
<i>Sample dilution</i>	17
<b>7. ASSAY PROCEDURE</b>	<b>17</b>
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
<b>8. RESULTS</b>	<b>19</b>
<b>9. LIMITATIONS</b>	<b>20</b>
<b>10. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Accuracy – Precision</i>	20
<i>Accuracy – Trueness</i>	21
<i>Linearity</i>	21
<i>Analytical sensitivity</i>	22
<b>12. PRECAUTIONS</b>	<b>22</b>
<b>13. TECHNICAL HINTS</b>	<b>23</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>15. REFERENCES</b>	<b>24</b>

## 1. INTENDED USE

This IDK® assay is an enzyme immunoassay intended for the quantitative determination of S100A12 in serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Alternative names: Calgranulin C, EN-RAGE

Members of the S100 protein family are low molecular mass acidic proteins characterised by cell-type-specific expression and the presence of 2 EF-hand calcium-binding domains. The S100 proteins or calgranulins are expressed in activated granulocytes under conditions of chronic inflammation. Hofmann et al. (1999) reported that RAGE is a central cell surface receptor for S100A12, which they referred to as EN-RAGE (Extracellular Newly identified RAGE-binding protein), and related members of the S100/calgranulin super family. Interaction of S100A12 with cellular RAGE on endothelium, mononuclear phagocytes, and lymphocytes triggered cellular activation, with generation of key pro-inflammatory mediators. In murine models, blockade of S100A12 quenched delayed-type hypersensitivity and inflammatory colitis by arresting activation of central signaling pathways and expression of inflammatory mediators.

Since S100A12 is secreted by activated neutrophilic granulocytes and acts as a pro-inflammatory mediator, it can be used as marker for diagnosis and therapeutic control of various inflammatory diseases (see review by Meijer et al., 2012).

### Indications

- Marker for acute inflammation
- Estimation of inflammation degree

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 6939	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 6939	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 100 ml
K 6939	STD	S100A12 standards, lyophilised (0; 0.66; 2.0; 6.0; 18.0; 54 ng/ml)	4 x 6 vials
K 6939	CTRL1	Control, lyophilised (see specification for range)	4 x 1 vial

Cat. No.	Label	Kit components	Quantity
K 6939	CTRL2	Control, lyophilised (see specification for range)	4 x 1 vial
K 6939	CONJ	Conjugate, peroxidase-labelled, concentrate, 100x	1 x 200 µl
K 0002.15	SUB	Substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. The crystals must be redissolved at room tem-

perature or in a water bath at 37°C before dilution of the concentrate. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask **for 1 month at 2–8 °C**.

- The **lyophilised standards** (STD) and **controls** (CTRL) are stable at **2–8 °C** until the expiry date stated on the label. **Reconstitution** details are given in the **specification data sheet**. **Standards and controls** (reconstituted STD and CTRL) **are not stable and cannot be stored**.
- **Preparation of the conjugate:** Before use, the **conjugate concentrate** (CONJ) has to be diluted **1:101** in **wash buffer** (100 µl CONJ + 10ml wash buffer). The **CONJ** is stable at **2–8 °C** until the expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored**.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *Sample storage*

S100A12 is stable in undiluted serum for 9 months at -20°C as well as for 7 days at 2–8°C and room temperature.

### *Sample dilution*

The sample is diluted **1:40** with **sample dilution buffer** (SAMPLEBUF). For example:

**20 µl** sample + **780 µl** SAMPLEBUF = 1:40

For analysis, pipet **100 µl** diluted sample per well.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of S100A12. A monoclonal antibody specific for S100A12 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and the immobilised antibody binds any S100A12 present. After washing away any unbound substances, an HRP-conjugated polyclonal antibody specific for S100A12 is added to the wells. Following a wash to remove any unbound antibody HRP-conjugate, the remaining conjugate is allowed to react with the HRP substrate tetramethylbenzidine. The reaction is stopped by ad-

dition of acidic solution and absorbance of the resulting yellow product is measured at 450 nm. The absorbance is proportional to the concentration of S100A12. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. S100A12, present in the samples, is determined directly from this curve.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	<b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Add each <b>100 µl standards/controls/samples</b> into the respective wells.
3.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal mixer*</b> .
4.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5.	Add <b>100 µl conjugate</b> (CONJ) in each well.
6.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal mixer*</b> .
7.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add <b>100 µl substrate</b> (SUB) in each well.
9.	Incubate for <b>10–20 min**</b> at room temperature (15–30 °C) <b>in the dark</b> .

10.	Add <b>100 µl stop solution</b> (STOP) and mix well.
11.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### Serum samples

The obtained S100A12 levels of samples have to be multiplied by the dilution factor of **40**.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## 9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard can be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

Serum (n = 63): 203 ng/ml (95th percentile)

Based on an Immundiagnostik AG study of evidently healthy persons (n = 63), a mean value of 101 ng/ml (standard deviation: 52 ng/ml) and a median value of 89 ng/ml were estimated. The 95th percentile (203 ng/ml) should be considered as the preliminary upper reference value of the test.

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n = 80**

The repeatability was assessed with 2 serum samples under **constant** parameters (same operator, measurement system, day and kit lot).

Sample	Mean value [ng/ml]	CV [%]
1	215.46	5.8
2	421.53	5.0

### Reproducibility (Inter-Assay); n = 21

The reproducibility was assessed with 2 serum samples under **varying** parameters (different operators, measurement systems, days and kit lots).

Sample	Mean value [ng/ml]	CV [%]
1	261.38	11.4
2	1012.62	8.0

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, S100A12-spikes with known concentrations were added to 2 different serum samples.

Sample [ng/ml]	Spike [ng/ml]	Obtained [ng/ml]	Expected [ng/ml]	Recovery [%]
1.25	26.21	29.10	27.46	106.0
	13.10	14.68	14.35	102.3
	6.55	7.73	7.80	99.0
	3.28	5.01	4.53	110.7
	1.64	3.00	2.89	103.8
2.51	26.21	25.56	28.72	89.0
	13.10	13.74	15.61	88.0
	6.55	8.91	9.06	98.3
	3.28	5.79	5.79	100.1
	1.64	4.92	4.15	118.7

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A with a serial dilution of 2 different serum samples.

For S100A12 in serum, the method has been demonstrated to be linear from 24.88–3.22 ng/ml, showing a non-linear behaviour of less than  $\pm 20\%$  in this interval.

Sample	Dilution	Obtained [ng/ml]	Expected [ng/ml]	Recovery [%]
A	1:40	24.88	24.88	100.0
	1:50	19.00	19.90	95.5
	1:60	16.26	16.58	98.0
	1:70	13.78	14.21	97.0
	1:80	12.70	12.44	102.1
	1:160	5.82	6.22	93.7
B	1:40	13.89	13.89	100.0
	1:60	8.48	9.26	91.6
	1:80	6.31	6.94	90.8
	1:100	4.86	5.56	87.5
	1:120	3.78	4.63	81.6
	1:140	3.64	3.97	91.6
	1:160	3.22	3.47	92.8

### Analytical sensitivity

The following values have been estimated based on the concentrations of the standard/calibration curve without considering possibly used sample dilution factors

Limit of blank, LoB 0.157 ng/ml

Limit of detection, LoD 0.400 ng/ml

Limit of quantitation, LoQ 0.605 ng/ml

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20 % CV.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still should be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK*® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure,


which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Hofmann MA, Drury S; Fu C, Qu W, Taguchi A, Lu Y, Avila C, Kambham N, Bierhaus A, Nawroth P, Neurath MF, Slattery T, Beach, D, McClary J, Nagashima M, Morser J, Stern D, Schmidt AM. RAGE mediates a novel proinflammatory axis: a central cell surface receptor for S100/calgranulin polypeptides. *Cell*. **97**: 889-901, 1999
2. Meijer, B., Geary, R.B. & Day, A.S., 2012. The role of S100A12 as a systemic marker of inflammation. *International journal of inflammation*, **2012**, p.907078.

### Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant