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
Manual*For professional use only*

SDMA ELISA

For the in vitro determination of SDMA in human, canine and feline serum, EDTA and Li-heparin plasma

Valid from 2025-04-03

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: + 49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Safety information

The assay has to be performed exclusively according to the instructions for use enclosed with the kit.

Important safety information for this product can be found in the chapter WARNINGS AND PRECAUTIONS.

Table of Contents

1. INTENDED USE	19
2. INTRODUCTION	19
3. MATERIAL SUPPLIED	19
4. MATERIAL REQUIRED BUT NOT SUPPLIED	20
5. STORAGE AND PREPARATION OF REAGENTS	20
6. STORAGE AND PREPARATION OF SAMPLES	21
7. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Derivatisation procedure</i>	22
<i>Test procedure</i>	22
8. RESULTS	23
9. LIMITATIONS	24
<i>Biotin interference</i>	24
10. QUALITY CONTROL	25
<i>Reference range</i>	25
11. PERFORMANCE CHARACTERISTICS	25
<i>Precision and reproducibility</i>	25
<i>Spiking recovery</i>	25
<i>Dilution recovery</i>	26
<i>Analytical sensitivity</i>	28
<i>Specificity</i>	28
<i>Correlation with HPLC-MS</i>	28
12. WARNINGS AND PRECAUTIONS	29
13. TECHNICAL HINTS	29
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	30
15. REFERENCES	30
<i>General literature</i>	30
<i>Literature using Immundiagnostik SDMA ELISA</i>	31
16. SYMBOLS	31

1. INTENDED USE

This Immundiagnostik AG assay is intended for the quantitative determination of SDMA in human, canine and feline serum, EDTA and Li-heparin plasma. For *in vitro* diagnostic use only.

2. INTRODUCTION

The dosage of most drugs must be adapted in renal insufficiency, making accurate assessment of renal function a prerequisite in clinical medicine. Furthermore, even a modest decline in renal function has been recognized as a cardiovascular risk.

In clinical practice serum creatinine is typically used to assess renal function, but this serum creatinine does not increase at modest decline in renal function. Consequently, there is an ongoing search for suitable endogenous markers of renal function.

SDMA is a methylated derivative of L-arginine which is strictly eliminated by renal extraction, thus the SDMA plasma level is strongly correlated to renal function. In 18 studies with more than 2136 patients systemic SDMA concentrations correlated highly with inulin clearance, as well as with various clearance estimates combined and serum creatinine. With respect to this, SDMA exhibits properties of a reliable marker of renal dysfunction.

Moreover, there are hints that increased SDMA correlates with total sequential organ failure, indicating both renal and hepatic failure, and with an increased cardiovascular risk.

Indication

- Renal failure
- Cardiovascular risk in renal dysfunction
- Hypertension in renal dysfunction

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 7780	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 7780	STD	Standards, ready-to-use (0, 0.1, 0.3, 0.6, 1.5, 4.0 µM)	6 x 500 µl
K 7780	CTRL 1	Control, ready-to-use (see specification for range)	1 x 500 µl
K 7780	CTRL 2	Control, ready-to-use (see specification for range)	1 x 500 µl

K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 7780	AB	SDMA antibody, ready-to-use	1 x 6 ml
K 7780	CONJ	Conjugate, ready-to-use	1 x 12 ml
K 0012.15	DERBUF	Reaction buffer, ready-to-use	1 x 15 ml
K 7780	DER	Derivatisation reagent, ready-to-use	1 x 6 ml
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Calibrated precision pipets and 10-1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** can be used until the expiry date stated on the label when stored at **2-8 °C**. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.

- Store **standards and controls (STD/CTRL)** frozen at **-20 °C**. They are stable at -20 °C until the expiry date stated on the label. Thaw before use in the test and mix well. Re-freeze standards and controls after use.
- The ready-to-use **derivatisation reagent (DER)** is already dissolved in DMSO and is stable at **2-8 °C** until the expiry date stated on the label. DMSO crystallises at 2-8 °C. Before opening the DER, bring to room temperature and ensure that all crystals are dissolved. Please note: DMSO attacks plastics but not polypropylene products and laboratory glass.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2-8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Serum, EDTA and Li-heparin plasma

- Freshly collected serum, EDTA and Li-heparin plasma samples can be stored for up to 48 hours at room temperature or for 3 days at 2-8 °C. For longer storage keep samples frozen at -20 °C.
- The samples are analysed **undiluted**.
If the sample volume is less than 50 µl, we recommend a 1:2 dilution in reaction buffer (25 µl sample + 25 µl DERBUF). This dilution factor must be considered in data evaluation.
- For sample preparation, a derivatisation reagent is added (see derivatisation procedure).

7. ASSAY PROCEDURE

Principle of the test

This ELISA is designed for the quantitative determination of SDMA. The assay is based on the method of competitive enzyme linked immunoassays.

The sample preparation includes the addition of a derivatisation reagent for SDMA derivatisation. Afterwards, the treated samples and the polyclonal SDMA antiserum are incubated in wells of a microtiter plate coated with SDMA derivative (tracer). During the incubation period, the target SDMA in the sample competes with the tracer, immobilised on the wall of the microtiter wells, for the binding of the polyclonal antibodies.

During the second incubation step, a peroxidase conjugated antibody is added to detect the anti-SDMA antibodies. After washing away the unbound components,

tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The colour changes from blue to yellow and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow colour is inverse proportional to the SDMA concentration in the sample; this means high SDMA concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. SDMA, present in the patient samples, is determined directly from this curve.

Derivatisation procedure

Bring **all reagents and samples to room temperature** (15-30 °C) and mix well.

Derivatisation of standards, controls and samples is carried out in single analysis in vials (e.g. 1.5 ml polypropylene vials).

We recommend preparing one derivatisation per standard, control and sample and transferring it in duplicate determinations into the wells of the microtiter plate.

1.	Add 200 µl standard (STD), 200 µl control (CTRL) and 50 µl sample in the corresponding vials.
2.	Add 150 µl reaction buffer (DERBUF) only to the samples .
3.	Add 50 µl derivatisation reagent (DER) into each vial (STD, CTRL, sample), mix thoroughly by repeated inversion or several seconds on a vortex mixer. Incubate for 45 min at room temperature (15-30 °C) on a horizontal shaker ..
4.	Incubate for 45 min at room temperature (15-30 °C) on a horizontal shaker .

2 x 50 µl of the derivatised standards, controls and samples are used in the ELISA as duplicates.

Test procedure

Mark the positions of standards/controls/samples in duplicate on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips covered with foil at 2-8 °C. Strips are stable until expiry date stated on the label.

5.	For the analysis in duplicate, take 2 x 50 µl of the derivatised standards/controls/samples out of the vials and add into the respective wells of the microtiter plate.
6.	Add 50 µl SDMA antibody (AB) into each well of the microtiter plate.

7.	Cover the strips and incubate for 2 hours at room temperature (15-30 °C) on a horizontal shaker .
8.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
9.	Add 100 µl conjugate (CONJ) into each well.
10.	Cover the strips tightly with foil and incubate overnight at 2-8°C .
11.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
12.	Add 100 µl substrate (SUB) into each well.
13.	Incubate for 10-15 min* at room temperature (15-30 °C) in the dark .
14.	Add 100 µl stop solution (STOP) into each well and mix well.
15.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference.

* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the program used, the duplicate values should be evaluated manually.

Serum, EDTA and Li-heparin plasma

Since the sample dilution is already considered in the standard curve, the dilution factor is 1.

In case an additional dilution factor is used, multiply the obtained result by the additionally used dilution factor.

9. LIMITATIONS

Samples with concentrations above the measurement range can be diluted with reaction buffer and re-assayed. Please consider this dilution factor when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

analytical sensitivity × sample dilution factor to be used

Analytical sensitivity see chapter "Performance Characteristics".

Biotin interference

Samples containing a biotin concentration of < 40 ng/ml show a change of the results of ≤ 25 %. Higher concentrations of biotin can lead to false results. Patients taking > 5 mg biotin per day should wait at least 24 hours after taking biotin to have their samples collected. Results of patients taking biotin supplements or receiving a high-dose biotin therapy should generally be interpreted along with the total clinical picture.

10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference range

Based on internal studies with human serum samples from apparently healthy persons (n = 40), a mean value of 0.47 µmol/l was estimated. The standard deviation (SD) was 0.09 µmol/l. From mean value $\pm$ 2 SD a normal range of 0.29 – 0.65 µmol/l was estimated.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-assay (n = 10)

sample	SDMA [µmol/l]	CV [%]
1	0.27	7.5
2	0.67	4.8

Inter-assay (n = 6)

sample	SDMA [µmol/l]	CV [%]
1	0.22	9.9
2	0.63	7.4

Spiking recovery

Human

One serum sample was spiked with different SDMA concentrations and measured in this assay. The mean recovery rate was 100.4 % (n = 5).

spike [μmol/l]	expected [μmol/l]	measured [μmol/l]	recovery [%]
0		0.736	
0.5	1.236	1.220	98.7
1	1.736	1.774	102.2

Dog

Two serum samples were spiked with different SDMA concentrations. The mean recovery rate was 94.1 % (n = 2).

spike [μmol/l]	expected [μmol/l]	measured [μmol/l]	recovery [%]
0		0.526	
0.5	1.026	0.937	91.3
1	1.526	1.561	102.3
0		0.454	
0.5	0.954	0.897	94.0
1	1.454	1.288	88.6

Cat

Two serum samples were spiked with different SDMA concentrations. The mean recovery rate was 95.4 % (n = 2).

spike [μmol/l]	expected [μmol/l]	measured [μmol/l]	recovery [%]
0		0.395	
0.5	0.895	0.833	93.1
1	1.395	1.252	89.8
0		0.326	
0.5	0.826	0.818	99.0
1	1.326	1.321	99.6

Dilution recovery

Human

One spiked sample was diluted with reaction buffer. The mean recovery rate was 90.4 % (n = 5).

dilution	expected [μmol/l]	measured [μmol/l]	recovery [%]
		1.774	
1:2	0.887	0.876	98.8
1:4	0.444	0.364	82.1

Dog

Two samples were diluted with serum containing 0.126 μmol/l SDMA. The mean recovery rate was 109.0 % (n = 2).

dilution	expected [μmol/l]	measured [μmol/l]	recovery [%]
		0.526	
1:2	0.326	0.334	102.5
1:4	0.226	0.241	106.6
		0.454	
1:2	0.290	0.340	117.2
1:4	0.208	0.228	109.6

Cat

Two samples were diluted with serum containing 0.126 μmol/l SDMA. The mean recovery rate was 97.1 % (n = 2).

dilution	expected [μmol/l]	measured [μmol/l]	recovery [%]
		0.395	
1:2	0.261	0.247	94.8
1:4	0.193	0.194	100.4
		0.326	
1:2	0.226	0.231	102.2
1:4	0.176	0.160	90,9

Analytical sensitivity

The zero-standard (STD 1) was measured 40 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be $0.05 \mu\text{mol/l}$.


Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to SDMA. The specificity is calculated in percent in relation to the SDMA-binding activity.

ADMA	< 0.1 %
L-arginine	< 0.001 %

Correlation with HPLC-MS

16 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.99$.

12. WARNINGS AND PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are harmful to health and the environment. Substrates for enzymatic colour reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found on the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

Warning: Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention..

- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

General literature

1. Fleck C, Schweitzer F, Karge E, Busch M, Stein G. Serum concentrations of asymmetric (ADMA) and symmetric (SDMA) dimethylarginine in patients with chronic kidney diseases. *Clinical Chimica Acta* **336**: 1 – 12 (2003).
2. D'Apolito O, Paglia G, Tricarico F, Garofalo D, Pilotti A, Lamacchia O, Cignarelli M, Corso G. Development and validation of a fast quantitative method for plasma dimethylarginines analysis using liquid chromatography-tandem mass spectrometry. *Clinical Biochemistry* **41**: 1391 – 1395 (2008).
3. Kielstein JT, Salpeter SR, Bode-Böger SM, Cooke JP, Fliser D. Symmetric dimethylarginine (SDMA) as endogenous marker of renal function – a meta-analysis. *Nephrol. Dial. Transplant* **21**: 2446 - 2451 (2006).
4. Bode-Böger SM, Scalera F, Kielstein JT, Martens-Lobenhoffer J, Breithardt G, Fobker M, Reinecke H. Symmetrical Dimethylarginine: A new combined parameter for renal function and extent of coronary artery disease. *J. Am. Soc. Nephrol.* **17**: 1128 – 1134 (2006).

Literature using Immundiagnostik SDMA ELISA

5. Koch A, Weiskirchen R, Bruensing J, Dückers H, Buendgens L, Kunze J, Matthes M, Luedde T, Trautwein C, Tacke F. Regulation and prognostic relevance of symmetric dimethylarginine serum concentrations in critical illness and sepsis. *Mediators of Inflammation* **2013**: 413826 (2013).
6. Tenderenda-Banasiuk E, Wasilewska A, Taranta-Janusz K, Korzeniecka-Kozerska A. Asymmetric and symmetric dimethylarginine in adolescents with hyperuricemia. *Disease Markers* **35**: 407-412 (2013).

16. SYMBOLS

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin