# Intaktes Proinsulin ELISA

*Zur in-vitro-Bestimmung von intaktem Proinsulin  
in humanem Plasma*

# Intact Proinsulin ELISA

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
# Intact Proinsulin

*For the in vitro determination of intact proinsulin  
in human plasma*

Valid from 2023-11-14

**REF** K 7821

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

|  |           |
|--|-----------|
| <b>1. INTENDED USE</b>                                     | <b>19</b> |
| <b>2. INTRODUCTION</b>                                     | <b>19</b> |
| <b>3. MATERIAL SUPPLIED</b>                                | <b>19</b> |
| <b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>               | <b>20</b> |
| <b>5. STORAGE AND PREPARATION OF REAGENTS</b>              | <b>20</b> |
| <b>6. STORAGE AND PREPARATION OF SAMPLES</b>               | <b>21</b> |
| <i>Dilution of samples</i>                                 | 21        |
| <i>Sample matrices</i>                                     | 21        |
| <i>Sample collection</i>                                   | 21        |
| <i>Specimen storage</i>                                    | 21        |
| <i>Sample stability</i>                                    | 22        |
| <b>7. ASSAY PROCEDURE</b>                                  | <b>23</b> |
| <i>Principle of the test</i>                               | 23        |
| <i>Test procedure</i>                                      | 23        |
| <b>8. RESULTS</b>  | <b>25</b> |
| <b>9. LIMITATIONS</b>                                      | <b>25</b> |
| <b>10. QUALITY CONTROL</b>                                 | <b>26</b> |
| <i>Expected values</i>                                     | 26        |
| <i>Limitations</i>   | 27        |
| <i>Interfering substances</i>                              | 27        |
| <b>11. PERFORMANCE CHARACTERISTICS</b>                     | <b>27</b> |
| <i>Precision and reproducibility</i>                       | 27        |
| <i>Spiking recovery</i>                                    | 28        |
| <i>Linearity</i>   | 28        |
| <i>Cross reactivity</i>                                    | 28        |
| <i>Sensitivity</i>   | 29        |
| <i>High dose hook effect</i>                               | 29        |
| <i>Comparison with the gold standard chemoluminescence</i> | 29        |
| <b>12. PRECAUTIONS</b>                                     | <b>30</b> |
| <b>13. TECHNICAL HINTS</b>                                 | <b>30</b> |
| <b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>    | <b>31</b> |
| <b>15. REFERENCES</b>                                      | <b>31</b> |

## 1. INTENDED USE

This Immundiagnostik AG assay is an enzyme immunoassay intended for the quantitative determination of intact proinsulin in serum and plasma. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Measurements of proinsulin are used in the diagnosis and treatment of patients with type 2 diabetes

Proinsulin is a precursor molecule for insulin and is synthesised by the pancreatic  $\beta$ -cells. Under normal circumstances, virtually all proinsulin is cleaved at residues 32–33 and 65–66 to produce insulin during the formation of secretory granules. Some unmodified proinsulin is released into the circulation, though it is believed to have little or no biological activity. Increased concentrations of circulating proinsulin may occur in insulin-resistant syndromes such as type 2 diabetes and in patients with insulinoma. When used in conjunction with a highly specific insulin assay, it may provide useful information on changes in the processing of insulin in such situations.

### Indications

- Diabetes type 1 and 2
- Insulinomas, functional hypoglycaemia, hyperinsulinemia
- Kidney dysfunction, cirrhosis and hyperthyreosis
- Cardiovascular diseases

## 3. MATERIAL SUPPLIED

| Cat. No.     | Label     | Kit components   | Quantity     |
|--------------|-----------|--|--------------|
| K 7821       | PLATE     | Microtiter plate, pre-coated                                     | 12 x 8 wells |
| K 0001.C.100 | WASHBUF   | Wash buffer concentrate, 10x                                     | 2 x 100 ml   |
| K 7821       | AB        | Biotin-conjugated proinsulin antibody, ready-to-use              | 1 x 12 ml    |
| K 7821       | SAMPLEBUF | Sample buffer, ready-to-use                                      | 1 x 12 ml    |
| K 7821       | CONJ      | HRP-conjugate, ready-to-use                                      | 1 x 12 ml    |
| K 7821       | STD       | Standards, lyophilised<br>(see specification for concentrations) | 5 x 1 ml     |

| Cat. No.  | Label  | Kit components  | Quantity      |
|-----------|--------|---|---------------|
| K 7821    | CTRL A | Control, lyophilised<br>(see specification for range) | 1 x 1 ml vial |
| K 7821    | CTRL B | Control, lyophilised<br>(see specification for range) | 1 x 1 ml vial |
| K 0002.15 | SUB    | Substrate (tetramethylbenzidine),<br>ready-to-use     | 1 x 15 ml     |
| K 0003.15 | STOP   | Stop solution, ready-to-use                           | 1 x 15 ml     |

The used standards have been calibrated on the WHO reference material IRR 09/296.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

- The **lyophilised standards** (STD) and **controls** (CTRL) are stable at **2–8°C** until the expiry date stated on the label. Before use, the STD and CTRL have to be reconstituted with **1000 µl of ultrapure water**. Allow the vial content to dissolve for 5 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. **Standards and controls** (reconstituted STD and CTRL) **can be stored at 2–8°C for 1 month or at least for 8 weeks at -20°C with up to 3 freeze-thaw-cycles**.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8°C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *Dilution of samples*

No sample dilution is necessary.

### *Sample matrices*

It is recommended to use heparin or EDTA Plasma for intact proinsulin measurements. Full recovery of intact proinsulin cannot be achieved from serum samples. Do not use severely hemolysed specimens.

### *Sample collection*

Whole blood should be collected into a tube containing EDTA or heparin as anticoagulant, allowed to stand for maximal 20 min and then centrifuged for at least 15 min at 2000 to 3000 rpm. The clear cell-free supernatant should be immediately collected and used in the assay.


### *Specimen storage*

After collection, the specimens should be capped and may be stored for up to 24 hours at 2–8°C prior to assaying. Specimens held for a longer time should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.


We recommend duplicate analysis for each sample.

### Sample stability

Samples can be stored for maximal 24 hours at 2–8°C prior to assaying. For a longer storage, the samples should be frozen at -20°C.

## long term stability at -20°C

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of intact proinsulin.

The intact proinsulin assay is a two-site immunoassay, employing a specific solid phase antibody immobilised on microtiter wells and a soluble antibody labelled with biotin. The sample is incubated in the microtiter well together with a buffer and, after a wash step, the labelled antibody solution is added. After a second incubation and wash step, HRP-labelled streptavidin is added. A third incubation and wash is followed by the addition of substrate solution. Following colour development, stop reagent is added and the colour intensity measured in a 96-well ELISA reader. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Intact proinsulin, present in the samples, is determined directly from this curve.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Mark the positions of standards/samples/controls on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8°C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

|     |  |
|-----|--|
| 1.  | Add <b>50 µl sample buffer</b> (SAMPLEBUF) into each well.   |
| 2.  | Add each <b>50 µl standards/samples/controls</b> into the respective wells.  |
| 3.  | Cover the strips tightly and incubate for <b>2 hours at room temperature</b> (15–30 °C).   |
| 4.  | Discard the content of each well. Remove residual fluids by tapping the inverted microtiter plate on absorbent paper. Fill wells with <b>300 µl wash buffer</b> . Empty wells and remove residual buffer by tapping the inverted microtiter plate on absorbent paper.<br><b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the inverted microtiter plate on absorbent paper after each washing step. |
| 5.  | Add <b>100 µl biotin conjugated antibody</b> (AB) into each well.  |
| 6.  | Cover the strips tightly and incubate for <b>1 hour at room temperature</b> (15–30 °C).  |
| 7.  | Discard the content of each well. Remove residual fluids by tapping the inverted microtiter plate on absorbent paper. Fill wells with <b>300 µl wash buffer</b> . Empty wells and remove residual buffer by tapping the inverted microtiter plate on absorbent paper.<br><b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the inverted microtiter plate on absorbent paper after each washing step. |
| 8.  | Add <b>100 µl HRP-conjugate</b> (CONJ) into each well.   |
| 9.  | Cover the strips tightly and incubate for <b>30 minutes at room temperature</b> (15–30 °C).  |
| 10. | Discard the content of each well. Remove residual fluids by tapping the inverted microtiter plate on absorbent paper. Fill wells with <b>300 µl wash buffer</b> . Empty wells and remove residual buffer by tapping the inverted microtiter plate on absorbent paper.<br><b>Repeat this step 2x</b> and make sure to remove residual buffer by tapping the inverted microtiter plate on absorbent paper after each washing step. |
| 11. | Add <b>100 µl TMB substrate</b> (SUB) into each well.  |

|     |  |
|-----|--|
| 12. | Incubate for <b>10–20 minutes*</b> at room temperature (15–30°C) <b>in the dark.</b>   |
| 13. | Add <b>100 µl</b> ELISA stop solution (STOP) and mix well.   |
| 14. | Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. |

\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

We recommend using the cubic spline algorithm. The following algorithms can be used alternatively to calculate the results.

### 1. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 2. 4-Parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

### 3. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

## 9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard can be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

It is strongly recommended that each laboratory determines its own normal and abnormal values.

Studies have been performed with the Intact Proinsulin Kit with adult males and females that had been diagnosed as having type 2 diabetes previously and were being treated with oral anti-diabetes drugs (1-2).

Samples from patients with type 2 diabetes with oral medication or dietary treatment were collected from 149 sites that participated in the IRIS-II study. In total, 2,146 male and 2,124 female patients with type 2 diabetes without insulin therapy participated in the study.

In an additional study 10 groups of 50 patients, each with incremental homeostasis model assessment (HOMA) scores, were randomly chosen out of a 4,265-person cohort in order to investigate intact proinsulin and adiponectin over a wide range of insulin resistance.

Another study evaluated 48 patients with type 2 diabetes and on oral anti-diabetic treatment. Twenty women and 28 men, aged 60 ( $\pm 9$  years), were studied by means of an intravenous glucose tolerance test. Determinations of fasting values of intact proinsulin, insulin, resistin, adiponectin, and glucose were performed. The results of these studies showed that a fasting intact proinsulin concentration of  $> 10$  pmol/l predicts the presence of insulin resistance in patients with type 2 diabetes mellitus at a very high specificity and high sensitivity. Fasting proinsulin levels in normal subjects were found to be  $< 10$  pmol/l.

These studies were conducted prior to WHO International Reference Reagent IRR 84/611 being replaced by IRR 09/296. Based on the value assignment of the new preparation samples measured against IRR 09/296 read 72.5% of the value measured against 84/611. The equivalent fasting plasma concentration in normal subjects when measured against the new reference reagent (IRR 09/296) used in this kit would be  $< 7.25$  pmol/l.

### Limitations

- The values obtained from this assay are intended to aid in diagnosis only. As with all serological tests, interpretation of results obtained with this test must be used in conjunction with the patient's clinical symptoms, medical history and other clinical and/or laboratory findings.
- Only if test instructions are rigidly followed will optimum results be achieved.
- Use fresh plasma or specimens frozen and thawed no more than twice. Specimens that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.
- Reproducible results depend on careful pipetting, observation of incubation periods and temperature, as well as thorough mixing of all prepared solutions.
- While rinsing, check that all wells are filled evenly with washing solution, and that there are no residues in the wells.
- Instructions for using appropriate microtiter plater reader are to be observed. Check that the instrument has the correct measurement protocol installed.

### Interfering substances

Interferences were studied in accordance with CLSI recommendations (CLSI EP7-A2). To study the effect of lipaemia, test pools were prepared by spiking plasma samples with a commercial lipid emulsion (Intralipid Sigma). Test samples for investigating the effect of haemolysis were obtained by osmotic shock. Icteric samples were prepared by spiking plasma samples with commercial bilirubin (Sigma).

No effect of lipaemia was observed at a lipaemic index up to 975. Interference due to haemolysis was not apparent at a haemolysis index up to 467. Bilirubin produced no apparent interference up to an icterus index of 1065.

## 11. PERFORMANCE CHARACTERISTICS

### Precision and reproducibility

A study was performed using 2 control samples assayed in 20 individual assays in duplicate. The following results were obtained:

#### Intra-Assay

| Intact Proinsulin<br>[pmol/l] | SD<br>[pmol/l] | CV<br>[%] |
|-------------------------------|----------------|-----------|
| 3.8                           | 0.19           | 4.9       |
| 70                            | 1.84           | 2.6       |

**Inter-Assay**

| <b>Intact Proinsulin<br/>[pmol/l]</b> | <b>SD<br/>[pmol/l]</b> | <b>CV<br/>[%]</b> |
|---------------------------------------|------------------------|-------------------|
| 3.8                                   | 0.33                   | 8.8               |
| 70                                    | 3.82                   | 5.5               |

*Spiking recovery*

Five plasma samples containing low endogenous intact proinsulin were spiked with recombinant proinsulin at 3 concentration levels. Recoveries are shown as percentages of the expected result for samples within the range of 9 to 22 pmol/l.

| <b>Sample</b> | <b>1</b> | <b>2</b> | <b>3</b> | <b>4</b> | <b>5</b> |
|---------------|----------|----------|----------|----------|----------|
| Spike 5 %     | 102.4    | 107.5    | 100.4    | 98.8     | 97.6     |
| Spike 10 %    | 105.1    | 107.1    | 101.9    | 101.9    | 96.1     |
| Spike 15 %    | 104.4    | 107.5    | 101.3    | 101.3    | 100.4    |

Mean spiking recovery was 102.4%.

*Linearity*

Five patient samples containing elevated proinsulin concentrations were diluted in Sample dilution buffer. The following table shows the measured intact proinsulin concentrations of the undiluted and diluted specimens in pmol/l.

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*Cross reactivity*

Cross reactivities of related proteins were investigated at concentrations of 100 pmol/l. Results are expressed as percentages of the reactivity of an identical concentration of intact proinsulin.

- Intact Proinsulin 100.0%
- Insulin 0.0%
- 32-33 split proinsulin 5.6%
- Des 31-32 split proinsulin 1.4%
- 65-66 split proinsulin 37%
- Des 65-66 split proinsulin 63%
- C-peptide 0.0%

### Sensitivity

Limit of blank, LoB

0,02 pmol/l

Limit of detection, LoD


0,36 pmol/l

### High dose hook effect

Because of the assay architecture, which employs separate incubations with solid phase and labeled antibodies, no high dose hook effect is experienced.

### Comparison with the gold standard chemoluminescence

#### Correlation Proinsulin intact ELISA / LUMI

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide or ProClin are hazardous to health and the environment. Substrates for enzymatic colour reactions may also cause skin and/or respiratory irritation. Any contact with the substances must be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.  
**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE


- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Langenfeld MR, et al. (2004) IRIS II Study: Sensitivity and specificity of intact proinsulin, adiponectin and the proinsulin/adiponectin ratio as markers for insulin resistance. *Diabetes Technology & Therapeutics*. **6**:836-843
2. Pfützner A, et al. (2005) IRIS II Study: Intact proinsulin is confirmed as a highly specific indicator for insulin resistance in a large cross-sectional study design. *Diabetes Technology & Therapeutics*. **7**:478-486

**Used symbols:**

|   |                                    |   |                                   |
|---|------------------------------------|---|-----------------------------------|
|  | Temperature limitation             |  | Catalogue number                  |
|  | In Vitro Diagnostic Medical Device |  | To be used with                   |
|  | Manufacturer                       |  | Contains sufficient for <n> tests |
|  | Lot number                         |  | Use by                            |
|  | Attention                          |  | Consult instructions for use      |
|  | Consult specification data sheet   |  | Irritant                          |