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
### 15. LITERATUR

1. Clement O, Dewachter P, Mouton-Faivre C, Nevoret C, Guilloux L, et al. Immediate Hypersensitivity to Contrast Agents: The French 5-year CIRTACI Study. *EClinical*

- Medicine*. 2018 Jul 28;1:51-61. doi: 10.1016/ j.eclinm.2018.07. 002. PMID: 31193689; PMCID: PMC6537532.
2. Kakolyri M, Strikou F, Kavallari A, Chliva C, Makris M, Tiligada E. Increased Basal Blood Histamine Levels in Patients with Self-Reported Hypersensitivity to Non-Steroidal Anti-Inflammatory Drugs. *Int Arch Allergy Immunol*. 2020;181(1):24-30. doi: 10.1159/000503968. Epub 2019 Nov 21. PMID: 31752003.
  3. Kohno R, Cannom DS, Olshansky , XiSC, Krishnappa D, et al. Mast Cell Activation Disorder and Postural Orthostatic Tachycardia Syndrome: A Clinical Association. *J Am Heart Assoc*. 2021 Sep 7;10(17):e021002. doi: 10.1161/JAHA.121.021002. Epub 2021 Aug 16. PMID: 34398691; PMCID: PMC8649306
  4. Kovacova-Hanuszkova E, Buday T, Gavliakova S, Plevkova J. Histamine, histamine intoxication and intolerance. *Allergol Immunopathol (Madr)*. 2015;43(5):498-506. doi:10.1016/j.aller.2015.05.001
  5. Li H, Xiao Y, Li Q, Yee C, Lee JJ, Yu D, et al. The allergy mediator histamine confers resistance to immunotherapy in cancer patients via activation of the macrophage histamine receptor H1. *Cancer Cell*. 2022 Jan 10;40(1):36-52.e9. doi: 10.1016/j.ccell.2021.11.002. Epub 2021 Nov 24. PMID: 34822775; PMCID: PMC8779329
  6. Nuñez-Borque E, Fernandez-Bravo S, Yuste-Montalvo A, Esteban V. Pathophysiological, Cellular, and Molecular Events of the Vascular System in Anaphylaxis. *Front Immunol*. 2022 Mar 8;13:836222. doi: 10.3389/fimmu.2022.836222. PMID: 35371072; PMCID: PMC8965328
  7. Shibata Y, Hirota S, Saito I, Asahina A. Diffuse cutaneous mastocytosis: Identification of KIT mutation and long-term follow-up with serum tryptase level. *J Dermatol*. 2021 May;48(5):672-675. doi: 10.1111/1346-8138.15764. Epub 2021 Jan 31. PMID: 33521998
  8. Zou W, Yang S, Chen L, Hu S, Hao G, Hu C. Iodixanol activation of mast cells: Implications in the pathogenesis of iodixanol-induced delayed cutaneous adverse reactions. *Toxicology*. 2022 Jan 15;465:153034. doi: 10.1016/j.tox.2021.153034. Epub 2021 Nov 11. PMID: 34774977

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
**Manual***For professional use only*

# Histamine ELISA

*For the in vitro determination of histamine  
in EDTA plasma and urine*

Valid from 2025-01-21

**REF** K 8212**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: + 49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)[www.immundiagnostik.com](http://www.immundiagnostik.com)

## **Safety information**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit.

Important safety information for this product can be found in the chapter **WARNINGS AND PRECAUTIONS**.

# Table of Contents

<b>1. INTENDED USE</b>	<b>21</b>
<b>2. INTRODUCTION</b>	<b>21</b>
<b>3. MATERIAL SUPPLIED</b>	<b>21</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>22</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>22</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>23</b>
<b>7. ASSAY PROCEDURE</b>	<b>24</b>
<i>Principle of the test</i>	24
<i>Derivatisation procedure</i>	24
<i>Test procedure</i>	25
<b>8. RESULTS</b>	<b>26</b>
<b>9. LIMITATIONS</b>	<b>27</b>
<i>Biotin interference</i>	27
<b>10. QUALITY CONTROL</b>	<b>27</b>
<i>Reference range</i>	28
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>28</b>
<i>Precision and reproducibility</i>	28
<i>Spiking recovery</i>	29
<i>Dilution recovery</i>	29
<i>Analytical sensitivity</i>	30
<i>Specificity</i>	30
<i>Comparison of ELISAs</i>	31
<b>12. WARNINGS AND PRECAUTIONS</b>	<b>31</b>
<b>13. TECHNICAL HINTS</b>	<b>32</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>32</b>
<b>15. REFERENCES</b>	<b>33</b>
<b>16. SYMBOLS</b>	<b>34</b>

## 1. INTENDED USE

This Immundiagnostik AG assay is intended for the quantitative determination of histamine in EDTA plasma and urine. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Histamine is a biogenic amine that derives from the decarboxylation of histidine. It is synthesised in mast cells, basophils, platelets, histaminergic neurons and enterochromaffine cells, where it is stored in vesicles. After stimulation and release, Histamine acts by binding to its 4 receptors (H1R, H2R, H3R and H4R) on target cells in various tissues.

It causes smooth muscle cell contraction, vasodilation, increased vascular permeability, increased mucous secretion, tachycardia, alterations of blood pressure and arrhythmias, among other effects.

In humans, histamine is one of the most important mediators and takes part in the initial phase of an anaphylactic reaction ("immediate type" allergy).

The quantification of histamine release from basophilic leucocytes in allergies is also of clinical interest.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 8212	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 8212	STD	Standards, ready-to-use (0, 0.3, 0.6, 1.2, 5, 25 ng/ml)	6 x 2 ml
K 8212	CTRL 1	Control, ready-to-use (see specification for range)	1 x 2 ml
K 8212	CTRL 2	Control, ready-to-use (see specification for range)	1 x 2 ml
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 8212	AB	Histamine antibody, peroxidase-labelled, ready-to-use	1 x 6 ml
K 8212	REABUF	Reaction buffer, ready-to-use	2 x 70 ml
K 8212	DER	Derivatisation reagent, lyophilised	1 vial
K 0008.10	DMSO	Dimethylsulfoxide (DMSO)	1 x 10 ml

K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipets and 10-1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once ensure that reagents are stored at the conditions stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** can be used until the expiry date stated on the label when stored at **2-8 °C**. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.
- **DMSO** crystallises at 2-8 °C. Before use, dissolve the crystals at room temperature or in a water bath.
- The **lyophilised derivatisation reagent (DER)** can be used until the expiry date stated on the label when stored at **2-8 °C**. Bring to room temperature before opening and dissolve the content of the vial in **DMSO** as stated on the label. Allow to dissolve for **15 min** and mix thoroughly with a vortex-mixer. Ensure that it is completely dissolved. **The derivatisation reagent** (reconstituted DER) can

be stored at **2-8 °C for 2 months**. Bring to room temperature before reuse. Please note: DMSO attacks plastics but not polypropylene products and laboratory glass.

- All other test reagents are ready-to-use. Test reagents can be used until the expiry date stated on the label when stored at **2-8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### EDTA plasma

- Collect blood into a chilled tube containing EDTA only, mix gently by inversion and chill immediately on ice. Centrifuge 10 minutes at 900 *g* at 4 °C within 20 minutes of sample collection.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- EDTA plasma samples are used **undiluted**.

### Urine

- It is possible to use spontaneous as well as 24-h urine. 24-h urine: The total volume of urine excreted during 24 h is collected in a bottle containing 10-15 ml of 6 M HCl as preservative. Avoid exposure to direct sunlight. Determine total volume for calculation of results.
- Urine samples are **diluted 1:51** for derivatisation (see derivatisation procedure).
- Samples with visible amounts of precipitates should be centrifuged.

For sample preparation, a derivatisation reagent is added (see derivatisation procedure).

### Sample Storage

Freshly collected EDTA plasma can be stored for up to 6 hours at 2-8 °C. For longer storage keep plasma samples frozen at -20 °C. Avoid repeated thawing and freezing.

Spontaneous urine is stable for up to 6 hours at 2-8 °C. Acidified urine samples can be stored at 2-8 °C up to 72 h. For longer periods, store urine samples frozen at -20 °C. Avoid repeated thawing and freezing.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of histamine in EDTA plasma or urine. This assay is based on the method of competitive enzyme linked immunoassays.

The sample preparation includes the addition of a derivatisation reagent for histamine derivatisation. Afterwards, the treated samples and a peroxidase-conjugated polyclonal histamine antibody are incubated in wells of a microtiter plate coated with histamine derivative (tracer). During the incubation period, the target histamine in the sample competes with the tracer, immobilised on the wall of the microtiter wells, for the binding of the polyclonal antibodies.

After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The colour changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow colour is inverse proportional to the histamine concentration in the sample; this means, high histamine concentration in the sample reduces the concentration of tracer-bound antibody and lowers the photometric signal. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standards. Histamine, present in the patient samples, is determined directly from this curve.

### *Derivatisation procedure*

Bring **all reagents and samples to room temperature** (15-30 °C) and mix well.

For **EDTA plasma** samples **no dilution** is required.

Dilute **urine samples** with reaction buffer by **factor 1:51**, as follows:

**20 µl** urine sample + **1000 µl** reaction buffer (REABUF).

Derivatisation of standards, controls, plasma samples and diluted urine samples is carried out in vials (e.g. 1.5 ml polypropylene vials).

We recommend preparing one derivatisation per standard, control and sample and transferring it in duplicate determinations into the wells of the microtiter plate.

1.	Add <b>75 µl standard</b> (STD)/ <b>control</b> (CTRL)/ <b>plasma sample</b> or <b>diluted urine sample</b> in the corresponding vials.
2.	Add <b>250 µl reaction buffer</b> (REABUF) into each vial (STD, CTRL, sample).
3.	Add <b>75 µl derivatisation reagent</b> into each vial (STD, CTRL, sample), and <b>mix thoroughly</b> by repeated inversion or several seconds on a vortex mixer.
4.	Incubate for <b>1 hour</b> at room temperature (15-30 °C) on a horizontal <b>shaker</b> .

2 x 50 µl of the derivatised standards/controls/samples are used in the ELISA as duplicates.

### *Test procedure*

Mark the positions of standards/controls/samples in duplicate on a protocol sheet. Take as many microtiter strips (PLATE) as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2-8 °C. Strips are stable until the expiry date stated on the label.

5.	For the analysis in duplicate take <b>2 x 50 µl</b> of the <b>derivatised standards/controls/ samples</b> out of the vials and add into the respective wells of the microtiter plate.
6.	Add <b>50 µl histamine antibody</b> (AB) into each well of the microtiter plate.
7.	Cover the plate tightly with foil and incubate for <b>1 hour</b> at room temperature (15-30 °C) on a <b>horizontal shaker</b> .
8.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
9.	Add <b>100 µl substrate</b> (SUB) into each well.
10.	Incubate for <b>12-18 min*</b> at room temperature (15-30 °C) in the <b>dark</b> .
11.	Add <b>100 µl stop solution</b> (STOP) into each well and mix well.

- |     |  |
|-----|--|
| 12. | Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm (690 nm) as a reference. |
|-----|--|

\* The intensity of the colour change is temperature sensitive. We recommend to observe the colour change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. Attention should be paid to the following points:

Derivatisation of the samples is possible in 60 minutes on a horizontal shaker or, after thorough mixing, in 90 minutes without shaking in an automated processor. In automated processing, the volumes of samples and diluents may be scaled while maintaining the respective dilutions. Device-specific minimum and maximum volumes must be taken into account. The homogeneity of the resulting dilution must be ensured.

For further details please contact your supplier or Immundiagnostik AG.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the program used, the duplicate values should be evaluated manually.

### EDTA plasma

**No factor** is required for the calculation of results from plasma samples.

### Urine

The obtained histamine levels of urine samples have to be multiplied by the **dilution factor of 51** to get the actual concentrations (ng/ml = µg/l).

Calculation of 24 h excretion:  $\mu\text{g}/24\text{h} = \mu\text{g}/\text{l} \times \text{l}/24\text{h}$

Conversion:  $\text{histamine} [\text{ng}/\text{ml}] \times 8.997 = \text{histamine} [\text{nmol}/\text{l}]$

## 9. LIMITATIONS

Samples with concentrations above the measurement range can be diluted with reaction buffer and re-assayed. Please consider this dilution factor when calculating the results.

Samples with concentrations lower than the measurement range cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

*highest concentration of the standard curve × sample dilution factor to be used*

The lower limit of the measurement range can be calculated as:

*analytical sensitivity × sample dilution factor to be used*

Analytical sensitivity see chapter "Performance Characteristics".

### *Biotin interference*

Samples containing a biotin concentration of  $\leq 44$  ng/ml show a change of the results of  $< 25$  %. Higher concentrations of biotin can lead to false results. Patients taking  $> 5$  mg biotin per day should wait at least 24 hours after taking biotin to have their samples collected. Results of patients taking biotin supplements or receiving a high-dose biotin therapy should generally be interpreted along with the total clinical picture.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within

the same assay one or more values of the quality control samples are outside of the acceptable limits.

### *Reference range*

Based on internal studies with plasma samples of apparently healthy persons (n = 36) a mean value of 0.67 ng/ml was calculated. The standard deviation (SD) was 0.18 ng/ml. The following normal range was calculated from mean value + 2 SD:

**EDTA plasma: < 1 ng/ml**

For urine, a mean value of 18.5 µg/g creatinine was calculated based on internal studies with samples of apparently healthy persons (n = 35). The standard deviation (SD) was 10.6. The following normal range was calculated from mean value + 2 SD:

**Urine: < 40 µg histamine / g creatinine**

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

#### **EDTA plasma**

Intra assay (n = 6)

sample	histamine [ng/ml]	CV [%]
1	1.36	8.9
2	3.06	6.2

Inter assay (n = 6)

sample	histamine [ng/ml]	CV [%]
1	0.84	10.0
2	0.54	10.1

#### **Urine**

Intra-Assay (n = 6)

sample	histamine [ng/ml]	CV [%]
1	10.10	10.5
2	16.02	8.7

Inter-Assay (n = 6)

sample	histamine [ng/ml]	CV [%]
1	8.49	12.1
2	12.80	8.1

### *Spiking recovery*

Two samples were spiked with different histamine concentrations and measured in this assay. The mean recovery rate was 105.4 % for plasma and 98.9% for urine (n = 2).

#### **EDTA plasma**

sample [ng/ml]	spike [ng/ml]	expected [ng/ml]	measured [ng/ml]	recovery [%]
0.96	3	3.96	4.46	112.6
	5	5.96	6.19	103.9
1.01	3	4.01	3.90	97.3
	5	6.01	6.46	107.6

#### **Urine**

sample [ng/ml]	spike [ng/ml]	expected [ng/ml]	measured [ng/ml]	recovery [%]
10.10	15	25.10	27.51	109.6
	30	40.10	39.88	99.5
16.03	15	31.03	28.02	90.3
	30	46.03	44.28	96.2

### *Dilution recovery*

Two spiked samples were diluted with reaction buffer. The mean recovery rate was 102.3 % for plasma and 102.4 % for urine (n = 2).

**EDTA plasma**

sample [ng/ml]	dilution	expected [ng/ml]	measured [ng/ml]	recovery [%]
6.19	1:2	3.09	3.16	102.3
	1:4	1.55	1.70	109.7
6.46	1:2	3.23	2.96	91.6
	1:4	1.62	1.83	113.0

**Urine**

sample [ng/ml]	dilution	expected [ng/ml]	measured [ng/ml]	recovery [%]
39.88	1:2	19.94	18.64	93.5
	1:4	9.97	9.99	100.2
44.28	1:2	22.14	21.94	99.1
	1:4	11.07	11.70	105.7

*Analytical sensitivity*

The zero standard was measured 75 times. The detection limit was set as  $B_0 - 2 SD$  and estimated to be 0.11 ng/ml.


*Specificity*

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to histamine. The specificity is calculated in percent in relation to the histamine binding activity:

3-methylhistamine	<0.1 %
tyramine	<0.001 %
L-phenylalanine	<0.0002 %
L-histidine	<0.0002 %
L-tyrosine	<0.0002 %
tryptamine	<0.0002 %
5-hydroxyindoleacetic acid	<0.0002 %
serotonin (5-hydroxytryptamine)	<0.0002 %

## Comparison of ELISAs

24 EDTA plasma samples and 16 urine samples were measured with this ELISA and with a commercially available ELISA. The correlation was  $r = 0.99$  for plasma and for urine.

## 12. WARNINGS AND PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Kit reagents contain ProClin or thimerosal as bactericides. ProClin and thimerosal are harmful to health and the environment. Substrates for enzymatic colour reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.

- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

### 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE


- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immunodiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Clement O, Dewachter P, Mouton-Faivre C, Nevoret C, Guilloux L, et al. Immediate Hypersensitivity to Contrast Agents: The French 5-year CIRTACI Study. *EClinical Medicine*. 2018 Jul 28;1:51-61. doi: 10.1016/j.eclinm.2018.07.002. PMID: 31193689; PMCID: PMC6537532.
2. Kakolyri M, Strikou F, Kavallari A, Chliva C, Makris M, Tiligada E. Increased Basal Blood Histamine Levels in Patients with Self-Reported Hypersensitivity to Non-Steroidal Anti-Inflammatory Drugs. *Int Arch Allergy Immunol*. 2020;181(1):24-30. doi: 10.1159/000503968. Epub 2019 Nov 21. PMID: 31752003.
3. Kohno R, Cannom DS, Olshansky , XiSC, Krishnappa D, et al. Mast Cell Activation Disorder and Postural Orthostatic Tachycardia Syndrome: A Clinical Association. *J Am Heart Assoc*. 2021 Sep 7;10(17):e021002. doi: 10.1161/JAHA.121.021002. Epub 2021 Aug 16. PMID: 34398691; PMCID: PMC8649306
4. Kovacova-Hanuszkova E, Buday T, Gavliakova S, Plevkova J. Histamine, histamine intoxication and intolerance. *Allergol Immunopathol (Madr)*. 2015;43(5):498-506. doi:10.1016/j.aller.2015.05.001
5. Li H, Xiao Y, Li Q, Yee C, Lee JJ, Yu D, et al. The allergy mediator histamine confers resistance to immunotherapy in cancer patients via activation of the macrophage histamine receptor H1. *Cancer Cell*. 2022 Jan 10;40(1):36-52.e9. doi: 10.1016/j.ccell.2021.11.002. Epub 2021 Nov 24. PMID: 34822775; PMCID: PMC8779329
6. Nuñez-Borque E, Fernandez-Bravo S, Yuste-Montalvo A, Esteban V. Pathophysiological, Cellular, and Molecular Events of the Vascular System in Anaphylaxis. *Front Immunol*. 2022 Mar 8;13:836222. doi: 10.3389/fimmu.2022.836222. PMID: 35371072; PMCID: PMC8965328
7. Shibata Y, Hirota S, Saito I, Asahina A. Diffuse cutaneous mastocytosis: Identification of KIT mutation and long-term follow-up with serum tryptase level. *J Dermatol*. 2021 May;48(5):672-675. doi: 10.1111/1346-8138.15764. Epub 2021 Jan 31. PMID: 33521998
8. Zou W, Yang S, Chen L, Hu S, Hao G, Hu C. Iodixanol activation of mast cells: Implications in the pathogenesis of iodixanol-induced delayed cutaneous adverse reactions. *Toxicology*. 2022 Jan 15;465:153034. doi: 10.1016/j.tox.2021.153034. Epub 2021 Nov 11. PMID: 34774977

## 16. SYMBOLS

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin