
# IDKmonitor<sup>®</sup> Infliximab free ADA ELISA

***Zur in-vitro-Bestimmung von freien humanen Antikörpern gegen Infliximab (z. B. REMICADE<sup>®</sup>) in EDTA-Plasma und Serum***  
***For the in vitro determination of free human antibodies against infliximab (e. g. REMICADE<sup>®</sup>) in EDTA plasma and serum***

Gültig ab / Valid from 2024-02-26

**REF** K 9650

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Inhalt

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
# **IDKmonitor<sup>®</sup> Infliximab free ADA ELISA**

***For the in vitro determination of free human antibodies  
against infliximab (e.g. REMICADE<sup>®</sup>)  
in EDTA plasma and serum***

Valid from 2024-02-26

**REF** **K 9650**

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>6. SAMPLE PREPARATION AND STORAGE</b>	<b>17</b>
<i>Storage of samples</i>	17
<i>Sample preparation</i>	17
<b>7. ASSAY PROCEDURE</b>	<b>18</b>
<i>Principle of the test</i>	18
<i>Test procedure</i>	18
<b>8. RESULTS</b>	<b>19</b>
<b>9. LIMITATIONS</b>	<b>20</b>
<b>10. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Accuracy – Precision</i>	20
<i>Analytical sensitivity</i>	21
<i>Linearity</i>	21
<b>12. PRECAUTIONS</b>	<b>22</b>
<b>13. TECHNICAL HINTS</b>	<b>22</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>15. REFERENCES</b>	<b>23</b>

## 1. INTENDED USE

This Immundiagnostik AG assay is intended for the qualitative determination of free human antibodies against Infliximab (e.g. REMICADE®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis or rheumatoid arthritis are often treated with anti-TNF $\alpha$  antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF $\alpha$  therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF $\alpha$  antibody [1-3]. Several factors influence the trough level, among them dosage and frequency of anti-TNF $\alpha$  blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation anti-drug antibodies, ADA). It is thought that ADA functionally neutralise the therapeutic antibodies or induce their rapid elimination [4]. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF $\alpha$  antibody application [1, 5].

The IDKmonitor® Infliximab free ADA ELISA for the detection of antibodies against Infliximab (e.g. REMICADE®) measures free antibodies against Infliximab. A co-determination of rheuma factors or irregular antibodies can be excluded. In combination with the drug level determination the IDKmonitor® Infliximab free ADA ELISA is an opportunity for the treating physician to monitor and optimise the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9650	PLATE	Microtiter plate, pre-coated with (F(ab) <sub>2</sub> )	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate 10x	1 x 100 ml
K 9650	CONJ	Conjugate concentrate, (therapy antibody, peroxidase labelled)	1 x 200 $\mu$ l
K 9650	CTRL POS	Positive control, lyophilised (see specification for range)	4 x 1 vial
K 9650	CTRL NEG	Negative control, lyophilised (see specification for range)	4 x 1 vial

Cat. No.	Label	Kit components	Quantity
K 9650	CTRL CUT-OFF	Cut-off control, lyophilised	4 x 1 vial
K 9650	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 30 ml
K 9650	ASYBUF	Assay buffer, ready-to-use	1 x 10 ml
K 0002.15	SUB	Substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Absorbent paper
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redis-

solved at room temperature or in a water bath at 37°C. The **WASHBUF** is stable at **2–8°C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8°C for 1 month**.

- The **lyophilised controls** (CTRL NEG, CTRL POS, CTRL CUT-OFF) are stable at **2–8°C** until the expiry date stated on the label. Before use, the controls have to be reconstituted with **300 µl ultrapure water** and mixed by gentle inversion to ensure complete reconstitution. Allow the vial content to dissolve for 10 minutes and then mix thoroughly. **Reconstituted controls are not stable and cannot be stored.**
- **Preparation of the conjugate:** The **conjugate concentrate (CONJ)** has to be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The CONJ is stable at **2–8°C** until expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8°C**.

## 6. SAMPLE PREPARATION AND STORAGE

### *Storage of samples*

Undiluted samples can be stored 2 months at -20 °C and for 7 days at 2–8°C or room temperature. Avoid more than 3 freeze-thaw-cycles.

Diluted samples are not stable and cannot be stored.

### *Sample preparation*

1.	Transfer <b>50 µl of each sample</b> in a 1.5 ml reaction tube and add <b>250 µl of sample dilution buffer</b> (SAMPLEBUF). Vortex well.
2.	Incubate for <b>15 min</b> at room temperature (15–30 °C) with gentle <b>shaking</b> .
3.	Add <b>50 µl of assay buffer</b> (ASYBUF) to each sample. Vortex well.
4.	Incubate for <b>15 min</b> at room temperature (15–30 °C) with gentle <b>shaking</b> .

For analysis, pipet 100 µl of each prepared sample per well. For the recommended analysis in duplicate, 2 x 100 µl are required.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the determination of free antibodies against infliximab (e.g. REMICADE®). In a first incubation step, the free anti-infliximab antibodies from the sample are bound to the infliximab F(ab)<sub>2</sub> fragments coated on the plate. To remove all unbound substances, a washing step is carried out. In a further incubation step, peroxidase labelled therapy antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic stop solution is then added. The colour converts to yellow. The absorbance of the colour compound is determined photometrically at 450 nm. The intensity of the colour is directly proportional to the amount of bound anti-infliximab antibodies (e.g. REMICADE®) from the sample. The results are evaluated by a cut-off control.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	<b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Add each <b>100 µl controls/prepared samples</b> into the respective wells.
3.	Seal the strips with foil and incubate <b>over night</b> (16–20 h) at <b>2–8 °C</b> on a <b>horizontal shaker</b> .*
4.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5.	Add <b>100 µl conjugate</b> (diluted CONJ) into each well.

6.	Seal the strips with foil and incubate for <b>1 hour</b> at room temperature (15–30°C) on a <b>horizontal shaker**</b> .
7.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8.	Add <b>100 µl substrate</b> (SUB) into each well.
9.	Incubate for <b>10–20 min***</b> at room temperature (15–30°C) in the <b>dark</b> .
10.	Add <b>100 µl stop solution</b> (STOP) into each well and mix well.
11.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference.

\* The above incubation step at 2–8°C and 550 rpm with an orbit of 2 mm is recommended by the producer. If there is no possibility to incubate at 2–8°C, while shaking, we recommend to incubate at 2–8°C without any shaking.

\*\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The analysis of the results is done using the cut-off control. Samples with a higher optical density (OD) as the OD of the cut-off control are positive. Samples with an OD lower than the OD of the cut-off control are negative.

$$\text{Cut-off} = 10 \text{ AU/ml} = \text{OD}_{\text{cut-off control}}$$

For the calculation of the sample concentrations, linear regression using a linear ordinate and abscissa is recommended.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### Sample calculation for a positive sample

average OD of patient's sample	0.735
average OD of cut-off control	0.065 = 10 AU/ml
Concentration of the patient's sample	$\frac{0,735 \times 10 \text{ AU/ml}}{0,065} = 113 \text{ AU/ml}$

## 9. LIMITATIONS

The lower limit of the measurement range is the LoB.

LoB see chapter "Performance Characteristics".

Samples with concentrations lower than the measurement range cannot be clearly quantified.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n=40**

The repeatability was assessed with 3 serum-samples under constant parameters (same operator, measurement system, day and kit lot).

Sample	Mean value [AU/ml]	CV [%]
1	69.45	3.0
2	167.58	3.3
3	118.86	3.2

### Reproducibility (Inter-Assay); n=11

The reproducibility was assessed with 5 serum-samples under varying parameters (different operators, measurement systems, days and kit lots).

Sample	Mean value [AU/ml]	CV [%]
1	174.56	10.2
2	107.64	12.5
3	47.28	10.5
4	67.24	10.4
5	73.07	11.4

### Analytical sensitivity

The LoB (limit of blank) was evaluated according to the CLSI guideline EP17-A2 and resulted in 5.751 AU/ml

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP06-A by serial dilution of 2 serum samples.

For anti-Infliximab antibodies in serum and plasma, the method has been demonstrated to be linear from 19.83–353.83 AU/ml, showing a non-linear behaviour of less than  $\pm 20\%$  in this interval.

Lot	Dilution	Expected [AU/ml]	Obtained [AU/ml]	Recovery [%]
1	1:20	353.83	353.83	100.00
	1:40	176.92	213.08	120.44
	1:80	88.46	98.50	111.35
	1:160	44.23	40.17	90.81
	1:320	22.11	19.83	89.68
2	1:40	324.00	324.00	100.00
	1:80	162.00	162.33	100.21
	1:160	81.00	68.83	84.98
	1:320	40.50	30.83	76.13

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide or ProClin are hazardous to health and the environment. Substrates for enzymatic colour reactions may also cause skin and/or respiratory irritation. Any contact with the substances must be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE




- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDKmonitor® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Afif, W. et al., 2010. Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *The American journal of gastroenterology*, **105**(5), pp.1133–9.
2. Kopylov, U. et al., 2012. Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-infliximab antibodies. *Inflammatory bowel diseases*, **18**(9), pp.1628–33.
3. Tak, P.P., 2012. A personalized medicine approach to biologic treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology (Oxford, England)*, **51**(4), pp.600–9.
4. Ordás, I. et al., 2012. Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clinical pharmacology and therapeutics*, **91**(4), pp.635–46.
5. Bender, N.K. et al., 2007. Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatology international*, **27**(3), pp.269–74.

**Used symbols:**

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant