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
# Vitamin A/E HPLC Kit

*For the determination of vitamin A/E in plasma and serum*

Valid from 2024-05-13

**REF** **KC1650**

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

**Safety information**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit. Important safety information for this product can be found in the chapter 16.

## Table of Contents

1. INTENDED PURPOSE	11
2. INTRODUCTION	11
3. PRINCIPLE OF THE TEST	11
4. MATERIAL SUPPLIED	12
5. MATERIAL REQUIRED BUT NOT SUPPLIED	12
6. PREPARATION AND STORAGE OF REAGENTS	12
7. SPECIMEN COLLECTION AND PREPARATION	13
8. PRECAUTIONS	13
9. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	13
10. REFERENCES	14
11. SYMBOLS	15

## 1. INTENDED PURPOSE

This HPLC application is intended for the quantitative determination of vitamin A/E in serum and plasma. This assay is designed for *in vitro* diagnostic use only.

## 2. INTRODUCTION

Vitamin A and E belong to the fat-soluble vitamins and can be stored in the fatty tissue over a longer period of time. Both an over- and an undersupply of vitamins A and E can manifest themselves in complaints.

Vitamin A (retinol) is essential for the visual process and keeps the skin and mucous membranes healthy. A deficiency of vitamin A can therefore have an impact on vision, especially during transitions from light to dark. A severe vitamin A deficiency can lead to blindness. Too high doses however, can also have a pathological effect and cause headaches, skin changes, liver damage, painful skeletal changes and damage to the foetus.

Vitamin E (tocopherol) is a natural antioxidant that protects unsaturated fatty acids from oxidation. It intercepts radicals before they can have a destructive effect on the cell. A deficiency of vitamin E showed damage to the musculature, nervous system and heart, liver and reproduction in animal experiments. In humans such effects have not been observed. Vitamin E can be stored in large quantities in fatty tissue. An undersupply can occur due to impaired fat digestion or absorption.

## 3. PRINCIPLE OF THE TEST

For the determination of vitamins A and E, a very simple sample preparation with a precipitation step is carried out

The separation via HPLC follows an isocratic method at 30 °C using a “reversed phase” column; one run lasts approx. 3 minutes. The detection is performed by an UV detector at two different wavelengths (Vitamin A: 325 nm; Vitamin E: 290 nm). The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas/heights in the external standard calibration mode.

### Summary

The application of vitamin A and E for HPLC allows the determination of both vitamins in an easy, fast, and precise method. The kit includes reagents in ready-to-use form for preparation and separation of the samples with exception of the columns.

#### 4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KC1650	MOPHA	Mobile phase, ready-to-use (can be recirculated)	1 000 ml
	PREC	Precipitation reagent, ready-to-use	50 ml
	CAL	Calibrator, lyophilised	4 x
	CTRL1	Control 1, lyophilised	4 x
	CTRL2	Control 2, lyophilised	4 x
K 0005.15	RECSOL	Reconstitution solution	15 ml

The HPLC column (KC1650RP), can be ordered separately from Immundiagnostik AG. To extend the lifetime of your HPLC column, pre-columns (KC1650VS) are highly recommended. These and also the pre-column holders (KC1650VH) can also be ordered from Immundiagnostik AG.

In addition to the complete kits, all components can also be ordered separately. Please ask for our single component price list.

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes (100 µl, 1 000 µl)
- Vortex
- HPLC system with UV-detector
- Reversed phase C18-column

#### 6. PREPARATION AND STORAGE OF REAGENTS

- The lyophilised controls (**CTRL1 and CTRL2**) are stable at -20°C until the expiry date stated on the label. Before use, they have to be reconstituted each with x µl reconstitution solution (**RECSOL**) (x= please see product specification). The concentration of vitamin A/E slightly changes from lot to lot, the exact concentration is stated in the product specification.

- The lyophilised calibrator (**CAL**) is stable at -20°C until the expiry date stated on the label. Before use, it has to be reconstituted with xµl reconstitution solution (**RECSOL**) (x= please see product specification). The concentration of vitamin A/E slightly changes from lot to lot, the exact concentration is stated in the product specification.
- All other test reagents are stable until the expiry date (see label of test package) when stored at 2–8°C.

## 7. SPECIMEN COLLECTION AND PREPARATION

Serum and EDTA plasma are suited as samples. The sample is light and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately. The samples are stable in the dark at 2–8°C for 12 h (vitamin A) and 3 days (vitamin E). Vitamin A is stable for 1 month, vitamin E for at least 3 months when stored at -20°C.

m should be equilibrated with ~ 30 ml mobile phase (**MOPHA**).

## 8. PRECAUTIONS

- Calibrator (**CAL**) and controls (**CTRL1** and **CTRL2**) are based on human plasma. They have been tested for HIV and hepatitis B and found to be negative. However, as a precaution, the test components should always be handled as potentially infectious material.
- The test components contain organic solvents. Avoid contact with skin or mucous membranes.

## 9. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Do not use reagents beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Do not mix plugs and caps from different reagents.
- Follow the guidelines for medical laboratories.



- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be lodged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Analyse controls with each run.
- Always perform the assay according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

## 10. REFERENCES

1. Sushil K.J., Mc Coy B., Wise R. (1994). Vitamin E and the hypercoagulability of neonatal blood. *Clin Cim Acta* **225**; 97-103.
2. Comstock G.W., Alberg A.J., Helzlsouer K.J. (1993). Reported effects of long-term freezer storage on concentrations of retinol,  $\beta$ -carotene, and  $\alpha$ -tocopherol in serum or plasma summarised. *Clin Chem* **39/6**; 1075-1078

## 11. SYMBOLS

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)