**Kurzanleitung / Short manual**

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
*For professional use only*

# **IDK<sup>®</sup> Vitamin K<sub>1</sub>/K<sub>2</sub> (SPE) HPLC Kit**

***For the in vitro determination of vitamin K<sub>1</sub> and K<sub>2</sub>  
in plasma and serum***

Valid from 2024-04-02

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany  
Tel.: +49 6251 70190-0                      Fax: + 49 6251 70190-363  
e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)    [www.immundiagnostik.com](http://www.immundiagnostik.com)

## **Safety information**

These accessories are to be used exclusively in accordance with the enclosed instructions for use. Important safety information for this product can be found in chapter 7.

# Table of Contents

1. INTENDED PURPOSE	12
2. INTRODUCTION	12
3. MATERIAL SUPPLIED	13
4. MATERIAL REQUIRED BUT NOT SUPPLIED	14
5. STORAGE AND PREPARATION OF REAGENTS	14
6. STORAGE AND PREPARATION OF SAMPLES	16
7. PRECAUTIONS	16
8. GENERAL NOTES ON THE TEST	16
9. REFERENCES	17
10. SYMBOL EXPLANATION	18

## 1. INTENDED PURPOSE

The Immundiagnostik AG assay is intended for the quantitative determination of vitamin K<sub>1</sub> and K<sub>2</sub> in plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Vitamin K is metabolized rapidly and is stored in the organism only in small quantities, so that deficiency symptoms become apparent after only a few days. The clinical symptoms are mainly coagulation disorders and are manifested by bleeding in the skin, mucous membranes, muscle tissue and internal organs. In addition, bone metabolism disorders due to deficient modification of osteocalcin have also been described.

The deficiency symptoms are consequences of a reduced effect of the vitamin as a cofactor of a microsomal, O<sub>2</sub>-dependent carboxylase in the carboxylation of specific glutamic acid residues of the coagulation factors, as well as of osteocalcin. The resulting reduced formation of  $\gamma$ -carboxyglutamic acid residues, which give the coagulation factors the ability to form complexes with Ca<sup>2+</sup> and membrane phospholipids, from which the active coagulation factors are then released after proteolytic cleavage, ultimately leads to the observed symptoms. Osteocalcin loses the ability to bind to the calcium hydroxyapatite of the bone matrix due to deficient  $\gamma$ -carboxylation. In addition, severe anemia may develop as a result of hypoplastic bone marrow changes in cases of severe vitamin K deficiency.

### Indications

- Blood clotting disorders
- Bone metabolism disorders
- Vitamin K deficiency induced by:
  - obstructive liver disease
  - obstructive icterus
  - malabsorption due to celiac disease, pancreatitis, diarrhea
  - long-term antibiotic therapy: destruction of the intestinal flora
- Haemorrhagic disorders of newborns

### 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KC2400	MOPHA	Mobile phase, ready-to-use (important: do not recirculate)	3 x 1 000 ml
	INTSTD	Internal Standard, ready-to-use	1 ml
	STD	Isopropanolic standard, ready-to-use	2 ml
	ZINC	Zinc	12 g
	ZUB	Accessories for post-column reduction reactor	1 x
KC2412	CAL	Calibrator, lyophilised	4 x
	CTRL1	Control 1, lyophilised	4 x
	CTRL2	Control 2, lyophilised	4 x
	ELUSOL	Elution solution, ready-to-use	110 ml
	PREC	Precipitating reagent, ready-to-use	525 ml
	WASHSOL	Wash solution, ready-to-use	315 ml
KC0002.50	SPE	Solid phase extraction cartridges	2 x 50
K 0005.15	RECSOL	Reconstitution solution, ready-to-use	1 x 15 ml

The HPLC column (KC2400RP) and the post-column reduction reactor (KC2400NR) can be ordered separately from Immundiagnostik AG. To extend the lifetime of your HPLC column, pre-columns (KC2400VS) should ideally be used. These and also the corresponding pre-column holders (KC2400VH) can also be ordered from Immundiagnostik AG. In addition to the complete kit, all components as well as required materials (e.g. borosilicate glass tubes (KC0001V), evaporator (KC0001SV)) can be ordered separately. Please ask for our single component price list.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Centrifuge
- Vortex
- 10 ml borosilicate glass tube, V-bottom
- 15 ml centrifuge tubes (e.g. Falcon tubes)
- HPLC-system with fluorescence detector
- Post-column reduction reactor
- Reversed phase C18-column
- Various precision pipettes
- Device for concentrating the sample (e.g. vacuum centrifuge or evaporator)

#### 5. STORAGE AND PREPARATION OF REAGENTS

##### Preparation of the calibrator and controls

All reagents are ready-to-use in dissolved form, except for the calibrator (**CAL**) and controls (**CTRL1** and **CTRL2**). Reconstitute calibrator (**CAL**) in  $x\ \mu\text{l}$  ( $x$  = see the enclosed product specification sheet) reconstitution solution (**RECSOL**). Unused reconstituted material must be disposed of, as refreezing partially destroys the calibrator. The content of vitamin K<sub>1</sub> and K<sub>2</sub> varies slightly from batch to batch, please refer to the enclosed product specification sheet for the exact amount.

Reconstitute controls (**CTRL1** and **CTRL2**) in  $x\ \mu\text{l}$  ( $x$  = see the enclosed product specification sheet) reconstitution solution (**RECSOL**).

The calibrator (**CAL**), internal standard (**INTSTD**), controls (**CTRL1** and **CTRL2**) and standard (**STD**) can be used up to the date of expiry stated on the label when stored at -20 °C. The other components, with the exception of zinc and the filters/seals for the reactor, must be stored at 2-8 °C.

##### Preparation of the mobile phase

If the used HPLC-system does not have an in-line degasser, the mobile phase is degassed in an ultrasonic bath for 10–15 min. Since the mobile phase is an organic solvent, care must be taken to seal the autosampler vials for injection into the HPLC-system with a septum in order to avoid evaporation and to ensure the consistent quality of the measurement.


### Preparation of the post-column reduction reactor

Vitamin K<sub>1</sub> and K<sub>2</sub> are reduced in the post-column reduction reactor and thus made accessible for detection in the fluorescence detector. In order to obtain a constant detector signal quality, it is necessary to freshly fill the post-column reduction reactor before each series of samples, as the surface of the zinc particles is oxidised after 48 h of running time. The filling process is very simple and takes only about 5 min.

For assembling the column see the following figure.

1. Cap nut
2. Stainless steel inlet
3. PTFE seal
4. Stainless steel sieve (grey)
5. Glass fiber sieve (3 pieces, white)
6. Stainless steel sieve (grey)
7. Column tube

1. Close one side of the column according to the figure above.
2. Fill in the zinc particles with a funnel or spatula while knocking the column slightly on the table, so that the packing will not show any cavities.
3. Close the upper side of the column.

The post-column reduction reactor should be mounted in the HPLC-system as described by the following picture:

**Note:** The performance of the reactor must be checked (see chapter 7. "Assay procedure" - Test procedure).

## 6. STORAGE AND PREPARATION OF SAMPLES

Plasma and serum can both be used for analysis. The samples are stable for 7 days at room temperature in the absence of light and can also be centrifuged 72 hours later in the absence of light at room temperature. For longer storage, samples should be frozen at -20°C.

## 7. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- The supplied reagents contain organic solvents. Although diluted, they still should be handled with care. It can cause burns and should be handled with gloves, eye protection and appropriate protective clothing. Do not breathe vapour and avoid inhalation.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.

## 8. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Do not use reagents beyond the expiration date stated on the kit label.
- Do not mix different lot numbers of any kit component within the same assay.
- Do not mix plugs and caps from different reagents.
- Follow the guidelines for medical laboratories.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be made within 14 days after receiving the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Analyse controls with each run.
- Always perform the assay according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- *IDK*<sup>®</sup> is a trademark of Immundiagnostik AG.

## 9. REFERENCES

1. Saupe J. (1992). Vitamin K und Indikationen zur Bestimmung im Plasma. Klin. Lab. 38: 51-54
2. Szulc P., Arlot M., Chapuy M.C., Duboef F., Meunier P.J., Delmas P.D. (1994). Serum undercarboxylated osteocalcin correlates with hip bone mineral density in elderly women. J Bone Miner Res 9:1591-95
3. Hodges S.J., Akesson K., Vergnaud P., Obrant K., Delmas P.D. (1993). Circulating levels of vitamin K1 and K2 decreased in elderly women with hip fracture. J Bone Miner Res 8:1241-45
4. Haroon Y, Bacon DS, Sadowski JA. (1986) Liquid-chromatographic determination of vitamin K1 in plasma, with fluorometric detection. Clin Chem. Oct;32(10):1925-9

## 10. SYMBOL EXPLANATION

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)