MutaPLATE[®] Faktor II (TM)

Real-Time-PCR-Kit

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MutaPLATE[®] Faktor II (TM)

Real-Time-PCR Kit

*For the analysis of the G20210A mutation
within the 3' UTR region of the prothrombin gene
(human coagulation factor II)*

Valid from 2022-03-04

KF1901632
KF1901696

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1	INTENDED USE	15
2	INTRODUCTION	15
3	PRINCIPLE OF THE TEST	15
4	PACKAGE CONTENTS	16
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	16
6	TRANSPORT, STORAGE AND STABILITY	16
7	WARNINGS AND PRECAUTIONS	16
8	SAMPLE MATERIAL	17
9	REAL-TIME-PCR	17
	9.1 <i>Important points before starting</i>	17
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	19
10	DATA ANALYSIS	19
11	TROUBLESHOOTING	20
12	LIMITATIONS OF THE METHOD	21
13	ABBREVIATIONS AND SYMBOLS	21
14	LITERATURE	22

1 INTENDED USE

The MutaPLATE® Faktor II (TM) Real-Time PCR Kit is a molecular biological test for the detection of the G20210A mutation within the 3' UTR region of the prothrombin gene (human coagulation factor II) in open real time PCR systems (e. g. RotorGene, SmartCycler, Light Cycler, ABI, Stratagene, Amplifa) by Taq-Man technology. This mutation is related to a big amount of prothrombin in plasma.

2 INTRODUCTION

Prothrombin (Factor II) is a protein in the blood required for the blood to clot. Blood clots are composed of a combination of blood platelets and a meshwork of the blood clotting protein fibrin. Prothrombin as blood clotting protein is needed to form fibrin. Persons with too little prothrombin have a bleeding tendency. If an individual has too much prothrombin, blood clots may form when they shouldn't. [1], [2], [3]

Poly-morphism	RsNumber	Effect		
		Wildtype	Heterozygous mutation	Homozygous mutation
G20210A	rs1799963	Normal levels of prothrombin	Intermediate increase in circulating prothrombin levels associated with higher thrombin formation. Prone to venous thrombosis.	Strong increase in circulating prothrombin levels associated with higher thrombin formation. Prone to venous thrombosis.

3 PRINCIPLE OF THE TEST

The MutaPLATE® Faktor II (TM) Real-Time PCR Kit contains specific primers and additional material for the detection of the G20210A polymorphism of the prothrombin gene. The variable area of the prothrombin gene is amplified by PCR using genomic DNA-template.

The standard PCR contains additionally two sequence specific oligonucleotides marked with fluorescence dye (TaqMan probes). Both probes bind at the amplified target-DNA which includes the single nucleotide polymorphism (SNP). Due to

this, a fluorescence signal is generated and detected by the optical unit of the used real time PCR instrument. The TaqMan probe for the G-allele (wildtype) is marked with **FAM (510 nm, green)** and the TaqMan probe for the A-allele (mutation) is marked with **YAK (555 nm, yellow)**.

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF1901632) or 96 (KF1901696) reactions.

Table 1: Components of the MutaPLATE® Faktor II (TM) Real-Time-PCR Kit .

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 438 µl	3 x 438 µl
Detection mix	yellow	1 x 368 µl	3 x 368 µl
Positive control	red	1 x 15 µl	1 x 45 µl
Negative control	green	1 x 150 µl	1 x 150 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/strips or tubes)
- Sterile PCR reaction tubes or 96-well plate/strips (white)
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaPLATE® Faktor II (TM) real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaPLATE® Faktor II (TM) real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.

- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, one reaction tube is needed per sample plus two additional reaction tubes for the negative and the positive control. The following tables show the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive control) (N) plus 10 % to compensate for inaccuracies. The master mixes are pipetted as described in tables 2 and 3:

Master mix

Table 2: Preparation of master mix

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (green).
- For the positive control add **2 µl** of the supplied positive control (red).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 4.

Table 4: Real-Time-PCR thermal profile

Description	Time	Temperature	Cycles	Acquisition
Initial Denaturation	120 s	94 °C	1	none
Denaturation	30 s	94 °C	45	none
Primer annealing and elongation	60 s	60 °C		single
Cooling	30 s	40 °C	1	-

10 DATA ANALYSIS

The results of the analysis for the **G20210A** polymorphisms are shown at **510-530 nm / green** and **550-560 nm / yellow** (corresponding to the Real-Time PCR instrument used).

The following figures show typical examples for both homozygous and heterozygous on a LightCycler® 2.0. With some real-time PCR devices (e.g. LightCycler® 2.0) a colour compensation file may be required.

G allele at 530 nm


Fig. 1: Evaluation of G allele

A allele at 560 nm


Fig. 2: Evaluation of A allele

The provided Positive Control (**red**) contains a template which is heterozygous for the G20210A-polymorphism (one allele carries the mutation, the other is wild type).

The following three discriminations are possible:

1. Homozygous **G/G**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe, no increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

2. Heterozygous **G/A**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe and increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

3. Homozygous **A/A**:

No increase of the fluorescent signal from the **FAM** labeled TaqMan probe, increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No fluorescence peak in the positive control or samples at about 510 - 530 nm or 550 - 560 nm.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8 °C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

Low fluorescence peak at about 510 - 530 nm or 555-560 nm

Mix individual components carefully before use (only by pipetting several times - do not vortex!).

Cool all stock solutions appropriately during the working steps and protect the detection mixes from light irradiation.

Work on ice or with a cooled block (4 °C).

12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix		Upper limit of temperature

	Positive control		Manufacturer
	Negative control		Lot number
	Content		Consult instructions for use
	Use by YYYY-MM-DD		<i>In vitro</i> diagnostic medical device

14 LITERATURE

1. V. D. Stefano, I. Martinelli, P. M. Mannucci, K. Paciaroni, P. Chiusolo, I. Casorelli, E. Rossi und G. Leone, „The risk of recurrent deep venous thrombosis among heterozygous carriers of both factor V Leiden and the G20210A prothrombin mutation,“ *N Engl J Med*, Bd. 341, Nr. 11, pp. 801-806, 1999.
2. A. Gerhardt, R. E. Scharf, M. W. Beckmann, S. Struve, H. G. Bender, M. Pillny, W. Sandmann und R. B. Zotz, „Prothrombin and factor V mutations in women with a history of thrombosis during pregnancy and the puerperium,“ *N Engl J Med*, Bd. 342, Nr. 6, pp. 374-380, 2000.
3. R. F. Franco, M. D. Trip, H. t. Cate, A. v. den, M. H. Prins, J. J. Kastelein und P. H. Reitsma, „The 20210 G-->A mutation in the 3'-untranslated region of the prothrombin gene and the risk for arterial thrombotic disease,“ *Br J Haematol*, Bd. 104, Nr. 1, pp. 50-54, 1999.

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com