# MutaPLATE<sup>®</sup> Aldolase B (TM)

## Real-Time-PCR-Kit

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# MutaPLATE<sup>®</sup> Aldolase B (TM)

## Real-Time-PCR Kit

*For the analysis of the mutations A149P und N334K  
within the Aldolase B gene*

Valid from 2023-12-06

**KF190432**  
**KF190496**

32/96

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1</b>	<b>INTENDED USE</b>	<b>19</b>
<b>2</b>	<b>INTRODUCTION</b>	<b>19</b>
<b>3</b>	<b>PRINCIPLE OF THE TEST</b>	<b>19</b>
<b>4</b>	<b>PACKAGE CONTENTS</b>	<b>19</b>
<b>5</b>	<b>EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER</b>	<b>20</b>
<b>6</b>	<b>TRANSPORT, STORAGE AND STABILITY</b>	<b>20</b>
<b>7</b>	<b>WARNINGS AND PRECAUTIONS</b>	<b>20</b>
<b>8</b>	<b>SAMPLE MATERIAL</b>	<b>21</b>
<b>9</b>	<b>REAL-TIME-PCR</b>	<b>21</b>
	9.1 <i>Important points before starting</i>	21
	9.2 <i>Procedure</i>	22
	9.3 <i>Instrument settings</i>	23
<b>10</b>	<b>DATA ANALYSIS</b>	<b>23</b>
<b>11</b>	<b>TROUBLESHOOTING</b>	<b>28</b>
<b>12</b>	<b>LIMITATIONS OF THE METHOD</b>	<b>29</b>
<b>13</b>	<b>ABBREVIATIONS AND SYMBOLS</b>	<b>29</b>
<b>14</b>	<b>LITERATURE</b>	<b>30</b>

## 1 INTENDED USE

The MutaPLATE® Aldolase B (TM) Real-Time PCR Kit is a molecular biological assay for the detection of point mutations A149P and N334K in the aldolase B gene in open real-time PCR systems using Taq-Man technology.

## 2 INTRODUCTION

The enzyme aldolase B cleaves the hexose fructose within the liver into the two tri-saccharides glyceraldehyde and dihydroxyacetone phosphate, which can then be further metabolized. The hereditary fructose intolerance is based on a defect within the enzyme aldolase B, which disturbs the degradation of fructose. Therefore fructose is accumulated and can have a harmful effect on the organism. [1]

## 3 PRINCIPLE OF THE TEST

The MutaPLATE® Aldolase B (TM) Real-Time PCR Kit contains two sequence specific primers flanking the region of interest and two TaqMan probes specific to the region containing the mutation. The two TaqMan probes are labeled at the 5' end with different fluorophores (reporter dyes) which are used for the allelic discrimination. On the 3' end the TaqMan probes are labeled with a non-fluorescent quencher. The proximity of the reporter dye to the quencher inhibits the fluorescence of the reporter molecule.

During amplification the probes hybridize specifically to the DNA fragments. The 5' nuclease activity of the polymerase cleaves the hybridized probes releasing the reporter from the quencher generating a fluorescent signal.

## 4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF190432) or 96 (KF190496) reactions.

Table 1: Components of the MutaPLATE® Aldolase B (TM) Real-Time-PCR Kit .

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 875 µl	2 x 1315 µl
Detection mix 1 (A149P)	yellow	1 x 368 µl	2 x 552 µl
Detection mix 2 (N334K)	white	1 x 368 µl	2 x 552 µl
Positive control 1 (A149P)	red	1 x 15 µl	1 x 45 µl
Positive control 2 (N334K)	orange	1 x 15 µl	1 x 45 µl
Negative control	green	1 x 150 µl	1 x 150 µl

## 5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/stripes or tubes)
- Sterile PCR reaction tubes or 96-well plate/stripes (white)
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

## 6 TRANSPORT, STORAGE AND STABILITY

The MutaPLATE® Aldolase B (TM) real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20°C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

## 7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).

- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

## 8 SAMPLE MATERIAL

Starting material for the MutaPLATE® Aldolase B (TM) real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

## 9 REAL-TIME-PCR

### 9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.

- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

## 9.2 Procedure

For amplification, two reaction tubes are required per sample and two additional reaction tubes per master mix are required for the negative and the positive controls. The following tables show the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive control) (N) plus 10 % to compensate for inaccuracies. The master mixes are pipetted as described in tables 2 and 3:

### Master mix 1 (A149P)

Table 2: Preparation of master mix

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 1 (yellow)	10.5 µl	$10.5 \mu\text{l} * (N + (N * 0.1))$
Enzyme mix (blue)	12.5 µl	$12.5 \mu\text{l} * (N + (N * 0.1))$

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (green).
- For the positive control add **2 µl** of the supplied positive control 1 (red).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

### Master mix 2 (N334K)

Table 3: Preparation of master mix

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 2 (white)	10.5 µl	$10.5 \mu\text{l} * (N + (N * 0.1))$
Enzyme mix (blue)	12.5 µl	$12.5 \mu\text{l} * (N + (N * 0.1))$

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.

- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 2 (**orange**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

Close the reaction tubes and shortly spin down. Then transfer to the real-time device and start the PCR protocol described in 9.3.

### 9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 3.

Table 4: Real-Time-PCR thermal profile

Description	Time	Temperature	Cycles	Acquisition
Initial Denaturation	120 s	94 °C	1	none
Denaturation	10 s	94 °C	45	none
Primer annealing and elongation	30 s	66 °C		<b>single</b>

## 10 DATA ANALYSIS

The TaqMan probes for the A149 and the N334 allele (wild type) are labeled with **FAM (green)**. The TaqMan probes for the 149P and the 334K allele (mutation) are labeled with **YAK (yellow)**.

The evaluation of the amplification curves (determination of crossing points) is carried out with an analysis of the type "absolute quantification". With some real-time PCR devices, a colour compensation file may be required.

To evaluate a **Bio-Rad CFX 96 / CFX Opus 96** run, we recommend the following settings in the CFX Maestro software:

- Cq Determination Mode: Single Treshold (100 RFU for A149P / 200 RFU for N334K).
- Baseline Setting: Baseline Substracted Curve Fit + Apply Fluorescence Drift Correction
- Cycles to Analyze: 10 to 45

## Evaluation A149P of a Bio-Rad CFX 96 / CFX Opus 96 run

The following graphs show a representative result for the A149P mutation:

### FAM (green) - detection of the A allele (wild type)


Fig. 1: Evaluation of Aldolase B A149P A allele - Bio-Rad CFX 96 / CFX Opus 96.

### YAK (yellow) - Detection of the P allele (mutation)


Fig. 2: Evaluation of Aldolase B A149P P allele - Bio-Rad CFX 96 / CFX Opus 96.

### Evaluation N334K of a Bio-Rad CFX96 / CFX Opus 96 run

The following graphs show a representative result for the N334K mutation:

#### FAM (green) - detection of the N allele (wild type)


Fig. 3: Evaluation of Aldolase B N334K N allele - Bio-Rad CFX 96 / CFX Opus 96.

#### YAK (yellow) - Detection of the K allele (mutation)


Fig. 4: Evaluation of Aldolase B N334K K allele - Bio-Rad CFX 96 / CFX Opus 96.

To evaluate a **Roche LightCycler® 480 II** run, we recommend the following settings in the software:

- Abs Quant/2nd Derivative Max
- Cycles to Analyze: 10 to 45

### Evaluation A149P of a Roche LightCycler® 480 II run

The following graphs show a representative result for the A149P mutation:

**FAM (green)** - detection of the A allele (wild type)


Fig. 5: Evaluation of Aldolase B A149P A allele - Roche LightCycler® 480 II.

**YAK (yellow)** - Detection of the P allele (mutation)


Fig. 6: Evaluation of Aldolase B A149P P allele - Roche LightCycler® 480 II.

## Evaluation N334K of a Roche LightCycler® 480 II run

The following graphs show a representative result for the N334K mutation:

**FAM (green)** - detection of the N allele (wild type)


Fig. 7: Evaluation of Aldolase B N334K N allele - Roche LightCycler® 480 II.

**YAK (yellow)** - Detection of the K allele (mutation)


Fig. 8: Evaluation of Aldolase B N334K K allele - Roche LightCycler® 480 II.

The provided Positive Control 1 (A149P) (**red**) contains a template, which is heterozygous for the point mutation A149P. The provided Positive Control 2 (N334K) (**orange**) contains a template, which is heterozygous for the point mutation N334K.

The following three discriminations are possible:

1. Homozygous **wild type**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe, no increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

2. Heterozygous **mutation**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe and increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

3. Homozygous **mutation**:

No increase of the fluorescent signal from the **FAM** labeled TaqMan probe, increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

## 11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on [info@immundiagnostik.com](mailto:info@immundiagnostik.com).

### **No fluorescence peak in the positive control or samples in the FAM or YAK channel:**

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8°C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

### **Low fluorescence peak in the positive control or samples in the FAM or YAK channel:**

Mix individual components carefully before use (only by pipetting several times - do not vortex!).

Cool all stock solutions appropriately during the working steps and protect the detection mixes from light irradiation.

Work on ice or with a cooled block (4°C).


## 12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

## 13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix 1		Upper limit of temperature
	Detection mix 2		Manufacturer
	Positive control 1		Lot number
	Positive control 2		Consult instructions for use
	Negative control		Use by YYYY-MM-DD
	Content		<i>In vitro</i> diagnostic medical device

## **14 LITERATURE**

1. Sánchez-Guitérrez et al., J Med Genet, 2002, 39:e56

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)