MutaPLATE[®] HLA DQ 2+8 (TM)

Real-Time-PCR-Kit

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www.immundiagnostik.com

Table of Contents

1	INTENDED USE	15
2	INTRODUCTION	15
3	PRINCIPLE OF THE TEST	15
4	PACKAGE CONTENTS	16
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	16
6	TRANSPORT, STORAGE AND STABILITY	16
7	WARNINGS AND PRECAUTIONS	17
8	SAMPLE MATERIAL	18
9	REAL-TIME-PCR	18
	9.1 <i>Important points before starting</i>	18
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	19
10	DATA ANALYSIS	20
11	TROUBLESHOOTING	23
12	LIMITATIONS OF THE METHOD	23
13	ABBREVIATIONS AND SYMBOLS	24
14	LITERATURE	24

1 INTENDED USE

The MutaPLATE® HLA DQ 2+8 (TM) Real-Time PCR Kit is a molecular biological assay for the detection of the HLA alleles DQA1*05, DQA1*03, DQB1*02¹ and DQB1*03:02/*03:05^{2,3} based on the TaqMan-technology.

¹ The number of copies of DQB1*02 is not determined by this method.

² Due to sequence homologies this method cannot differentiate between DQB1*03:02 and DQB1*03:05. However, the frequency of DQB1*03:05 is very low, about 0.4% in the Caucasian population, whereas DQB1*03:02 has a frequency of about 15% (<http://www.allelefrequencies.net>, Nov. 2018). Worldwide the frequency of DQB1*03:05 is 0.07% and the frequency of DQB1*03:02 is 11.5% (Solberg et al., Hum Immunol, 2008).

³ The detection of other very rare alleles cannot be completely excluded.

2 INTRODUCTION

Celiac disease / gluten intolerance is one of the most frequent chronic gastrointestinal diseases. It is a genetic disease, in which the body is not able to process the gluten present in many cereals. Almost all celiac disease patients are carrier of HLA-DQ2 (HLADQA1*05 and HLA-DQB1*02) or HLA-DQ8 (HLA-DQB1*03:02, often in combination with DQA1*03). If DQ2 and DQ8 are not detected, celiac disease can be excluded with a probability of over 95%.

3 PRINCIPLE OF THE TEST

The MutaPLATE® HLA DQ 2+8 (TM) real-time PCR kit contains two specific primers that flank the target sequence and a hydrolysis probe (TaqMan probes) that binds allele specifically. In addition, the kit includes an internal amplification control to indicate correct PCR if no HLA allele is present in the sample.

The hydrolysis probes are labelled at the 5' end with different fluorophores (reporter dyes), which are used to distinguish the alleles. At the 3' end the probes are labelled with a non-fluorescent quencher. The proximity of the reporter dye to the quencher inhibits the fluorescence of the reporter molecule. During amplification, the probes bind specifically to the DNA fragments. The 5' nuclease activity of the polymerase cleaves the hybridised probes, separating the reporter from the quencher and generating a fluorescent signal.

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF190532+) or 96 (KF190596+) reactions.

Table 1: Components of the MutaPLATE® HLA DQ 2+8 (TM) Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Reaction mix 1 (DQA1*05 / DQB1*02)	yellow	1 x 805 µl	3 x 805 µl
Reaction mix 2 (DQA1*03 / DQB1*03:02)	brown	1 x 805 µl	3 x 805 µl
Positive control	red	1 x 30 µl	1 x 90 µl
Negative control	green	1 x 150 µl	1 x 150 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/strips or tubes), e.g.: Bio-Rad CFX 96 or CFX Opus 96 or Roche LightCycler® 480 II.
 - *The CE conformity is only given with these instruments.*
- Sterile PCR reaction tubes or 96-well plates/strips (white)
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaPLATE® HLA DQ 2+8 (TM) real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid more than eight freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the reaction mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit components.
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaPLATE® HLA DQ 2+8 (TM) real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 Important points before starting

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the reaction mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, two reaction tubes are required per sample and two additional reaction tubes per reaction mix are required for the negative and the positive control.

The Reaction Mixes contain all of the components needed for PCR except the sample:

Reaction Mix 1: DQA1*05 / DQB1*02

- Mix the Reaction Mix 1 carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Reaction Mix 1 to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Reaction Mix 2: DQA1*03 / DQB1*03:02

- Mix the Reaction Mix 2 carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Reaction Mix 2 to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Thoroughly mix the PCR reaction tubes, close the reaction tubes and briefly centrifuge. Then transfer to the real-time PCR device and start the PCR programme described in 9.3..

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 2.

Table 2: Real-Time-PCR thermal profile

Description	Time	Temperature	Ramp rate	Cycles	Acquisition
Initial Denaturation	120 s	94 °C	max.	1	none
Denaturation	10 s	94 °C	max.	45	none
Primer annealing and elongation	50 s	60 °C	max.		single
Cooling	30 s	40 °C	max.	1	-

10 DATA ANALYSIS

The MutaPLATE® HLA DQ 2+8 Real-Time PCR kit detects the presence of the HLA alleles DQA1*05 DQA1*03, DQB1*02 and DQB1*03:02. The corresponding Taq-Man probes for the four alleles are FAM or ROX labelled. If the HLA alleles are not present, no amplification takes place and thus no fluorescence is detected by the FAM or ROX labelled TaqMan probe. In order to guarantee a successful PCR in this case, an internal amplification control (IC) is incorporated into the PCR. The TaqMan probe for the internal amplification control is measured on the YAK channel.

Corresponding to the genotype the following results can be achieved:

1. HLA allele present:

Increase of the fluorescence signal of the **FAM** resp. **ROX**-labeled TaqMan probe and increase of the fluorescence signal of the **YAK**-labeled TaqMan probe.

2. HLA allele is not present:

No increase of the fluorescence signal of the **FAM** resp. **ROX**-labeled TaqMan probe and increase of the fluorescence signal of the **YAK**-labeled TaqMan probe.

For the **evaluation** of a **Bio-Rad CFX 96 / CFX Opus 96** run we recommend the following settings in the CFX Maestro software:

- Cq Determination Mode: Regression
- Baseline Setting: Baseline Substracted Curve Fit + Apply Fluorescence Drift Correction
- Cycles to Analyze: 10 to 45

Reaction Mix 1:

FAM - channel: Detection of the allele DQA1*05


Fig. 1: Evaluation for HLA DQA1*05

YAK - channel: Internal amplification control


Fig. 2: Evaluation for Reaction Mix 1 - Internal amplification control

ROX - channel: Detection of the allele DQB1*02


Fig. 3: Evaluation for HLA DQB1*02.

Reaction Mix 2:

FAM - channel: Detection of the allele DQB1*03:02


Fig. 4: Evaluation for HLA DQB1*03:02

YAK - channel - Internal amplification control


Fig. 5: Evaluation for Reaction Mix 2 - Internal amplification control

ROX - channel: Detection of the allele DQA1*03


Fig. 6: Evaluation for HLA DQA1*03.

The provided positive control (**red**) contains a template, which is positive for the alleles DQA1*05, DQA1*03, DQB1*02 and DQB1*03:02.

The negative control (**green**) has to be negative in all channels.

11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No fluorescence peak in the positive control or samples in the FAM, YAK or ROX channel:

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Reaction mixes have been subjected to more than eight freeze cycles. Repeat the analysis with a fresh aliquot or new reaction mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

Low fluorescence peak in the positive control or samples in the FAM, YAK or ROX channel:

Mix individual components carefully before use (only by pipetting several times - do not vortex!).

Cool all stock solutions appropriately during the working steps and protect the reaction mixes from light irradiation.

Work on ice or with a cooled block (4 °C).

12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Reaction mix 1		Contains sufficient for <n> test
	Reaction mix 2		Upper limit of temperature
	Positive control		Manufacturer
	Negative control		Lot number
	<i>In vitro</i> diagnostic medical device		Consult instructions for use
	Content		Use by YYYY-MM-DD

14 LITERATURE

1. Solberg et al., Hum Immunol, 2008

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com