# MutaPLATE® COMT

## Real-Time-PCR-Kit

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# MutaPLATE<sup>®</sup> COMT

## Real-Time-PCR Kit

*For the analysis of the Val158Met mutation in the  
COMT gene.*

Valid from 2024-01-31

**KF191532**  
**KF191596**

32/96

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1</b>	<b>INTENDED USE</b>	<b>15</b>
<b>2</b>	<b>INTRODUCTION</b>	<b>15</b>
<b>3</b>	<b>PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>4</b>	<b>PACKAGE CONTENTS</b>	<b>16</b>
<b>5</b>	<b>EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER</b>	<b>16</b>
<b>6</b>	<b>TRANSPORT, STORAGE AND STABILITY</b>	<b>16</b>
<b>7</b>	<b>WARNINGS AND PRECAUTIONS</b>	<b>16</b>
<b>8</b>	<b>SAMPLE MATERIAL</b>	<b>17</b>
<b>9</b>	<b>REAL-TIME-PCR</b>	<b>18</b>
	9.1 <i>Important points before starting</i>	18
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	19
<b>10</b>	<b>DATA ANALYSIS</b>	<b>19</b>
<b>11</b>	<b>TROUBLESHOOTING</b>	<b>21</b>
<b>12</b>	<b>LIMITATIONS OF THE METHOD</b>	<b>21</b>
<b>13</b>	<b>ABBREVIATIONS AND SYMBOLS</b>	<b>22</b>
<b>14</b>	<b>LITERATURE</b>	<b>22</b>

## 1 INTENDED USE

The MutaPLATE® COMT Real-Time PCR Kit is a molecular biology test for the investigation of the Val158Met mutation in the catechol-O-methyl transferase (COMT) gene in open real-time PCR systems using Taq-Man technology.

## 2 INTRODUCTION

As a phase II enzyme, catechol-O-methyltransferase catalyses the transfer of methyl groups to endogenous/exogenous catechols and is thus involved in their inactivation. For example, the O-methylation of neurotransmitters such as dopamine, adrenalin or noradrenaline is a main degradation pathway of these catecholamines. In addition to the role in the metabolism of endogenous substances, the enzyme is also important in the metabolism of catecholamines, which are used, for example, in the treatment of high blood pressure, asthma and Parkinson's disease.

The exchange of valine for methionine at position 158 of catechol-O-methyltransferase results in a 3-4-fold reduction in the activity (low activity allele) of the enzyme (in the case of a homozygous genotype). [1], [2]

## 3 PRINCIPLE OF THE TEST

The MutaPLATE® COMT real-time PCR kit contains sequence specific primers flanking the region of interest and TaqMan probes specific to the region containing the mutation. The two TaqMan probes are labeled at the 5' end with different fluorophores (reporter dyes) which are used for the allelic discrimination. On the 3' end the TaqMan probes are labeled with a non-fluorescent quencher. The proximity of the reporter dye to the quencher inhibits the fluorescence of the reporter molecule.

During amplification the probes hybridize specifically to the DNA fragments. The 5' nuclease activity of the polymerase cleaves the hybridized probes releasing the reporter from the quencher generating a fluorescent signal.

## 4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF191532) or 96 (KF191596) reactions.

Table 1: Components of the MutaPLATE® COMT Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 438 µl	3 x 438 µl
Detection mix	yellow	1 x 368 µl	3 x 368 µl
Positive control	red	1 x 15 µl	1 x 45 µl
Negative control	green	1 x 150 µl	1 x 150 µl

## 5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/strips or tubes)
- Sterile PCR reaction tubes or 96-well plate/strips (white)
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

## 6 TRANSPORT, STORAGE AND STABILITY

The MutaPLATE® COMT real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

## 7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.

- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

## 8 SAMPLE MATERIAL

Starting material for the MutaPLATE® COMT real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

**ATTENTION:** The concentration of the DNA used should not exceed 30 ng/μl. Otherwise, cross-talk between the FAM or YAK channel may occur and thus influence the interpretation.

## 9 REAL-TIME-PCR

### 9.1 Important points before starting

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run the Positive Control and the Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mix from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

### 9.2 Procedure

For amplification, one reaction tube is needed per sample plus two additional reaction tubes per for the negative and the positive control. The following tables show the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mix is pipetted as described in table 2:

#### Master mix - COMT

Table 2: Preparation of master mix - COMT

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (green).

- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Close the reaction tubes and shortly spin down. Then transfer to the real-time device and start the PCR protocol described in 9.3.

### 9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 3.

Table 3: Real-Time-PCR thermal profile

Description	Time	Temperature	Cycles	Acquisition
Initial Denaturation	180 s	95 °C	1	none
Denaturation	10 s	95 °C	45	none
Primer annealing and Elongation	30 s	63 °C		<b>single</b>

## 10 DATA ANALYSIS

The TaqMan probe for wild type allele is labeled with **FAM (green)**. The TaqMan probe for the mutation allele are labeled with **YAK (yellow)**.

The evaluation of the amplification curves (determination of crossing points) is carried out with an analysis of the type "absolute quantification". With some real-time PCR devices, a colour compensation file may be required.

To evaluate a **Bio-Rad CFX 96 / CFX Opus 96** run, we recommend the following settings in the CFX Maestro software:

- Cq Determination Mode: Single Treshold (210 RFU).
- Baseline Setting: Baseline Substracted Curve Fit + Apply Fluorescence Drift Correction
- Cycles to Analyze: 10 to 45

### COMT\*Val158Met (Bio-Rad CFX Opus 96)

The following figures show the typical results for this polymorphism.

#### FAM (green) - Detection of the wild type allele


Fig. 1: Evaluation of COMT\*Val158Met – wild typ allele - Bio-Rad CFX Opus 96

#### HEX (yellow) - Detection of the mutation allele


Fig. 2: Evaluation of COMT\*Val158Met – mutation allele - Bio-Rad CFX Opus 96

The positive control (**red**) provided contains a template that is **heterozygous** for the mutation.

## 11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on [info@immundiagnostik.com](mailto:info@immundiagnostik.com).

### **No or weak fluorescence in the positive control or samples.**

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mix has been subjected to more than two freeze cycles or has been stored at 2 - 8 °C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mix was not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

## 12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test analyses only a selection of markers. When alleles are detected, the investigated polymorphism is indicated. Other rare alleles may be present and are not covered by this method. Therefore, a negative test result for a patient does not completely exclude a risk of any kind.

### 13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix		Upper limit of temperature
	Positive control		Manufacturer
	Negative control		Lot number
	<i>In vitro</i> diagnostic medical device		Consult instructions for use
	Use by YYYY-MM-DD		Content

### 14 LITERATURE

1. Milikan et al., Carcinogenesis 199 (11), 1998
2. Egan et al., PANS 98, 18917-22, 1998

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)