MutaREAL® HFE

Real-Time-PCR-Kit

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14 LITERATUR

1. Robson et al. (2000) DIAGNOSIS AND MANAGEMENT OF HAEMOCHROMATOSIS SINCE THE DISCOVERY OF THE HFE GENE: A EUROPEAN EXPERIENCE, British Journal of Haematology, 2000, 108, 31±39

MutaREAL[®] HFE

Real-Time-PCR Kit

*For the analysis of the three most important mutations
(C282Y, H63D and S65C) of the haemochromatosis
(HFE) gene*

Valid from 2025-01-07

KF291032
KF291096

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1	INTENDED USE	17
2	INTRODUCTION	17
3	PRINCIPLE OF THE TEST	17
4	PACKAGE CONTENTS	17
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	18
6	TRANSPORT, STORAGE AND STABILITY	18
7	WARNINGS AND PRECAUTIONS	18
8	SAMPLE MATERIAL	19
9	REAL-TIME-PCR	19
	9.1 <i>Important points before starting</i>	19
	9.2 <i>Procedure</i>	20
	9.3 <i>Instrument settings</i>	21
10	DATA ANALYSIS	21
11	TROUBLESHOOTING	25
12	LIMITATIONS OF THE METHOD	25
13	ABBREVIATIONS AND SYMBOLS	26
14	LITERATURE	26

1 INTENDED USE

The MutaREAL® HFE Real-Time PCR Kit is a FRET-based molecular biology test for the investigation of the point mutations C282Y, H63D and S65C in the HFE gene.

2 INTRODUCTION

Primary haemochromatosis is an inherited disease in which there is excessive absorption of iron into the body. absorption of iron into the body. This can lead to liver damage, diabetes and heart problems, among other things, diabetes and heart problems. Mutations in the gene for HFE can lead to haemochromatosis. Wild-type HFE interacts with β 2-microglobulin and regulates iron uptake via the transferrin- receptor. In the case of an HFE mutation, iron uptake can be impaired, causing iron to accumulate in the body. iron accumulates in the body. Homozygous as well as, in rarer cases, also heterozygous presence of mutations can lead to an outbreak of the disease. [1]

3 PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® HFE real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay contains two specific primers that flank the target sequence and two hybridisation probes that bind to the target sequences. The amount of hybridised probe and thus the fluorescence signal increases with the amount of amplified PCR product. The fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of the amplification by means of a melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. The hybridisation probes bind to the part of the target sequence that is present in the wild type and the mutation. As the temperature increases, the mismatched and therefore less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and therefore the fluorescence signal only decreases at a higher temperature.

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF291032) or 96 (KF291096) reactions.

Table 1: Components of the MutaREAL® HFE Real-Time-PCR Kit .

Label	Lid Colour	Content	
		32	96
Detection mix 1 (C282Y)	yellow	1 x 805 µl	3 x 805 µl
Detection mix 2 (H63D & S65C)	white	1 x 805 µl	3 x 805 µl
Positive control 1 (H63D & C282Y)	red	1 x 30 µl	1 x 90 µl
Positive control 2 (S65C)	orange	1 x 15 µl	1 x 45 µl
Negative control	green	1 x 150 µl	1 x 150 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/strips or tubes)
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaREAL® HFE real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20°C immediately after receipt. Avoid more than 4 freeze-thaw cycles (if necessary, make aliquots). Do not use after the expiry date indicated on the package.

After opening, the detection mixes can be stored at +4 to +8°C for up to 3 months. Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.

- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaREAL® HFE real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.

- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, two reaction tubes are needed per sample plus two or three additional reaction tubes per detection mix for the negative and the positive controls.

Detection mix 1 (C282Y)

- Mix the Detection Mix 1 (**yellow**) carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Detection Mix 1 (**yellow**) to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 1 (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Detection mix 2 (H63D & S65C)

- Mix the Detection Mix 2 (**white**) carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Detection Mix 2 (**white**) to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** each of the supplied positive control 1 (**red**) and 2 (**orange**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Close the reaction tubes, mix carefully and shortly spin down. Then transfer to the real-time device and start the PCR protocol described in 9.3.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 2.

Table 2: Real-Time-PCR thermal profile

Description	Time	Temperature	Cycles	Acquisition
Initial Denaturation	180 s	95 °C	1	none
Denaturation	10 s	95 °C	45	none
Primer annealing	15 s	53 °C		single
Elongation	15 s	72 °C		none
Melting curve	20 s	94 °C	1	none
	20 s	40 °C	1	none
	0 s	75 °C	1	constant
Cooling	30 s	40 °C	1	-

10 DATA ANALYSIS

Add an analysis of the type „genotyping“ for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is 530 nm for all three mutations.

H63D and S65C

For H63D, the probe is designed in such a way that the H63D mutation exactly matches the mutation and has a mismatch with the wild type. In addition, the probe is designed to detect the S65C genotype at 47 °C.

In the very rare case of heterozygous mutations of H63D and S65C, there is no wild-type peak, because due to the Simple-Probe method used and the spatial proximity of these two mutations, only the exactly matching mutation probes can bind. The „access“ for the wild-type probe is denied. Nevertheless, the result is similar to that of a (not described so far) homozygous double mutation of H63D + S65C.

Due to the base difference of only 4 bp in the SNPs for H63D and S65C, a joint probe must be used for both SNPs. The very rare combination of a heterozygous H63D and S65C allele carrier leads to an omitted wild-type peak, because the FRET technology only generates two mutation peaks. Since this melting profile then looks like a homozygous H63D- and S65C- allele carrier a stepwise diagnosis of this sample is recommended, (e.g. through sequencing) to exclude homozygosity of the respective SNP.

On the other hand, such a double homozygous genotype has not yet been described.

Figures 1-3 show typical examples of homozygous and heterozygous samples. Temperatures should be in the range +/- 2 °C*:

H63D & S65C

Temperature mutation allele (S65C)*: 47.0 °C (+/-2 °C)

Temperature wild type allele*: 52.0 °C (+/-2 °C)

Temperature mutation allele (H63D)*: 61.0 °C (+/-2 °C)


Fig. 1: Example figure from *Clin Chem Lab Med* 2008;46(7):985-986 showing all known genotypes in H63D and S65C analysis using the SimpleProbe method.

The following figures show the typical results for some of the possible genotypes for H63D.


Fig. 2: Evaluation of the mutations H63D & S65C on a Roche LightCycler® 2.0.


Fig. 3: Evaluation of the mutations H63D & S65C on a Bio-Rad CFX Opus 96.

C282Y

Temperature wild type allele*: 60.0°C (+/-2°C)

Temperature mutation allele*: 49.0°C (+/-2°C)

The following figures show the typical results for some of the possible genotypes for C282Y.


Fig. 4: Evaluation of the mutation C282Y on a Roche LightCycler® 2.0.


Fig. 5: Evaluation of the mutation C282Y on a Bio-Rad CFX Opus 96.

** Depending on the PCR device used, the T_m values determined may deviate by more than +/- 2 °C. In this case use the positive controls (red and orange) as a reference and apply +/- 2 °C to these T_m values.*

The supplied positive control 1 (red) contains a template that is **heterozygous** for the mutations **H63D & C282Y**. The positive control 2 (orange) contains a template that is **heterozygous** for the mutation **S65C**.

11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than four freeze cycles or have been stored at 2-8 °C for more than 3 month. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.


12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. When alleles are detected, the investigated polymorphism is indicated. Other rare alleles may be present and are not covered by this method. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Detection mix 1		Contains sufficient for <n> test
	Detection mix 2		Upper limit of temperature
	Positive control 1		Manufacturer
	Positive control 2		Lot number
	Negative control		Consult instructions for use
	<i>In vitro</i> diagnostic medical device		Content
	Use by YYYY-MM-DD		

14 LITERATURE

1. Robson et al. (2000) DIAGNOSIS AND MANAGEMENT OF HAEMOCHROMATOSIS SINCE THE DISCOVERY OF THE HFE GENE: A EUROPEAN EXPERIENCE, British Journal of Haematology, 2000, 108, 31±39

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com