MutaREAL[®] DPD-2

Real-Time-PCR-Kit

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MutaREAL[®] DPD-2

Real-Time-PCR Kit

*For the analysis of the polymorphisms *2A, *13, *D949V and the haplotype B3 in the DPYD gene*

Valid from 2024-03-04

KF292032
KF292096

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1	INTENDED USE	15
2	INTRODUCTION	15
3	PRINCIPLE OF THE TEST	15
4	PACKAGE CONTENTS	16
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	16
6	TRANSPORT, STORAGE AND STABILITY	17
7	WARNINGS AND PRECAUTIONS	17
8	SAMPLE MATERIAL	18
9	REAL-TIME-PCR	18
	9.1 <i>Important points before starting</i>	18
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	20
10	DATA ANALYSIS	20
11	TROUBLESHOOTING	22
12	LIMITATIONS OF THE METHOD	23
13	ABBREVIATIONS AND SYMBOLS	23
14	LITERATURE	24

1 INTENDED USE

The MutaREAL® DPD-2 Real-Time PCR Kit is a molecular biological test for the investigation of mutations and polymorphisms of the dihydropyrimidine dehydrogenase gene (DPYD) from genomic DNA. The test is based on FRET-based real-time PCR. The gene variants investigated are associated with dihydropyrimidine dehydrogenase deficiency, which has an influence on the metabolism of 5-Fluorouracil (5-FU). The following alleles and mutations are detected by the test:

Table 1: Detected alleles and mutations.

Allele / Mutation		rs number	Influence on the 5-FU metabolism
*2A	IVS14+1G>A	rs3918290	very high
*13	1679T>G	rs55886062	very high
D949V	2846A>T	rs67376798	very high
Haplotype B3	1236G>A	rs56038477	high

2 INTRODUCTION

The enzyme dihydropyrimidine dehydrogenase (DPD) belongs to the group of oxidoreductases and catalyses the degradation of various cytostatic drugs such as 5-fluorouracil (5-FU). Mutations in the DPYD gene lead to reduced activity of the enzyme. This can lead to an accumulation of the cytostatic drug in the body during treatment with 5-FU. This can lead to severe side effects and life-threatening poisoning. [1]

3 PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® DPD-2 real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay contains specific primers that flank the respective target sequence and hybridisation probes that bind to the target sequence in the region of the SNP. One of the hybridisation probes is labelled with a donor fluorophore and, after appropriate excitation, transfers its energy to the acceptor fluorophore with which the other hybridisation probe is labelled if they are in close proximity. After the energy transfer, the acceptor dye emits light with a longer wavelength. Energy transfer can only occur when both hybridisation probes have bound to the target sequence. The amount of hybridised probe pairs and thus the fluorescence signal increases with the amount of amplified PCR product. Here, the fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of amplification by melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. One of the hybridisation probes binds to a part of the target sequence that is present in the wild type and the mutation. The second hybridisation probe spans the mutation site. As the temperature rises, the mismatched and thus less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and thus the fluorescence signal decreases only at a higher temperature.

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF292032) or 96 (KF292096) reactions.

Table 2: Components of the MutaREAL® DPD-2 Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Enzyme / buffer mix	blue	1 x 876 µl	2 x 1314 µl
Detection mix 1 (DPD*2A / *Haplo B3)	yellow	1 x 368 µl	2 x 552 µl
Detection mix 2 (DPD*13 / *D949V)	white	1 x 368 µl	2 x 552 µl
Positive control 1 (DPD*2A)	red	1 x 30 µl	1 x 90 µl
Positive control 2 (DPD*Haplo B3)	orange	1 x 30 µl	1 x 90 µl
Positive control 3 (DPD*D949V)	purple	1 x 30 µl	1 x 90 µl
Positive control 4 (DPD*13)	brown	1 x 30 µl	1 x 90 µl
Negative control	green	1 x 150 µl	1 x 150 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Roche LightCycler® 1.2, 1.5 or 2.0 real-time PCR system
- Roche LightCycler® capillaries
- Roche LightCycler® Cooling Block
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaREAL® DPD-2 real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

After opening, the detection mixes and the enzyme / buffer mix can be stored at +2 to +8 °C for up to 6 months.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.

- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaREAL® DPD-2 real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 Important points before starting

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8°C) or on ice while working.

9.2 Procedure

For amplification, two reaction tubes (LightCycler® capillaries) per sample and six additional reaction tubes per master mix are required for the negative (2x) and the positive (4x) controls. The following table shows the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mix is pipetted as described in tables 3 and 4:

Master mix 1 (*2A / *Haplo B3)

Table 3: Preparation of master mix 1 (*2A / *Haplo B3)

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 1 (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme / buffer mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control 1 add **2 µl** of the supplied positive control 1 (**red**).
- For the positive control 2, add **2 µl** of the supplied positive control 2 (**orange**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Master mix 2 (*13 / *D949V)

Table 4: Preparation of master mix 2 (*13 / *D949V)

Reagenz	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 2 (white)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme / buffer mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control 3 add **2 µl** of the supplied positive control 3 (**purple**).
- For the positive control 4, add **2 µl** of the supplied positive control 4 (**brown**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Close the capillaries with the lids, transfer them into the LightCycler® carousel and centrifuge them in the LightCycler® centrifuge (if a tabletop centrifuge is used, centrifuge the capillaries in the inserts of the cooling block at 3000 rpm for 15 s). Then transfer the carousel to the LightCycler® and start the PCR programme described in 9.3.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 5.

Table 5: Real-Time-PCR thermal profile

Description	Time	Temperature	Heating rate	Cycles	Acquisition
Initial Denaturation	15 min	95 °C	20 °C/s	1	none
Denaturation	10 s	94 °C	20 °C/s	50	none
Primer annealing	25 s	56 °C	20 °C/s		single
Elongation	25 s	72 °C	20 °C/s		none
Melting curve	30 s	95 °C	20 °C/s	1	none
	30 s	40 °C	20 °C/s	1	none
	0 s	80 °C	0.2 °C/s	1	constant
Cooling	30 s	40 °C	20 °C/s	1	-

10 DATA ANALYSIS

Add an analysis of the type „genotyping“ for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is to be selected at 640 nm (F2) or 705 nm (F3). **The MutaREAL® CC universal (KF29-4-CC) colour compensation kit is required for analysis.**

DPD*2A (F2 / 640 nm)

Temperature mutation allele: 50.0 °C (+/-2 °C)

Temperature wild type allele: 59.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **green curve** - homozygous wild type, **black curve** - heterozygous mutant, **red curve** - homozygous mutant.


Fig. 1: Evaluation of master mix 1 - F2 / 640 nm.

DPD*Haplotyp B3 (F3 / 705 nm)

Temperature mutation allele: 61.0 °C (+/-2 °C)

Temperature wild type allele: 67.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **pink curve** - homozygous wild type, **dark blue curve** - heterozygous mutant, **green curve** - homozygous mutant.


Fig. 2: Evaluation of master mix 1 - F3 / 705 nm.

DPD*D949V (F2 / 640 nm)

Temperature mutation allele: 55.0 °C (+/-2 °C)

Temperature wild type allele: 60.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **green curve** - homozygous wild type, **black curve** - heterozygous mutant, **red curve** - homozygous mutant.


Fig. 3: Evaluation of master mix 2 - F2 / 640 nm.

DPD*13 (F3 / 705 nm)

Temperature mutation allele: 66.5 °C (+/-2 °C)

Temperature wild type allele: 58.5 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **pink curve** - homozygous wild type, **dark blue curve** - heterozygous mutant, **green curve** - homozygous mutant.


Fig. 4: Evaluation of master mix 2 - F3 / 705 nm.

The positive controls 1, 2, 3 & 4 provided contain one DNA each which is heterozygous for the DPD*2A, DPD*Haplotype B3, DPD*D949V or DPD*13 mutation.

11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8 °C for more than 6 months. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. When alleles are detected, the polymorphism tested is indicated. Other rare alleles may be present and are not covered by this method. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme / buffer mix		Contains sufficient for <n> test
	Detection mix 1		Upper limit of temperature

DETECTION MIX 2	Detection mix 2		Manufacturer
CONTROL 1 +	Positive control 1	LOT	Lot number
CONTROL 2 +	Positive control 2		Consult instructions for use
CONTROL 3 +	Positive control 3	CONTENT	Content
CONTROL 4 +	Positive control 4	IVD	<i>In vitro</i> diagnostic medical device
CONTROL -	Negative control		Use by YYYY-MM-DD

14 LITERATURE

1. Del Re et al., EPMA J. 2010 Sep; 1(3): 495–502

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com