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14. LITERATUR

1. Ambrosch, A. et al., 2001. Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic medicine : a journal of the British Diabetic Association*, **18**(3), pp.185–92.
2. Dierkes, J. et al., 2001. Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *Journal of renal nutrition : the official journal of the Council on Renal Nutrition of the National Kidney Foundation*, **11**(2), pp.67–72.
3. Dierkes, J. et al., 2001. Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis*, **158**(1), pp.161–4.
4. Morris, M.C. et al., 2004. Dietary niacin and the risk of incident Alzheimer's disease and of cognitive decline. *Journal of neurology, neurosurgery, and psychiatry*, **75**(8), pp.1093–9.
5. Selhub, J. et al., 1995. Association between plasma homocysteine concentrations and extracranial carotid-artery stenosis. *The New England journal of medicine*, **332**(5), pp.286–91.
6. Basheer, M.P. et al., 2016. Vitamin B complex and homocysteine status and Cognitive impairment in the elderly among Indian population. *Journal of Neuroscience and Behavioral Health*, **8**(4), pp.20–26.
7. Lindqvist, M., Hellstrom, A. & Henriksson, A.E., 2012. Abdominal aortic aneurysm and the association with serum levels of Homocysteine, vitamins B₆, B₁₂ and Folate. *Am J Cardiovasc Dis*, **2**(4), pp.318–322.
8. Loitsch, S.M. et al., 2014. Comparison of a new microbiological assay with a standard high-performance liquid chromatographic method for determination of vitamin B₆ in serum. *Clinical laboratory*, **60**(6), pp.1035–41.

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For use in human and veterinary medicine and in research


Valid from 2024-05-07

REF KIF006

REF KIF006.2

Σ 96

Σ 2x 96

IVD

CE

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Safety information

These accessories are to be used exclusively in accordance with the enclosed instructions for use. Important safety information for this product can be found in chapter 6.

Table of Contents

1. INTENDED PURPOSE	17
2. INTRODUCTION	17
3. PRINCIPLE OF THE TEST	18
4. MATERIAL SUPPLIED	18
5. MATERIAL REQUIRED BUT NOT SUPPLIED	18
6. PRECAUTIONS	19
7. STORAGE AND PREPARATION OF REAGENTS	19
7.1 Water	20
7.2 Preparation of the sterile assay medium	20
7.3 Preparation of the enzyme solution	20
7.4 Preparation of the controls	20
7.5 Preparation of the standard curve	21
7.6 Microtiter plate (PLATE)	21
8. SAMPLE STORAGE AND PREPARATION	22
8.1 Sample pretreatment	22
8.2 Sample dilution	22
9. ASSAY PROCEDURE	22
9.1 Test preparations	22
9.2 Test procedure	22
9.3 Measurement	23
10. EVALUATION OF RESULTS	23
10.1 Calculation	23
10.2 Quality control	24
11. LIMITATIONS	24
12. PERFORMANCE CHARACTERISTICS	24
12.1 Precision and reproducibility	24
12.2 Recovery	25
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
14. REFERENCES	27
15. SYMBOLS	28

1. INTENDED PURPOSE

ID-Vit® Vitamin B₆ is a microtiter plate test kit based on a microbiological method which measures the vitamin B₆ content in serum. The test kit contains the standard and all reagents required to perform the test. An ELISA reader is required for the evaluation of the results. For use in human and veterinary medicine and in research. For *in vitro* diagnostic use only.

2. INTRODUCTION

The vitamin B₆ group comprises three natural forms: pyridoxine, pyridoxamine and pyridoxal, which, during metabolism, are converted to the enzymatically active form pyridoxal-5-phosphate (P5P). Pyridoxal-5-phosphate (here described as „vitamin B₆“) is a cofactor in more than a hundred enzyme reactions. One of its functions is to perform transamination, a key step in both breaking down and producing amino acids in the body. Vitamin B₆ is also necessary for the synthesis of neurotransmitters, as well as hemoglobin in red blood cells. And vitamin B₆ plays a central role in fat metabolism.

Vitamin B₆ deficiency symptoms

- Disturbance of protein biosynthesis
- Muscular weakness, loss of muscle control
- Skin disorders (dermatitis, pigment abnormality etc.)
- Nervous disorders (irritability, depression, palsy)
- Insomnia

Vitamin B₆ deficiency can be considered as a risk factor for myocardial infarction, peripheral vascular diseases and atherosclerosis, especially in connection with the regulation of homocysteine metabolism.

Indications

Vitamin B₆ deficiency can result from:

- Chronic inflammatory bowel diseases (colitis ulcerosa, Crohn´s disease, gluten sensitivity)
- Dialysis
- Liver diseases
- Alcohol abuse
- Homocysteinuria
- Pregnancy and lactation
- taking oral contraceptives or medication for epilepsy or asthma

3. PRINCIPLE OF THE TEST

The serum samples are enzymatically pre-treated and then transferred into the wells of a microtiter plate coated with *Saccharomyces cerevisiae*. The addition of vitamin B₆ in either standards or samples gives a vitamin B₆-dependent growth response until vitamin B₆ is consumed. After incubation at **30°C** for **46–50 h**, the growth of *Saccharomyces cerevisiae* is measured turbidimetrically at 610–630 nm (alternatively at 540–550 nm) in an ELISA reader and compared to a standard curve generated from the dilution series. The amount of vitamin B₆ is directly proportional to the turbidity.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity	
			KIF006	KIF006.2
KIF006/ KIF006.2	DIL	Water	4 x 30 ml	5 x 30 ml
	PLATE	<i>Saccharomyces cerevisiae</i> -precoated microtiter plate	1 x	2 x
	SOL	Sample preparation buffer	7 x 5 ml	13 x 5 ml
	ENZ	Enzyme, lyoph.	7 x	13 x
	ASYMED	Vitamin B ₆ assay medium	4 x	4 x
	STD	Vitamin B ₆ standard, lyoph.	4 x	3 x
	FOL	Adhesive cover foil	1 x whole 3 x half	3 x whole
	FRA	Replacement holder for microtiter strips	1 x	1 x
	CTRL1	Vitamin B ₆ control 1, lyoph.	4 x	3 x
	CTRL2	Vitamin B ₆ control 2, lyoph.	4 x	3 x

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber, 30–37°C
- Water bath (90–100°C)

- Microtiter plate reader 610–630 nm (540–550 nm)
- Calibrated precision pipettors and sterile single-use 20–1 000 µl tips
- 5 ml and 10 ml pipets
- 1.5–2 ml reaction vials
- 0.2 µm sterile polyethersulfone (PES) filter with a disposable syringe
- 15 ml centrifuge tubes, (e.g. Falcon tubes)
- Biocentrifuge (10 000g)
- Vortex

6. PRECAUTIONS

- The test is based on a microbiological method. Contaminations lead to erroneous results.
- Water quality is extremely important for the test. Use only the water delivered with the test kit (**DIL**).
- For sterile filtration, only a sterile polyethersulfone filter must be used.
- It is essential to run a standard curve for each separate assay.
- Measure controls with each assay.
- We recommend measurements in duplicate.
- Do not use reagents beyond the expiration date shown on the label.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- Used microtiter stripes and materials that have been in contact with patient samples must be handled and disposed of as potentially infectious.

7. STORAGE AND PREPARATION OF REAGENTS

- Store test kit and reagents at 2–8 °C.
- Prepare reagents freshly and use them immediately after preparation.
- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. Prepare only the appropriate amount necessary for each run. The kit can be used up to 3 x (KIF006.2) or 4 x (KIF006) within the expiry date stated on the label.

7.1 Water

- Water (**DIL**) for medium (**ASYMED**), standard (**STD**), controls (**CTRL1**, **CTRL2**) and dilutions.
- Push the lid up and pull it back to the rim of the glass, then twist the whole cap off.

7.2 Preparation of the sterile assay medium

- Fresh sterile assay medium has to be prepared each time before performing a test.
- Remove the desiccant bag from the lyophilised assay medium bottle by taking the bag with a forceps and shaking it whilst still inside the bottle. Then remove the clean desiccant bag and discard it.
- Add 10 ml water (**DIL**) to the assay medium bottle (**ASYMED**), close the bottle firmly and vortex well. This amount is sufficient for 6 microtiter stripes.
- Filter the medium using a disposable syringe (10 ml) and the 0.2 µm PES filter into a centrifuge tube (15 ml, e.g. Falcon).
- After this preparation, the sterile assay medium can be used in the test.

7.3 Preparation of the enzyme solution

- Add 4 ml sample preparation buffer (**SOL**) to a vial of lyophilised enzyme (**ENZ**), then homogenise using a vortex.
- Enzyme solution cannot be stored.

7.4 Preparation of the controls

- The lyophilised controls (**CTRL1**, **CTRL2**) have to be resuspended with each **750 µl** water (**DIL**) from the test kit, then homogenise using a vortex.
- After reconstitution, the controls are treated like samples.
- The concentration of the controls changes from lot to lot and is stated in the product specification.

7.5 Preparation of the standard curve

- For the preparation of the standard curve, standard concentrate is needed. To prepare standard concentrate, resuspend the lyophilised standard (**STD**) with x ml water (**DIL**) (x = see quality control protocol) supplied with the test kit, then homogenise using a vortex.
- Prepare a standard curve in 6 sterile reaction tubes (1.5–2 ml volume) from standard concentrate and water (**DIL**) following the scheme depicted in the table below:

Vitamin B ₆ [µg/l]	Water (DIL) [µl]	+	Standard concentrate [µl]	=	Total volume [µl]
Blank: 0	500	+	0	=	500
Standard 1: 0.09	985	+	15	=	1 000
Standard 2: 0.72	440	+	60	=	500
Standard 3: 1.5	375	+	125	=	500
Standard 4: 2.4	300	+	200	=	500
Standard 5: 3.6	200	+	300	=	500

7.6 Microtiter plate (PLATE)

- Store the microtiter plate (**PLATE**) in the aluminium packaging containing the desiccant bag at 2–8 °C.
- The microtiter plate (**PLATE**) has to be protected from humidity and contamination.
- Take care that the aluminium packaging is not damaged.
- Carefully close the aluminium packaging after opening.
- Take only the microtiter stripes needed directly before usage to avoid contamination.

8. SAMPLE STORAGE AND PREPARATION

- Use serum for analysis.
- Samples are stable at 2–8°C for one day in the dark. For longer storage, samples can be frozen and kept at -20°C for 5 months.
- Centrifuge samples prior to measurement (at least 5 min at 10 000 *g*). Use the resulting supernatant in the test.
- Hemolytic samples may give erroneous results and should not be used for analysis. Centrifuge lipemic samples at 13 000 *g* for 10 min before assaying to obtain a serum that is as fat free as possible.

8.1 Sample pretreatment

Add 300 µl serum/control to 300 µl of prepared enzyme solution (ratio 1:2), mix, and incubate at 37°C for 30 min in the dark. Then heat to 95°C for 30 min, cool at 2–8°C for 10 min and centrifuge for 10 min at 10 000 *g*.

8.2 Sample dilution

Take 100 µl from the supernatant of the prepared sample/control, add 400 µl water [DIL] and mix. The sample treatment and dilution result in a total dilution of 1:10 (= sample dilution factor).

9. ASSAY PROCEDURE

9.1 Test preparations

Take as many microtiter strips as needed from kit. Return unused strips and any unused test kit components to the original packaging, and store in the refrigerator. Bring all necessary reagents to room temperature.

9.2 Test procedure

- Take as many microtiter strips as needed from the kit and put them in the second microtiter strip holder (**FRA**).
- Put 150 µl sterile assay medium into each cavity.
- Add 150 µl of the prepared standard dilutions (blank, standard 1–5), samples and controls into the respective cavities. Pre-rinse each pipet tip with standard, control or sample solution, respectively.

- Carefully seal the plate with adhesive cover foil (**FOL**). Important: the cavities must be made airtight by pressing the foil down with the hand!
- Keep at **30 °C** for **46–50 h** in an incubator.

9.3 Measurement

- Press the adhesive cover foil (**FOL**) firmly down again with the hand.
- Turn the microtiter plate (**PLATE**) upside down, place it onto a tabletop and shake the microbes well.
- Turn the microtiter plate (**PLATE**) over again and carefully remove the adhesive cover foil (**FOL**). During this, fix the strips in the frame with your hand because the foil is highly adhesive.
- Remove air bubbles in the cavities using a pipet tip or a needle.
- Read turbidity in an ELISA reader at E 610–630 nm (alternatively at E 540–550 nm).

Please note

- After 46–50 h incubation time, the microtiter plate (**PLATE**) may be stored for a maximum of 48 h in the refrigerator before measuring the turbidity.

10. EVALUATION OF RESULTS

We recommend to use the 4 parameter algorithm to calculate the results. The sample dilution factor has to be considered for data evaluation.

The blank serves as a visual control to exclude contamination and is not taken into account in the calculation. The optical density must be < standard 1. If this is not the case, the analysis must be carried out again.

10.1 Calculation

Vitamin B₆ in µg/l = value from the standard curve × sample dilution factor (10).

Reference value for human serum

Based on studies of serum samples of apparently healthy persons (n = 74), the following values were estimated.

Vitamin B₆: 4.8–17.7 µg/l

Please note

A concentration range of 0.9–36 µg/l vitamin B₆ is covered at a sample dilution of 1:10.

We recommend each laboratory to develop its own normal range as normal ranges strongly depend on the choice of the patient collective. The values mentioned above are only for orientation and can deviate from other published data.

10.2 Quality control

The extinction of the highest standard has to be > 0.6.

Results, generated from the analysis of control samples, should be evaluated for acceptability. The results for the samples may not be valid if within the same assay one or more values of the quality control sample or the highest standard are outside the acceptable limits.

11. LIMITATIONS

Only serum can be used for the test.

12. PERFORMANCE CHARACTERISTICS

The following performance characteristics have been collected using human serum samples.

12.1 Precision and reproducibility**Intraassay (n = 16)**

	Vitamin B₆ [µg/l]	CV [%]
Sample	9.87	9.17

Interassay (n = 5)

	Vitamin B₆ [µg/l]	CV [%]
Sample	10.98	6.78

12.2 Recovery

Samples from 3 patients were spiked with vitamin B₆ and analysed. The mean values are shown below.

Sample (n=5)	Mean value original sample [µg/l]	Spike [µg/l]	Vitamin B ₆ expected [µg/l]	Vitamin B ₆ measured [µg/l]	Recovery Rate [%]
A	11.99	10	21.99	26.00	140
		20	31.99	35.60	118
Recovery rate in total [%]					129

Sample (n=5)	Mean value measured in original sample [µg/l]	Spike [µg/l]	Vitamin B ₆ expected [µg/l]	Vitamin B ₆ measured [µg/l]	Recovery Rate [%]
B	9.46	10	19.46	22.99	135
		20	29.46	32.50	115
Recovery rate in total [%]					125

Sample (n=5)	Mean value measured in original sample [µg/l]	Spike [µg/l]	Vitamin B ₆ expected [µg/l]	Vitamin B ₆ measured [µg/l]	Recovery Rate [%]
C	8.85	10	18.85	18.31	95
		20	28.85	32.36	118
Recovery rate in total [%]					107

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- ID-Vit® is a trademark of Immundiagnostik AG.
- Do not use reagents beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Follow the guidelines for medical laboratories.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be made within 14 days after reception of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Analyse controls with each run.
- Always perform assay according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

14. REFERENCES

1. Ambrosch, A. et al., 2001. Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic medicine : a journal of the British Diabetic Association*, **18**(3), pp.185–92.
2. Dierkes, J. et al., 2001. Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *Journal of renal nutrition : the official journal of the Council on Renal Nutrition of the National Kidney Foundation*, **11**(2), pp.67–72.
3. Dierkes, J. et al., 2001. Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis*, **158**(1), pp.161–4.
4. Morris, M.C. et al., 2004. Dietary niacin and the risk of incident Alzheimer's disease and of cognitive decline. *Journal of neurology, neurosurgery, and psychiatry*, **75**(8), pp.1093–9.
5. Selhub, J. et al., 1995. Association between plasma homocysteine concentrations and extracranial carotid-artery stenosis. *The New England journal of medicine*, **332**(5), pp.286–91.
6. Basheer, M.P. et al., 2016. Vitamin B complex and homocysteine status and Cognitive impairment in the elderly among Indian population. *Journal of Neuroscience and Behavioral Health*, **8**(4), pp.20–26.
7. Lindqvist, M., Hellstrom, A. & Henriksson, A.E., 2012. Abdominal aortic aneurysm and the association with serum levels of Homocysteine, vitamins B6, B12 and Folate. *Am J Cardiovasc Dis*, **2**(4), pp.318–322.
8. Loitsch, S.M. et al., 2014. Comparison of a new microbiological assay with a standard high-performance liquid chromatographic method for determination of vitamin B6 in serum. *Clinical laboratory*, **60**(6), pp.1035–41.

15. SYMBOLS

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com