# **IDKmagnetics® Vitamin D duo LC-MS/MS Extraktionskit / extraction kit**

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## 12. LITERATUR

1. Chapuy MC, Schott AM, Garnero P, Hans D, Delmas PD, Meunier PJ. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *J Clin Endocrinol Metab.* 1996;**81**(3):1129-1133.
2. Scharla SH. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporos Int.* 1998;**8** Suppl 2:S7-S12.
3. Grant WB, Holick MF. Benefits and requirements of vitamin D for optimal health: a review. *Altern Med Rev.* 2005;**10**(2):94-111.
4. Visser M, Deeg DJ, Puts MT, Seidell JC, Lips P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *Am J Clin Nutr.* 2006;**84**(3):616-672.
5. Wicherts IS, van Schoor NM, Boeke AJ, et al. Vitamin D status predicts physical performance and its decline in older persons. *J Clin Endocrinol Metab.* 2007;**92**(6):2058-2065.
6. Wielders JP, Wijnberg FA. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clin Chem.* 2009;**55**(8):1584-1585.
7. Soldin OP, Sharma H, Husted L, Soldin SJ. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clin Biochem.* 2009;**42**(9):823-827.
8. Dirks NF, Ackermans MT, de Jonge R, Heijboer AC. Reference values for 24,25-dihydroxyvitamin D and the 25-hydroxyvitamin D/24,25-dihydroxyvitamin D ratio. *Clin Chem Lab Med.* 2019;**57**(10):e259-e261.

9. Berg AH, Powe CE, Evans MK, *et al.* 24,25-Dihydroxyvitamin d3 and vitamin D status of community-dwelling black and white Americans. *Clin Chem.* 2015;**61**(6):877-884.
10. Cavalier E, Huyghebaert L, Rousselle O, *et al.* Simultaneous measurement of 25(OH)-vitamin D and 24,25(OH)<sub>2</sub>-vitamin D to define cut-offs for CYP24A1 mutation and vitamin D deficiency in a population of 1200 young subjects. *Clin Chem Lab Med.* 2020;**58**(2):197-201.
11. Shieh A, Ma C, Chun RF, *et al.* Associations Between Change in Total and Free 25-hydroxyvitamin D With 24,25-Dihydroxyvitamin D and Parathyroid Hormone. *J Clin Endocrinol Metab.* 2018;**103**(9):3368-3375.
12. Melamed ML, Chonchol M, Gutiérrez OM, *et al.* The Role of Vitamin D in CKD Stages 3 to 4: Report of a Scientific Workshop Sponsored by the National Kidney Foundation. *Am J Kidney Dis.* 2018;**72**(6):834-845.

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
# **IDKmagnetics® Vitamin D duo** **LC-MS/MS extraction kit**

*For the extraction of 24,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>  
and 25-OH-vitamin D<sub>3</sub>/D<sub>2</sub> from plasma and serum*

Valid from 2022-05-25

**KM1321**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>11</b>
<b>2. INTRODUCTION AND CLINICAL RELEVANCE</b>	<b>11</b>
<b>3. MATERIAL SUPPLIED</b>	<b>12</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>13</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>13</b>
<i>Storage</i>	13
<i>Preparation of test reagents</i>	13
<b>6. SAMPLE PREPARATION</b>	<b>13</b>
<b>7. LIMITATIONS</b>	<b>14</b>
<b>8. PRECAUTIONS</b>	<b>14</b>
<b>9. DISPOSAL</b>	<b>14</b>
<b>10. TECHNICAL HINTS</b>	<b>15</b>
<b>11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>15</b>
<b>12. REFERENCES</b>	<b>16</b>

## 1. INTENDED USE

IDKmagnetics® Vitamin D duo LC-MS/MS is a kit for the extraction of the vitamin D metabolites 24,25-dihydroxyvitamin D<sub>3</sub> and 25-hydroxyvitamin D<sub>3</sub>/D<sub>2</sub> from human serum and plasma via immunoaffinity enrichment. The assay is an *in vitro* diagnostic tool for manual and automated use by professional laboratory staff. It is used for the sample preparation for LC-MS based determination of the vitamin D metabolite composition in patients with suspected impaired vitamin D metabolism.

## 2. INTRODUCTION AND CLINICAL RELEVANCE

Vitamin D is either formed in the skin (under the influence of UV light) from precursor molecules or ingested from food. The 25-hydroxylated storage form of vitamin D, is formed in the liver while a second hydroxylation in the kidney can subsequently form the active 1,25-dihydroxyvitamin D or the inactive 24,25-dihydroxyvitamin D.

The enzyme responsible for the formation of the inactive 24,25-dihydroxyvitamin D, 24-hydroxylase, belongs to the cytochrom-P450 superfamily and is encoded by the gene *CYP24A1*. The 24-hydroxylase is regulated in the opposite direction to the 1 $\alpha$ -hydroxylase (*CYP27B1*) and therefore has an indirect regulatory effect on the concentration of the active form 1,25-dihydroxyvitamin D.

The 1 $\alpha$ -hydroxylase normally regulates the concentration of 1,25-dihydroxyvitamin D in the serum very strictly, adapted to the requirements of the metabolism. The concentration of 1,25-dihydroxyvitamin D in the serum will only decrease in the case of a pronounced vitamin D deficiency.

Reasons for an deficiency of 1,25-dihydroxyvitamin D can be metabolic disorders either due to genetic defects of the 1 $\alpha$ -hydroxylase (rare) or due to renal dysfunction (more frequent). Even a slightly impaired kidney function leads to a decrease of the 1,25-dihydroxyvitamin D and 24,25-dihydroxyvitamin D concentration.

Since 1,25-dihydroxy vitamin D has important functions in calcium metabolism as well as supplementing secretion of parathyroid hormone from the parathyroid glands, increasing kidney malfunctioning and an associated deficiency of 1,25-dihydroxy vitamin D leads to development of renal osteopathy, which is characterised by osteomalacia and osteitis fibrosa.

An unphysiological accumulation of 25-hydroxyvitamin D and 1,25-dihydroxyvitamin D occurs, among other things, in case of a mutation of the *CYP24A1* gene since an excess of 25-hydroxyvitamin D is no longer converted to 24,25-dihydroxyvitamin D.

Therefore, 24,25-dihydroxyvitamin D should be quantified beside 25-hydroxyvitamin D and 1,25-dihydroxyvitamin D for an evaluation of the vitamin D status.

The concentration of 24,25-dihydroxyvitamin D in the serum correlates closely with the concentration of 25-hydroxyvitamin D in the serum in case of a sufficient supply with vitamin D. In order to evaluate the response of a vitamin D supplementation with ergocalciferol (vitamin D<sub>2</sub>) or cholecalciferol (vitamin D<sub>3</sub>) it is useful to study the vitamin D metabolite ratio (VMR) of 24,25-dihydroxyvitamin D and 25-hydroxyvitamin D. The VMR remains unaffected by individual variations such as e.g. ethnicity and thus allows a more reliable evaluation of the current 25-hydroxyvitamin D concentration.

### Indications:

- Defect of kidney functions
- Chronic kidney failure
- Haemodialysis following kidney transplantation
- Renal osteopathy
- Osteomalacia from various types of vitamin D metabolism disturbances
- Kidney tubules function disturbances (diabetes insipidus, Fanconi-Syndrom)
- Monitoring of therapy with active vitamin D metabolites
- Ideopathic hypercalciuria
- Hypercalcaemia

## 3. MATERIAL SUPPLIED

Cat. no.	Label	Kit components	Amount
2xKM0003	2xWASHSOL	2xWashing solution	2x80 ml
KM1300	ELUSOL	Elution solution	11 ml
	PURSOL	Purification solution	4 ml
KM1320	M96DWP	96-deep well plate with magnetic beads for extraction of 25-OH-vitamin D <sub>3</sub> /D <sub>2</sub> and 24,25-(OH) <sub>2</sub> -vitamin D <sub>3</sub>	1 piece

For reorders of single components, please use the catalogue number followed by the label without space as product number.

The following accessories for the IDKmagnetics® Vitamin D duo-LC-MS/MS application can be ordered separately at Immundiagnostik AG:

- IDKmagnetics® Vitamin D duo LC-MS/MS kit (KM1320)
- tuning solution for the analytes (KM1300TU)
- tuning solution for the internal standard (KM1300TS)

- UPLC column (KM300SP)
- in-line filter (KM1300IF)
- in-line filter holder (KM1300IH)

Please ask for our single component price list.

#### **4. MATERIAL REQUIRED BUT NOT SUPPLIED**

- Vortex mixer
- 2 ml reaction tubes (e.g. Eppendorf)
- Precision pipettes and disposable tips to deliver 10–1 000  $\mu$ l
- Magnetic separation plate
- Horizontal microtiter plate shaker
- Ethanol, HPLC grade
- Centrifuge

#### **5. STORAGE AND PREPARATION OF REAGENTS**

##### *Storage*

The test reagents should be stored protected from light, dry and their specified storage temperature (2–8 °C). The test reagents stored in this way are usable until the indicated expiry date.

##### *Preparation of test reagents*

The purification solution (PURSOL) is mixed with with 25 % ethanol (v/v) immediately before use (e.g. 2 ml PURSOL with 0.5 ml ethanol). It is recommended to prepare only as much as is needed for the test approach (50  $\mu$ l per sample).

#### **6. SAMPLE PREPARATION**

Serum and plasma can be used as sample. The samples should be stored at -20 °C until testing. Samples are stable for at least four weeks when stored at -20 °C.

Prior to use in the assay, allow all samples and reagents to come to room temperature (18–26 °C).

Mix samples and reagents well before use.

1.	Put 50 µl sample in one well of the M96DWP.
2.	Add 50 µl diluted purification solution (PURSOL diluted with ethanol).
3.	Incubation on the microtiter plate shaker(1 000 rpm*) for 30 min at room temperature.
4.	Remove the supernatant with the help of the magnetic separation plate; keep the beads in the well.
5.	Add 500 µl of wash solution (WASHSOL), shake carefully and discard wash solution with the help of the magnetic separation plate; keep the beads in the well. Carry out this wash step three times in total.
6.	Add 100 µl elution solution (ELUSOL), shake carefully (approx. 1 min).
7.	Transfer approx. 70 µl of the eluate into a new reaction tube with the help of the magnetic separation plate. Note: the magnetic beads should be separated completely since they might harm the LC column.

\*The exact intensity of shaking depends on the microtiter plate shaker and should be adjusted that the magnetic beads stay in suspension during incubation.

## 7. LIMITATIONS

The trueness of the measurement cannot be guaranteed with lipaemic or haemolytic blood samples.

## 8. PRECAUTIONS

- Control samples should be carried along each run.
- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immunodiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.

## 9. DISPOSAL

Washing solution (WASHSOL) purification solution (PURSOL) and elution solution (ELUSOL) must be disposed as non-halogenated solvent.

## 10. TECHNICAL HINTS

- The assay should always be performed according to the enclosed manual.
- Components of the kit package should not be mixed with other kit packages.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Plugs and caps of different reagents should not be swapped.
- The individual components of the kit are designed for a maximum of the specified number of test runs. Any part of the components that has already been used must not be reused, but must be disposed of properly in accordance with local regulations.

## 11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE





- This product has been produced and placed on the market in accordance with the regulation (EU) 2017/746 (IVDR).
- IDKmagnetics® is a brand of Immundiagnostik AG.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- This instruction for use replaces the version dated 2021-11-18.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 12. REFERENCES

1. Chapuy MC, Schott AM, Garnero P, Hans D, Delmas PD, Meunier PJ. Healthy elderly French women living at home have secondary hyperparathyroidism and high bone turnover in winter. EPIDOS Study Group. *J Clin Endocrinol Metab.* 1996;**81**(3):1129-1133.
2. Scharla SH. Prevalence of subclinical vitamin D deficiency in different European countries. *Osteoporos Int.* 1998;**8** Suppl 2:S7-S12.
3. Grant WB, Holick MF. Benefits and requirements of vitamin D for optimal health: a review. *Altern Med Rev.* 2005;**10**(2):94-111.
4. Visser M, Deeg DJ, Puts MT, Seidell JC, Lips P. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *Am J Clin Nutr.* 2006;**84**(3):616-672.
5. Wicherts IS, van Schoor NM, Boeke AJ, et al. Vitamin D status predicts physical performance and its decline in older persons. *J Clin Endocrinol Metab.* 2007;**92**(6):2058-2065.
6. Wielders JP, Wijnberg FA. Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clin Chem.* 2009;**55**(8):1584-1585.
7. Soldin OP, Sharma H, Husted L, Soldin SJ. Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clin Biochem.* 2009;**42**(9):823-827.
8. Dirks NF, Ackermans MT, de Jonge R, Heijboer AC. Reference values for 24,25-dihydroxyvitamin D and the 25-hydroxyvitamin D/24,25-dihydroxyvitamin D ratio. *Clin Chem Lab Med.* 2019;**57**(10):e259-e261.
9. Berg AH, Powe CE, Evans MK, et al. 24,25-Dihydroxyvitamin d3 and vitamin D status of community-dwelling black and white Americans. *Clin Chem.* 2015;**61**(6):877-884.
10. Cavalier E, Huyghebaert L, Rousselle O, et al. Simultaneous measurement of 25(OH)-vitamin D and 24,25(OH)2-vitamin D to define cut-offs for CYP24A1 mutation and vitamin D deficiency in a population of 1200 young subjects. *Clin Chem Lab Med.* 2020;**58**(2):197-201.
11. Shieh A, Ma C, Chun RF, et al. Associations Between Change in Total and Free 25-hydroxyvitamin D With 24,25-Dihydroxyvitamin D and Parathyroid Hormone. *J Clin Endocrinol Metab.* 2018;**103**(9):3368-3375.

12. Melamed ML, Chonchol M, Gutiérrez OM, et al. The Role of Vitamin D in CKD Stages 3 to 4: Report of a Scientific Workshop Sponsored by the National Kidney Foundation. *Am J Kidney Dis.* 2018;**72**(6):834-845.

**Used symbols:**

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Unique Device Identification

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)