IDK[®] Cortisol & Cortisone LC-MS/MS Kit

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
IDK[®] Cortisol & Cortisone LC-MS/MS Kit

*For the in vitro determination of
cortisol & cortisone in urine*

Valid from 2023-10-27

REF KM2110

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	16
2. INTRODUCTION	16
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION, STORAGE AND STABILITY OF REAGENTS	19
<i>Calibrators and controls</i>	19
<i>Dilution solution</i>	19
<i>Internal standard</i>	20
<i>Mobile phases A and B</i>	20
<i>Autosampler wash solution</i>	20
6. PREPARATION OF SAMPLE	21
7. LC-MS/MS METHOD	21
8. EXAMPLES OF CHROMATOGRAMS	22
9. QUALITY CONTROL	23
<i>Reference range</i>	23
10. TEST CHARACTERISTICS	23
<i>Repeatability - intra-day precision</i>	23
<i>Reproducibility - inter-day precision</i>	23
<i>Measuring range with limit of blank (LOB)</i>	23
11. PRECAUTIONS	24
12. DISPOSAL	24
13. TECHNICAL HINTS	24
14. GENERAL NOTES ON THE TEST	25
15. REFERENCES	25
16. SYMBOL EXPLANATION	26

1. INTENDED USE

IDK® Cortisol & Cortisone LC-MS/MS assay is for the quantitative detection of cortisone and cortisol levels in 24-hour human urine samples. The assay is an *in vitro* diagnostic medical device for automated use by professional laboratory staff in a laboratory environment. The assay serves as an aid in the diagnosis and monitoring of functional disorders of the pituitary and adrenal glands (e.g. Cushing's syndrome, Addison disease, congenital adrenal hyperplasia or adrenal insufficiency). It may also serve as an aid in the assessment and monitoring of chronic stress and in the comprehensive understanding of the complex context of obesity and metabolic syndrome.

2. INTRODUCTION

Urinary free cortisol (UFC) analysis represents the first biochemical laboratory approach for the screening of endogenous Cushing's syndrome (CS). Endogenous CS is caused by prolonged exposure to elevated levels of endogenous cortisol that may occur from excess production by one or both adrenal glands, or from overproduction of the adrenocorticotropic hormone (ACTH), which normally regulates cortisol production. As symptoms are always non-specific, including hypertension, truncal obesity and mood disorders, specific biochemical tests are required for diagnosing CS.

One of the first-line tests for this diagnosis, the measurement of 24-h UFC, can also be useful in other clinical conditions characterized by a high urine cortisol level, such as in non-autonomous hypercortisolism (pseudo-CS), psychiatric disorders, morbid obesity, poorly controlled diabetes mellitus and alcoholism. Altered cortisol metabolism is also responsible for a condition called apparent mineralocorticoid excess (AME) syndrome.

While cortisol is mainly essentially secreted by the adrenal gland, cortisone is mainly produced by type 2 11 β -hydroxysteroid dehydrogenase (11 β -HSD) in multiple human tissues, which interconverts bioactive cortisol to hormonally inactive cortisone to prevent activation of the mineralocorticoid receptor by cortisol. Hence the simultaneous determination of cortisol and cortisone can help in the diagnosis of AME syndrome, but also in congenital adrenal hyperplasia and adrenal insufficiency.¹


Fig. 1: Cortisol.


Fig. 2: Cortisone.

3. MATERIAL SUPPLIED

Cat. No.	Identifier	Kit components	Quantity
KM2110	AUTOWASH	Autosampler wash solution	1 x 1 000 ml
	CAL1–6	Calibrators 1–6; lyophilised (see product specification for concentration)	2 vials (à 1 ml) per level
	CTRL1–3	Controls 1–3; lyophilised (see product specification for concentration)	3 vials (à 1 ml) per level
	DILSOL	Dilution solution	3 x 100 ml
	INTSTD	Internal standard; lyophilised	3 x 6 ml
	MOPHAA	Mobile phase A	1 x 500 ml
	MOPHAB	Mobile phase B	1 x 500 ml

For reorders of single components, please use the catalogue number followed by the identifier without space as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

The following accessories are required for the IDK® Cortisol & Cortisone LC-MS/MS application (not included in the kit):

- Ultrapure water*
- Precision pipettors and disposable tips to deliver 10–1 000 µl
- Centrifuge 10 000 *g* (at least)
- Vortex mixer or microtiter plate shaker
- Standard laboratory disposable plastic reagent tubes (e.g. 1.5 ml)
- LC-MS/MS system and LC-MS vials

*Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696/LC-MS grade), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

The following accessories for the IDK® Cortisol & Cortisone LC-MS/MS application can be ordered separately at Immundiagnostik AG:

- UPLC column (KM2110SP)
- all single components

Please ask for our single component price list. Please contact us for customized inquiries.

5. PREPARATION, STORAGE AND STABILITY OF REAGENTS

Note: Please unpack the kit components from the transport packaging immediately upon receipt and follow the instructions for storage conditions printed on the product labels. In general, all components should be stored protected from light, dry and at the specified storage temperature.

All components are for LC-MS/MS use only, components may also contain other ingredients than those listed as active ingredients below which might influence the measurement. All declared stabilities are only valid in case of no bacterial contamination.

Calibrators and controls

Handling:

Always remove the cap and rubber plug carefully (in order to avoid loss of content).

Reconstitute the calibrators and controls as follows:

- Reconstitute each calibrator and control with exactly 1 ml distilled or deionised water and incubate for 15 min at room temperature.
- Next, mix the component thoroughly to make sure that all dry material has dissolved; do not shake too vigorously to avoid foam formation.
- Handle the prepared component as a patient sample during the test procedure.

Stability and storage:

Before reconstitution: -20°C Until expiry date printed on the product label.

After reconstitution: -20°C 2 weeks

 2–8°C 1 week

Dilution solution

Handling:

The component is liquid and ready for use.

Stability and storage:

Store at 2–8°C

After first opening the component can be used for up to 6 weeks if closed and stored at 2–8°C or 2 weeks on the LC.

Internal standard

Handling:

Always remove the cap and rubber plug carefully (in order to avoid loss of content).

Reconstitute the internal standard as follows:

- Reconstitute the content of one vial INTSTD with exactly 6 ml DILSOL and incubate for 15 min at room temperature. Put aside the used vial of DILSOL with the remaining 94 ml.
- Next, mix the component thoroughly to make sure that all dry material has dissolved; do not shake too vigorously to avoid foam formation.
- Finally, return the 6 ml reconstituted INTSTD into the DILSOL vial (with the remaining 94 ml) and mix thoroughly.

Stability and storage:

Before reconstitution: 2–8 °C Until expiry date printed on the product label.

After reconstitution: -20 °C 2 weeks

2–8 °C 1 week

Mobile phases A and B

Handling:

The component is liquid and ready for use.

Stability and storage:

Store at 2–8 °C

After first opening the component can be used for up to 6 weeks if closed and stored at 2–8 °C or 2 weeks on the LC.

Autosampler wash solution

Handling:

The components are liquid and ready for use.

Stability and storage:

Store at 2–8 °C

After first opening the component can be used for up to 6 weeks if closed and stored at 2–8 °C or 2 weeks on the LC.

6. PREPARATION OF SAMPLE

Use human 24-hour urine.

Stability and storage:

Store at 2–8 °C 14 days

Store at -20 °C 2 months

Sample preparation

1.	Centrifuge sample, CAL or CTRL for 5 min at 10 000 <i>g</i> .
2.	Transfer 100 µl of the supernatant (see step 1.) into one reaction tube or well of a 96-well plate..
3.	Add 900 µl DILSOL (with added INTSTD, see "Preparation, storage and stability of reagents"). Mix well (avoid foam formation).
4.	Transfer the solution into a suitable LC-MS vial or 96-well plate for the LC-MS autosampler.
5.	Injection in the LC-MS/MS (see application note).

7. LC-MS/MS METHOD

Please refer to the application note or contact lcms@immundiagnostik.com for the LC-MS/MS method parameters.

8. EXAMPLES OF CHROMATOGRAMS


Fig. 3: Chromatogram (recorded with Waters TQS-micro LC-MS/MS) for the internal standard (top) and analyte (bottom) of cortisol.


Fig. 4: Chromatogram (recorded with Waters TQS-micro LC-MS/MS) for the internal standard (top) and analyte (bottom) of cortisone.

9. QUALITY CONTROL

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. Immundiagnostik AG cannot guarantee the trueness of the method if one or more values of the quality control samples are outside the acceptable range (see product specification).

Reference range

Cortisol: 5–133 nmol/l (1.8–48.2 ng/ml)

Cortisone: 74–328 nmol/l (26.7–118.2 ng/ml)

We recommend each laboratory to establish its own reference range.

10. TESTCHARACTERISTICS

Repeatability - intra-day precision

Analyte	CV		
	Control 1 [%]	Control 3 [%]	Sample [%]
Cortisol	2.0	1.0	1.4
Cortisone	1.7	1.2	0.9

Reproducibility - inter-day precision

Analyte	CV		
	Control 1 [%]	Control 3 [%]	Sample [%]
Cortisol	3.0	2.5	3.2
Cortisone	2.9	2.6	3.4

Measuring range with limit of blank (LOB)

Analyte	nmol/l (ng/ml)
Cortisol	0–1 200 (0–435)
Cortisone	0.218–3 100 (0.1–1 117)

11. PRECAUTIONS

- Human material used in the kit components was tested and found to be negative for HBsAg, HIV-1, HIV-2 and HIV p24 Ag antibodies, Anti-HTLV 1 and 2, HCV and HIV genome. However, for safety reasons, all kit components should be treated as potentially infectious.
- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immundiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.

12. DISPOSAL

Autosampler wash solution (AUTOWASH), mobile phase A (MOPHAA), mobile phase B (MOPHAB) and dilution solution (DILSOL) must be disposed as non-halogenated solvents.

Calibrators (CAL1–6) and controls (CTRL1–3) should be disposed as potentially infectious material in accordance with local regulations.

13. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according to the enclosed manual.
- Plugs and caps of different reagents should not be swapped.
- The individual components of the kit are designed for a maximum of the specified number of test runs. Any part of the components that has already been used must not be reused, but must be disposed properly in accordance with local regulations.


14. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- IDK® is a trademark of Immundiagnostik AG.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be sent to Immundiagnostik AG together with a written complaint.

15. REFERENCES

1. Antonelli, G. et al. (2014) Cortisol and cortisone ratio in urine: LC-MS/MS method validation and preliminary clinical application. *Clinical Chemistry and Laboratory Medicine*, 52(2), 213-220.

16. SYMBOL EXPLANATION

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com