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

**Manual**  
*For professional use only*

# 3-OH-Isovalerylcarnitine LC-MS/MS kit

*For the determination of 3-OH-isovalerylcarnitine  
in citrate plasma*

Gültig ab/Valid from 2023-05-15

**REF** KM3200

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com) [www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED PURPOSE</b>	<b>16</b>
<b>2. MATERIAL PROVIDED</b>	<b>16</b>
<b>3. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>4. REAGENTS: STORAGE AND PREPARATION</b>	<b>17</b>
<i>Storage</i>	17
<i>Preparation of mobile phases</i>	17
<i>Preparation of calibrators and controls</i>	18
<b>5. SAMPLES: STORAGE AND PREPARATION</b>	<b>18</b>
<i>Sampling</i>	18
<i>Sample storage</i>	18
<i>Sample preparation</i>	18
<b>6. TEST PROCEDURE</b>	<b>18</b>
<i>Test principle</i>	18
<i>Processing</i>	19
<i>Supporting information</i>	19
<b>7. RESULT EVALUATION</b>	<b>20</b>
<i>Quality control</i>	20
<i>Calculation</i>	20
<i>Metrological traceability of the calibrators and controls</i>	20
<i>Example of chromatograms</i>	20
<b>8. ANALYTICAL AND CLINICAL PERFORMANCE</b>	<b>22</b>
<i>Accuracy and precision</i>	22
<i>Measuring range with quantification limits (LLOQ &amp; ULOQ)</i>	22
<i>Matrix effect and interferences</i>	23
<i>Carry over</i>	23
<i>Recovery</i>	23
<i>Cut off value</i>	23
<b>9. WARNINGS AND PRECAUTIONS</b>	<b>23</b>
<b>10. DISPOSAL</b>	<b>24</b>
<b>11. SYMBOLS</b>	<b>25</b>

**SAFETY INFORMATION**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit. Important safety information for this product can be found in chapter 9.

## 1. INTENDED PURPOSE

The 3-OH-Isovalerylcarnitine LC-MS/MS kit is an LC-MS/MS assay for the quantification of 3-OH-isovalerylcarnitine by LC-MS/MS in citrate plasma after precipitation reaction. This can be used for the diagnosis of functional biotin deficiency. The assay is an in vitro diagnostic medical device for manual or automatic use by professional laboratory staff in a laboratory environment.

## 2. MATERIAL PROVIDED

Cat. no.	Identifier	Kit components	Amount
KM0001	ACTSOL	Activation solution	1 x 1,5 ml
KM0002	RECSOL	Reconstitution solution	1 x 15 ml
KM3200	96DWP	96-deep well plate	2 x 1 piece
	CAL1–6	Calibrators 1–6; lyophilised (see product specification for concentration)	2 vials (à 250 µl) per level each
	CTRL1–3	Control 1–3; lyophilised (see product specification for concentration)	2 vials (à 250 µl) per level each
	DILSOL	Dilution solution	1 x 11 ml
	MOPHAA	Mobile phase A	1 x 500 ml
	MOPHAB	Mobile phase B	1 x 500 ml
	PREC	Precipitation reagent (contains the internal standard)	2 x 17 ml

For reorders of single components, please use the catalogue number followed by the identifier without space as ordering number.

## 3. MATERIAL REQUIRED BUT NOT SUPPLIED

The following material is required for the 3-OH-Isovalerylcarnitine LC-MS/MS application (not included in the kit):

- Precision pipettors and disposable tips to deliver 10–1 000 µl
- Centrifuge, 150 g for microtiter plates
- Vortexer
- microtiter plate mixer
- LC-MS/MS system

- LC-MS vials or 96-well plate suitable for the autosampler
- for the procedure in reaction tubes: 1.5 ml reaction tubes

The following accessories for the 3-OH-Isovalerylcarnitine LC-MS/MS kit can be ordered separately at Immundiagnostik AG:

- Tuning solution for the analyte (KM3200TU)
- Tuning solution for the internal standard (KM3200TS)
- UPLC-column (KM3200SP)
- In-line filter (KM3200IF)
- In-line filter holder (KM3200IH)

Please ask for our single component price list and contact us for your individual solution.

## 4. REAGENTS: STORAGE AND PREPARATION

### *Storage*

The test reagents should be stored protected from light, dry and at their specified storage temperature (CAL1–6, CTRL1–3, PREC: -20°C; all others 2–8°C). The test reagents stored in this way are usable until the indicated expiry date.

**Note:** After preparation of the reagents their stability may vary (see respective preparation step).

### *Preparation of mobile phases*

Before use, the mobile phases (MOPHAA and MOPHAB) must be activated by adding activation solution (ACTSOL) according to the following chart:

Component			ACTSOL [µl]
Name	[ml]		
Mobile phase A (MOPHAA)	500	+	500
Mobile phase B (MOPHAB)	500		500

Mobile phases should be degassed prior to use.

**Note:** After activation with activation solution (ACTSOL), the components mobile phase A (MOPHAA) and mobile phase B (MOPHAB) can be stored at 2–8°C up to 2 weeks. It is therefore recommended to prepare only as much as is needed for the test approach.

**Attention:** The activation solution (ACTSOL) must be added under the fume hood. All used vials must be absolutely clean, detergent-free and preferably made of glass or plastic suitable for LC-MS/MS.

### *Preparation of calibrators and controls*

Calibrators (CAL1–6) and controls (CTRL1–3) have to be reconstituted with 250 µl reconstitution solution (RECSOL) each. Mix the vial content thoroughly for 30 s by using a vortex mixer.

**Note:** The reconstituted CAL1–6 and CTRL1–3 are stable for 20 days at 2–8 °C.

## **5. SAMPLES: STORAGE AND PREPARATION**

### *Sampling*

Citrate plasma samples are to be used with this assay.

The blood sampling must be performed with appropriate blood collection tubes according to the manufacturer's instructions and by qualified medical personnel.

### *Sample storage*

The samples can be stored for 4 months at -20 °C.

Repeated freeze and thaw cycles should be avoided.

### *Sample preparation*

The samples must be brought to room temperature (18–26 °C) prior the test procedure.

## **6. TEST PROCEDURE**

### *Test principle*

For the extraction of 3-OH-isovalerylcarnitine, a precipitation reaction is performed. After chromatographic separation, the analyte is measured via tandem mass spectrometry (MS/MS) in multiple reaction monitoring (MRM) mode.

A stable isotope-labelled version of the analyte is used as internal standard to correct sample loss during the sample preparation and matrix effects.

## Processing

Only reagents and samples that are at room temperature (18–26°C) should be used in the test (exception: PREC).

Vortex samples and reagents before use.

	96-deep well plate	1.5 ml reaction tube
1.	Add 50 µl sample, calibrator (CAL) or control (CTRL) into a well of the 96DWP respectively a reaction tube.	
2.	Add 150 µl ice-cold (-20°C) PREC.	
3.	Seal 96DWP and shake for 5 min at 1 000 rpm.	Seal tubes and vortex for 1 min.
4.	Centrifuge for 15 min at 150 g.	Centrifuge for 15 min at 8 900 g.
5.	Add 50 µl DILSOL into one well of a new 96 well plate* or LC-MS vial (analogous to the number of samples, CAL and CTRL used).	
6.	Transfer 50 µl of the supernatant (see step 4.) to the DILSOL presented.	
7.	Mix the sample by pipetting up and down.	-
8.	Seal plate and shake for 5 min at 1 000 rpm.	Seal LC-MS vials and vortex for 5 s.
9.	Injection into the LC-MS system (see application note).	

\*suitable for the autosampler of the LC-MS system.

## Supporting information

Please refer to the application note or contact [lcms@immundiagnostik.com](mailto:lcms@immundiagnostik.com) for the LC-MS/MS method parameters.

The exact concentrations of the controls (CTRL1–3) and calibratos (CAL1–6) may vary from lot to lot and can be found in the enclosed product specification.

## 7. RESULT EVALUATION

### *Quality control*

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results may not be valid if one or more values of the quality control sample are outside the acceptable limits (see product specification) within the same assay.

We furthermore recommend to monitor the intensity of the internal standard within one batch and, in the long term, the ion transition ratio in order to identify deviations in performance.

### *Calculation*

Linear regression is used as model for evaluation of the results. The six calibrator concentration points are connected by a straight line (weighting factor:  $1/x^2$ ). This is used to calculate the concentration of the samples and controls via the peak area.

When using citrate plasma, the dilution factor in the blood collection tube with the citrate solution (normally factor 1:10) has to be noted for the evaluation.

The results of the assay must be interpreted by medical professionals with experience in the intended field of application and taking into account reference values as well as patient-specific factors.

### *Metrological traceability of the calibrators and controls*

The calibrators (CAL1–6) and controls (1–3) are based on synthetic pure substance and every lot is compared to its previous lot as well as external calibration in multiple experiments.

### *Example of chromatograms*

It must be noted that the retention time and signal intensity may vary depending on the device. The chromatograms shown here were recorded by using a Waters TQ-XS LC-MS/MS.

## Blank


Fig. 1: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 3-OH-isovalerylcarnitine for a blank sample.

## Control


Fig. 2: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 3-OH-isovalerylcarnitine for a control (CTRL1).

## Sample


Fig. 3: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 3-OH-isovalerylcarnitine for a patient sample (citrate plasma).

## 8. ANALYTICAL AND CLINICAL PERFORMANCE

### Accuracy and precision

Sample [ng/ml]	Accuracy		Precision	
	<i>intra-day</i> (n=5)	<i>inter-day</i> (n=40)	<i>intra-day</i> (n=5)	<i>inter-day</i> (n=40)
0.47	96.1 – 114.7 %	101.9 %	≤ 11.4 %	9.3 %
3.43	100.3 – 109.8 %	104.2 %	≤ 9.4 %	5.2 %
10.72	99.6 – 104.5 %	101.2 %	≤ 9.6 %	5.4 %
59.49	102.6 – 109.2 %	105.1 %	≤ 8.3 %	4.8 %

### Measuring range with quantification limits (LLOQ & ULOQ)

The assay has been verified as linear in a range from 0.47–88.82 ng/ml ( $R^2 \geq 0.99$ ; 95 % CI of the slope is between 0.9 and 1.1; 95 % CI of the y-intercept brackets the origin).

### *Matrix effect and interferences*

Citrate plasma exhibits an ion suppressing effect which can be compensated by the internal standard. No interference was detected for the internal standard in citrate plasma.

### *Carry over*

Carry over can never be completely excluded due to the analytical method. However, this was verified as not significant for the analyte and the internal standard (< 15% deviation at the LLOQ after preceding injection of an ULOQ sample).

### *Recovery*

The recovery was verified to be stable ( $96.1 \pm 2.8\%$  and  $89.8 \pm 2.7\%$ ) across different concentration ranges (4.0 ng/ml; 10.0 ng/ml; 60.0 ng/ml) for the analyte and internal standard, respectively.

### *Cut off value*

Cut off value\* in citrate plasma: 12.5 ng/ml

Higher values of 3-OH-isovalerylcarnitine can be taken to represent a functional biotin deficiency.

\*97.5th percentile of healthy blood donors: N = 100, 40% female, aged 20–60.

## **9. WARNINGS AND PRECAUTIONS**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only and in accordance with the intended use.
- The assay should always be performed according to the enclosed manual. Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Human material used in the kit components was tested and found to be negative for Anti-HIV, Anti-HCV, HBsAg, Anti-HBc, Anti-*Treponema pall.*, HIV-1-RNA and HCV-RNA. However, for safety reasons, all kit components should be treated as potentially infectious.

- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immundiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- Any part of the components that has already been used must not be reused, but must be disposed in accordance with local regulations.
- Do not interchange different lot numbers of any kit component within the same assay. Do not mix plugs and caps from different reagents.
- Reagents should not be used beyond the expiry date stated on the kit label.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.


## 10. DISPOSAL

Mobiles phases (MOPHAA, MOPHAB) and precipitation reagent (PREC) have to be disposed as non-halogenated solvents.

Dilution solution (DILSOL) has to be disposed as halogenated solvent.

Calibrators (CAL1–6) and controls (CTRL1–3) should be disposed as potentially infectious material in accordance with local regulations.

## 11. SYMBOLS

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)