

IDK[®] Casomorphin/Gliadorphin- Peptide LC-MS/MS

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
IDK[®] Casomorphin/Gliadorphin peptides LC-MS/MS

For the determination in stabilised urine

Valid from 2021-10-04

REF **KM8000**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. INTRODUCTION	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. STORAGE AND PREPARATION OF REAGENTS	19
<i>Storage</i>	19
<i>Preparation of mobile phases and test reagents</i>	19
<i>Preparation of the calibrators and controls</i>	19
6. STORAGE OF SAMPLES	20
7. ASSAY PROCEDURE	20
<i>Preparation of the cartridges</i>	20
<i>Test procedure</i>	20
8. LC-MS/MS METHOD	21
9. CALCULATION	21
10. CHROMATOGRAM EXAMPLES	21
11. QUALITY CONTROL	24
<i>Reference range</i>	24
12. PERFORMANCE CHARACTERISTICS	25
<i>Precision - reproducibility and repeatability</i>	25
<i>Sensitivity</i>	26
13. PRECAUTIONS	27
14. DISPOSAL	27
15. TECHNICAL HINTS	27
16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	27

1. INTENDED USE

The *IDK*[®] Casomorphin/Gliadorphin peptides LC-MS/MS kit is an *in vitro* diagnostic tool for the quantitative determination of casomorphin/gliadorphin peptides in stabilised urine by LC-MS/MS after solid phase extraction. The kit is intended for manual use by professional laboratory staff.

2. INTRODUCTION

Exorphins are bioactive peptides which can react with opiate receptors. These peptides are derived from e.g. the incomplete hydrolysis of the milk protein casein and the wheat protein gluten and gliadin. Casomorphin (BCM7 from casein) and gliadorphin-7 (from gluten) have been widely studied and are linked to gastrointestinal, neurological, and neuro-developmental disorders.

These short peptides are usually broken down by the enzyme dipeptidyl-peptidase IV (DPPIV). In people with autism, schizophrenia and other neurological disorders, either a genetic deficiency or inhibition of this enzyme, e.g. by drugs treating diabetes type 2, mercury or yeast, can occur.

Both ways, this proteolytic enzyme can't act properly and leaves incompletely digested peptides behind which are able to pass the blood-brain barrier after the uptake from the intestinal lining into the bloodstream. In the brain, these bioactive peptides can bind to opiate receptors and mimic the effects of opiate drugs like heroine and morphine. This can result in change of behavior, lack of focus and attention, sleepiness and even aggression and self-abuse.

Wheat and milk allergies are associated casomorphin/gliadorphin problems. Even if no milk or wheat allergies occur some people can still react negatively to these peptides.

Children with autism and ADHS can benefit from an elimination diet – gluten-free, casein-free (GfCf) diet. Physical and mental health can be improved.

This test for casomorphin/gliadorphin peptides is superior to any other tests because no nutritional regimen, e.g. abandonment of soy products, prior to sample collection is required.

3. MATERIAL SUPPLIED

Art. no.	Label	Kit components	Amount
2 x KM0001	2 x ACTSOL	2 x Activation solution	2 x 1.5 ml
4 x KM0002	4 x RECSOL	4 x Reconstitution solution	4 x 15 ml
6 x KM0003	6 x WASHSOL	6 x Wash solution	6 x 80 ml
KM8000	CAL1–2	Calibrators 1 und 2, lyophilised (for concentration see product specification)	5 vials (á 2.2 ml) per level each
	CTRL1–2	Controls 1 und 2, lyophilised (for concentration see product specification)	5 vials (á 2.2 ml) per level each
	ELUSOL	Elution solution	2 x 27 ml
	INTSTD	Internal standard	1.5 ml
	MOPHAA	Mobile phase A	1 000 ml
	MOPHAB	Mobile phase B	1 000 ml
	SOLA	Solution A	22 ml

For reorders of single components, please use the catalogue number followed by the label without space as product number.

The following accessories for the *IDK*[®] Casomorphin/Gliadorphin peptides LC-MS/MS kit can be ordered separately at Immundiagnostik AG:

- tuning solution for the analytes (KM8000TU)
- tuning solution for the internal standard (KM8000TS)
- UPLC column (KM8000SP)
- in-line filter (KM8000IF)
- in-line filter holder (KM8000IH)
- SPE cartridge (KM8000CK)

Please ask for our single component price list.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Standard laboratory glass vials, suitable for LC-MS/MS
- Calibrated precision pipettors and 10–1 000 µl tips
- Solid phase extraction unit
- Evaporation unit (e.g. vacuum centrifuge or nitrogen distributor)
- Vortex
- LC-MS/MS equipment
- Borosilicate tubes

- Methanol p.a.
- Ultra pure water*

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩcm).

5. STORAGE AND PREPARATION OF REAGENTS

Storage

The test reagents should be stored protected from light, dry and their specified storage temperature (CAL1–2, CTRL1–2: -20 °C; all others 2–8 °C). The test reagents stored in this way are usable until the indicated expiry date.

Note: After preparation of the test reagents for the test procedure other stabilities might apply (see respective preparation step).

Preparation of mobile phases and test reagents

Before use, the mobile phases (MOPHAA and MOPHAB) and solution A (SOLA) must be activated by adding activation solution (ACTSOL) according to the following chart:

Component			ACTSOL [µl]
Name	[ml]		
Mobile phase A (MOPHAA)	500	+	500
Mobile phase B (MOPHAB)	500		500
Solution A (SOLA)	25		25

Prior use mobile phases should be degassed.

Note: After activation with activation solution (ACTSOL), the components mobile phase A (MOPHAA), mobile phase B (MOPHAB), sample buffer (SAMPLEBUF), wash solution 2 (WASHSOL2), dilution solution (DILSOL) and elution solution (ELUSOL) can be stored at 2–8 °C up to 2 weeks. It is therefore recommended to prepare only as much as is needed for the test approach.

Attention: The activation solution (ACTSOL) must be added under the fume cupboard. All vessels to be used must be absolutely clean, free of detergents and preferably made of LC-MS/MS suitable glass.

Preparation of the calibrators and controls

The calibrators (CAL1–2) and the controls (CTRL1–2) are dissolved in 2.2 ml reconstitution solution (RECSOL) each while 30 s vortexing.

Note: After reconstitution with the reconstitution solution (RECSOL), the calibrators (CAL1–2) and the controls (CTRL1–2) are not stable and can not be stored.

6. STORAGE OF SAMPLES

In unstabilised urine, the concentration of casomorphin/gliadorphin peptides declines quickly, therefore only stabilised urine is suited as sample. Immundiagnostik offers urine sample stabilising tubes (order no DZ9030UT) specifically designed for this purpose. Samples stabilised using those tubes are stable for 48h at room temperature. For longer storage, the sample must be frozen at $< -17^{\circ}\text{C}$.

7. ASSAY PROCEDURE

Control samples should be analysed with each run.

Prior to use in the assay, allow all samples and reagents to come to room temperature ($18\text{--}26^{\circ}\text{C}$).

Mix samples and reagents well before use.

Preparation of the cartridges

Conditioning of the cartridges with 1 ml methanol, followed by equilibration with 1 ml ultra pure water.

Test procedure

1.	Mix 1 ml reconstituted CAL/reconstituted CTRL/sample and 10 μl INTSTD in a reaction tube.
2.	Add solution onto the SPE cartridges and aspirate.
3.	Wash with 1 ml wash solution (WASHSOL), aspirate. Carry out this wash step four times in total.
4.	Eluate with 500 μl elution solution (ELUSOL) into borosilicate tubes.
5.	Dry the eluate, e.g. in a vacuum centrifuge.
6.	Reconstitute the dried sample with 200 μl activated solution A (SOLA).
7.	Injection into the LC-MS/MS system (see applicatio note).

8. LC-MS/MS METHOD

Please refer to the application note or contact lcms@immundiagnostik.com for the parameters for setting the LC-MS/MS method.

9. CALCULATION


The linear regression is used as model for evaluation of the results. The two calibrator concentration points are connected by a straight line. The samples can be calculated using the obtained line.

10. CHROMATOGRAM EXAMPLES


β -Casomorphin 7 (molecular weight: 789.95 g/mol)

Sample Name: Casoglia150616_001

F7:MRM of 2 channels,ES+
790.43 > 136.1
5.926e+004


F7:MRM of 2 channels,ES+
790.43 > 70.05
2.175e+005


INTSTD (β -Casomorphin 7 heavy labelled, molecular weight: 795.95 g/mol)

Sample Name: Casoglia150616_001

F8:MRM of 2 channels,ES+
796.63 > 74.95
1.517e+005


F8:MRM of 2 channels,ES+
796.63 > 70.36
2.215e+004


Casomorphin 1-3 (molecular weight: 425.50 g/mol)**Sample Name: Casoglia150616_001**


F1:MRM of 2 channels,ES+
426.82 > 136.1
9.308e+003


F1:MRM of 2 channels,ES+
426.82 > 70.32
6.050e+003

**Casomorphin 1-4** (molecular weight: 522.62 g/mol)**Sample Name: Casoglia150616_001**


F3:MRM of 3 channels,ES+
523.13 > 70.08
1.668e+005


F3:MRM of 3 channels,ES+
523.13 > 116.07
4.734e+004







**Casomorphin 1-4 amide** (molecular weight: 521.63 g/mol)**Sample Name: Casoglia150616_001**



F2:MRM of 2 channels,ES+
522.27 > 136.07
1.409e+005

F2:MRM of 2 channels,ES+
522.27 > 70.05
2.280e+005

Exorphin A5 (molecular weight: 599.65 g/mol)**Sample Name: Casoglia150616_001****Exorphin B5** (molecular weight: 594.67 g/mol)**Sample Name: Casoglia150616_001****Exorphin C** (molecular weight: 591.72 g/mol)**Sample Name: Casoglia150616_001**

Gliadorphin (molecular weight: 876.00 g/mol)**Sample Name: Casoglia150616_001****11. QUALITY CONTROL**

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on Immundiagnostik studies of samples of apparently healthy persons (n = 18), the following mean values were estimated

Analyte	ng/ μ mol Creatinine
Casomorphin 1–3	1.406
Exorphin A5	2.033
β -Casomorphin 7	0.138
Gliadorphin	0.187
Casomorphin 1–4 amide	0.234
Casomorphin 1–4	0.188
Exorphin C	0.076
Exorphin B5	0.611

We recommend each laboratory to establish its own reference range.

12. PERFORMANCE CHARACTERISTICS

Precision - reproducibility and repeatability

Inter-Assay (n = 10)

Analyte	Sample	ng/ml	CV [%]
Casomorphin 1–3	1	1.05	10.79
	2	4.55	2.41
Exorphin A	1	3.75	24.35
	2	11.37	9.27
β-Casomorphin 7	1	1.56	12.32
	2	6.93	4.66
Gliadorphin	1	1.96	7.38
	2	7.80	7.14
Casomorphin 1–4 amide	1	1.44	11.40
	2	5.57	6.07
Casomorphin 1–4	1	2.09	5.13
	2	8.71	5.54
Exorphin C	1	2.84	5.16
	2	11.76	8.53
Exorphin B	1	2.88	5.37
	2	10.91	4.31

Intra-Assay (n = 4)

Analyte	Sample	ng/ml	CV [%]
Casomorphin 1–3	1	3.44	12.92
	2	1.14	3.63
Exorphin A	1	6.29	26.98
	2	3.05	39.92
β-Casomorphin 7	1	3.19	8.42
	2	1.52	11.88

Analyte	Sample	ng/ml	CV [%]
Gliadorphin	1	3.17	8.03
	2	2.14	7.59
Casomorphin 1-4 amide	1	2.33	17.08
	2	1.34	15.58
Casomorphin 1-4	1	2.83	7.79
	2	2.06	7.77
Exorphin C	1	3.91	4.11
	2	3.05	4.69
Exorphin B	1	5.24	13.38
	2	2.72	8.84

Sensitivity

The detection limit (LLOD) designates the lowest concentration of the analyte that can still be detected.

Analyte	pg/ml Urine
Casomorphin 1-3	1 132
Exorphin A	1 872
β -Casomorphin 7	420
Gliadorphin	520
Casomorphin 1-4 amide	624
Casomorphin 1-4	307
Exorphin C	473
Exorphin B	976

It should be noted that the detection limit depends not only on the application method but also on the instrument.

13. PRECAUTIONS

- Control samples should be analysed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immunodiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.

14. DISPOSAL

Mobile phase B (MOPHAB), elution solution (ELUSOL), solution A (SOLA) and activation solution (ACTSOL) must be disposed as non-halogenated solvents. The calibrators (CAL1–2) and controls (CTRL1–2) should be disposed due to their treatment as potentially infectious material in accordance with local regulations.

15. TECHNICAL HINTS


- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according to the enclosed manual.
- Plugs and caps of different reagents should not be swapped.

16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- *IDK*® is a trademark of Immunodiagnostik AG.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com