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

[www.immundiagnostik.com](http://www.immundiagnostik.com)

**Safety information**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit. Important safety information for this product can be found in chapter 9.

## Table of Contents

<b>1. INTENDED PURPOSE</b>	<b>12</b>
<b>2. MATERIALS PROVIDED</b>	<b>12</b>
<b>3. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>13</b>
<b>4. PREPARATION, STORAGE AND STABILITY OF REAGENTS</b>	<b>13</b>
<i>Preparation of mobile phases</i>	13
<i>Preparation of calibrators and controls</i>	14
<i>Preparation of enzyme</i>	14
<i>Preparation of internal standard</i>	14
<b>5. SAMPLE PREPARATION</b>	<b>14</b>
<b>6. LC-MS/MS METHOD</b>	<b>16</b>
<b>7. EVALUATION</b>	<b>16</b>
<b>8. QUALITY CONTROL</b>	<b>16</b>
<b>9. PRECAUTIONS</b>	<b>16</b>
<b>10. DISPOSAL</b>	<b>16</b>
<b>11. TECHNICAL HINTS</b>	<b>17</b>
<b>12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>17</b>
<b>13. SYMBOL EXPLANATION</b>	<b>18</b>

## 1. INTENDED PURPOSE

This LC-MS/MS application is intended for the quantitative determination of human  $\alpha_4\beta_7$  integrin therapeutic antibody vedolizumab (e.g. Entyvio®) in serum. This kit is designed for laboratory professional and research use only. Not for use in diagnostic procedures.

## 2. MATERIALS PROVIDED

Cat. No.	Label	Kit components	Quantity
KMR0001	ACTSOL	Activation solution	1 x 1.5 ml
KMR0002	RECSOL	Reconstitution solution	1 x 15 ml
KMR9000	COLUMNS	ImmuTube®-columns for extraction of vedolizumab	50 pieces
	ELUREAG	Elution reagent	1 x 3 ml
	ENZ	Enzyme	5 x 21 µl
	REABUF	Reaction buffer	1 x 2 ml
	STOP	Stop solution	1 x 1 ml
	WASHBUF	Wash buffer	1 x 55 ml
	WASHSOL	Wash solution	1 x 30 ml
KMR9600	CAL1–6	Calibrators 1–6; lyophilised (see product specification for concentration)	1 vial (à 200 µl) per level
	CTRL1–3	Controls 1–3; lyophilised (see product specification for concentration)	1 vial (à 200 µl) per level
	INTSTD	Internal standard	1 x 11 µl
	MOPHAA	Mobile phase A	1 x 250 ml
	MOPHAB	Mobile phase B	1 x 250 ml

For reorders of single components, please use the catalogue number followed by the label without space as product number.

The following accessories for the IDKmonitor® Vedolizumab drug level LC-MS/MS extraction kit can be ordered separately at Immundiagnostik AG:

- tuning solution for the analyte (KMR9600TU)
- tuning solution for the internal standard (KMR9600TS)

- UPLC column (KMR9600SP)
- in-line filter (KMR9600IF)
- in-line filter holder (KMR9600IH)

Please ask for our single component price list. Feel free to contact us for customized inquiries.

### 3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Calibrated precision pipettes and 10–1 000 µl tips
- 1.5 ml and 2 ml reaction tubes (Immundiagnostik AG recommends Protein LoBind Tubes from Eppendorf)
- Multistepper
- Centrifuge, 1 000 g
- Vortex
- Standard single-use laboratory plastic vials, cups, etc.
- Thermal block, heatable to at least 70 °C
- Overhead rotator for reaction tubes
- LC-MS/MS equipment

### 4. PREPARATION, STORAGE AND STABILITY OF REAGENTS

To run the assay multiple times, ensure that the reagents are stored at the conditions stated on the label. Prepare only the amount required for each run.

Reagents with a volume less than 100 µl should be centrifuged before use to avoid loss of volume.

#### *Storage*

The test reagents should be stored protected from light, dry and at their specified storage temperature (INTSTD, CAL1–6, CTRL1–3, ENZ: -20 °C; all others: 2–8 °C). The test reagents stored in this way are usable until the specified expiry date.

**Attention:** The matrix of the COLUMNS must not run dry.

**Note:** After preparation of the reagents their stability may vary (see respective preparation step).

#### *Preparation of mobile phases*

MOPHAA and MOPHAB have to be activated with 0.5 ml ACTSOL before use and sonicated in the ultrasonic bath for 15 min.

**Attention:** If possible, use glass pipettes to pipet ACTSOL. The ACTSOL must be added under the fume hood. All used vials must be absolutely clean, detergent-free and preferably made of glass or plastic suitable for LC-MS/MS.

**Note:** Mobile phases spiked with activation reagent can be stored up to 2 weeks in a closed flask at 2–8°C.

### *Preparation of calibrators and controls*

CAL1–6 and CTRL1–3 have to be reconstituted with 200 µl RECSOL each. Mix the vial content thoroughly and allow it to dissolve for 10 minutes to ensure complete reconstitution.

Reconstituted CAL1–6 and CTRL1–3 have to be diluted 1:10 with WASHSOL prior to use (e.g. 10 µl CAL and 90 µl WASHSOL). 50 µl diluted sample will be needed per sample preparation.

**Note:** The reconstituted CAL1–6 and CTRL1–3 are stable up to 10 days at 4°C.

### *Preparation of enzyme*

ENZ should be thawed and diluted with 84 µl WASHSOL just before the assay is performed (see step 8.). One aliquot (à 21 µl) is sufficient for 10 samples.

**Note:** The diluted ENZ is not stable and cannot be stored.

### *Preparation of internal standard*

The INTSTD is briefly centrifuged and diluted with 1 089 µl WASHBUF in the reaction tube of the INTSTD.

**Note:** The diluted INTSTD is stable for 7 days at 4°C.

## 5. SAMPLE PREPARATION

This assay is intended for serum samples only.

The samples must be diluted 1:10 with WASHSOL before use (e.g. 10 µl sample and 90 µl WASHSOL). For the test procedure 50 µl per diluted sample are required.

Quality control samples should be analysed with each run.

Reagents and samples should be used at room temperature (18–26 °C) - Exception: the enzyme (ENZ). Vortex samples, solutions and reagents before use.

- |    |   |
|----|---|
| 1. | Place the COLUMNS into an appropriate rack. Assure that the included suspension does not stick to the lid. Label the COLUMNS. |
|----|---|

2.	<p>Unscrew the lids of the COLUMNS and open the outlet by breaking the endpiece at the predetermined breaking point. Do <u>not</u> discard the endpiece. Put the COLUMNS into 2 ml reaction tubes.</p> <p>Add 150 µl WASHBUF and centrifuge for 30s at 1000g; discard the flowthrough. Carry out this wash step (2.) two times in total.</p>
3.	<p>Add 100 µl WASHSOL into the COLUMNS. Close the COLUMNS firmly by recapping the outlet with the endpiece in reversed orientation.</p> <p>Add 20 µl diluted INTSTD.</p> <p>Add 50 µl diluted CAL, diluted CTRL or diluted sample. Close the COLUMNS by screwing the lid.</p>
4.	<p>Mix for 1 h by slow end-over-end rotation at room temperature. After the incubation, shortly centrifuge the columns to collect matrix remains in the lid.</p>
5.	<p>Unscrew the lid of the COLUMNS and open the outlet. Put the COLUMNS into plastic reaction tubes and centrifuge for 30s at 1000g.</p> <p>Add 150 µl WASHBUF and centrifuge for 30s at 1000g; discard the flowthrough. Carry out this wash step with WASHBUF (5.) two times in total.</p>
6.	<p>Add 150 µl WASHSOL and centrifuge for 30s at 1000g; discard the flowthrough. The COLUMNS should not stand dry for too long.</p> <p>Carry out this wash step with WASHSOL (6.) two times in total.</p>
7.	<p>Place COLUMNS in already labelled new 1.5 ml reaction tubes.</p> <p>Add 40 µl ELUREAG to the matrix, centrifuge for 30–60s at 1000g until the matrix is completely dry. Discard the COLUMNS.</p>
8.	<p>Add 20 µl REABUF to the eluate.</p>
9.	<p>Thaw the required amount of ENZ aliquots and dilute with 84 µl WASHSOL. Add 10 µl diluted ENZ to the reaction tubes.</p>
10.	<p>Incubation for 90 min at 70 °C and 1 000 rpm agitating in a thermoblock.</p>
11.	<p>Allow reaction tubes to cool down to room temperature and centrifuge briefly before the next step.</p>
12.	<p>Add 10 µl STOP.</p>
13.	<p>Transfer to an LC-MS vial/96-well plate suitable for the autosampler and inject into the LC-MS/MS system (see application note).</p>

## 6. LC-MS/MS METHOD

Please refer to the application note or contact [lcms@immundiagnostik.com](mailto:lcms@immundiagnostik.com) for the LC-MS/MS method parameters.

## 7. EVALUATION

Linear regression is used as model for the evaluation of the results. The six calibrator concentration points are connected by a straight line (weighting factor:  $1/x$ ). The peak area is used to determine the concentration of the samples and controls.

## 8. QUALITY CONTROL

Control samples should be analysed with each run. The results of the control samples are used to confirm the accuracy of the method. The test results may not be valid if one or more values of the quality control sample are outside the acceptable limits (see product specification) within the same assay.

## 9. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. Still, all kit components should be treated as potentially infectious.
- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immundiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.
- The test components contain organic solvents. Avoid contact with skin or mucous membranes.

## 10. DISPOSAL

- MOPHAA, MOPHAB and ACTSOL have to be disposed as non-halogenated solvents, ELUREAG and STOP as halogenated solvents. CAL1–6 and CTRL1–3 should be disposed as potentially infectious material in accordance with local regulations.


## 11. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date stated on the kit label.
- The assay should always be performed according to the enclosed manual.
- Do not mix plugs and caps from different reagents.
- The individual components of the kit are designed for the specified number of test runs. Any part of the components that has already been used must not be reused.

## 12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the kit package are for research use only. Not for use in diagnostic procedures.
- The guidelines for laboratories should be followed.
- *IDKmonitor*® and *ImmuTube*® are trademarks of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints regarding deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

### 13. SYMBOL EXPLANATION

	Temperature limitation		Catalogue number
	For research use only		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)