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
*For professional use only*

# Coenzyme Q<sub>10</sub> UHPLC Kit

*For the determination of coenzyme Q<sub>10</sub> (ubiquinone) in  
EDTA-whole blood, serum and EDTA-plasma*

Valid from 2024-11-26

**REF** KU1700

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## **Safety information**

The assay has to be performed exclusively according to the instructions for use enclosed with the kit. Important safety information for this product can be found in the chapter WARNINGS AND PRECAUTIONS.

## Table of Contents

1. INTENDED USE _____	11
2. SUMMARY AND EXPLANATION OF THE TEST _____	11
3. PRINCIPLE OF THE TEST _____	12
4. MATERIAL SUPPLIED _____	12
5. MATERIAL REQUIRED BUT NOT SUPPLIED _____	13
6. STORAGE AND PREPARATION OF REAGENTS _____	13
7. SPECIMEN COLLECTION AND PREPARATION _____	14
8. WARNINGS AND PRECAUTIONS _____	14
9. GENERAL NOTES ON THE TEST AND TEST PROCEDURE _____	14
10. SYMBOL EXPLANATION _____	15

## 1. INTENDED USE

This HPLC application is intended for the quantitative determination of coenzyme Q<sub>10</sub> (ubiquinone) in EDTA-whole blood, serum and EDTA-plasma. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Ubiquinone was first isolated in the 1950s from the group of Prof. Green (Wisconsin). The function was investigated by Prof. Mitchel, who received the Nobel price for his research on the oxydative phosphorylation pathway.

Ubiquinone is a coenzyme which is represented in every cell and the whole metabolism. It is built up by a chinonic ring and an isoprenic side chain. In humans, ubiquinone can be synthesized and taken up by nutrition.

Ubiquinone has two different pysiological functions:

- 1) Component of the energy metabolism
- 2) Radical scavanger

During the reduction of oxygen in the oxydative phosphorylation 3 mol ATP are generated. In this reduction electrons are transferred from NADPH to oxygen via 6 different redox systems. Ubiquinone is the least concentrated and thus speed-determining in the redox system in the mitochondrial membranes. It therefore has a control function in the energy metabolism and, in the case of a deficit on the entire metabolism. With normal nutrition enough amounts are taken up and the human organism is in the able to provide coenzyme Q<sub>10</sub> itself. However, with increasing age and under exposure to sunlight the content of coenzyme Q<sub>10</sub> is reduced up to 50%.

Ubiquinone has a high amount of carbon doublebonds and therefore a higher potential of reduction than vitamin C or vitamin E. It will therefore be consumed first when a radical attack occurs and thus represents a protection for particularly sensitive systems. Q<sub>10</sub> reacts very quickly and sensitively to free radicals, even if other protective systems such as vitamin E still behave inertly. Q<sub>10</sub> thus becomes an optimal membrane stabilizer, which stabilizes the ion channels and mantains optimal metabolism.

### Indications

- Determination of ubiquinone status
- Cardiovascular disease
- Carcinogenesis
- Aging
- Fatigue symptoms

### 3. PRINCIPLE OF THE TEST

The first step in determining coenzyme Q<sub>10</sub> is a sample preparation. The internal standard, the precipitation reagent and the extraction solution are added to the calibrator, controls and samples. The supernatant is then completely dried. The residue is resuspended with dilution solution and then injected into the UHPLC system.

The separation via UHPLC follows an isocratic method using a reversed phase column. One run lasts 1.5 minutes. The chromatograms are recorded by a UV detector. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas/heights by the internal standard method.

#### Summary

The kit includes all reagents for preparation and separation of the samples, except the column.

As for many other parameters, the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical chemical routines quickly and precisely. Mostly, a one-point calibration is sufficient for calibrating the test system – unlike immunoassays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. With short test series, the one-point calibration is much more economic than 6-point calibration for immunoassays.

### 4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 0005.15	RECSOL	Reconstitution solution	15 ml
KU1700	CAL	Calibrator; lyophilised	4 x
	CTRL1	Control 1; lyophilised	4 x
	CTRL2	Control 2; lyophilised	4 x
	INTSTD	Internal Standard	1 ml
	EXTSOL	Extraction solution	110 ml
	PREC	Precipitation reagent	40 ml
	SOL	Dilution solution	25 ml
	MOPHA	Mobile phase	1 000 ml

The HPLC column (KU1700RP), can be ordered separately from Immundiagnostik. In addition to the complete kits, all components can also be ordered separately. Please ask for our single component price list.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Glass reaction tubes (plastic-coated vessels can lead to interference signals)
- Centrifuge
- Various pipettes
- UHPLC with UV detector
- Evaporation unit
- Reversed phase C<sub>18</sub> column
- Vortex

## 6. STORAGE AND PREPARATION OF REAGENTS

- The lyophilised calibrator (**CAL**) is stable at **-20 °C** until the expiry date stated on the label. Before use, it has to be reconstituted in **x µl** (x = see the enclosed product specification for the volume needed) reconstitution solution (**RECSOL**). The concentration of coenzyme Q<sub>10</sub> slightly changes from lot to lot, please refer to the enclosed product specification sheet for the exact amount.
- The lyophilised controls (**CTRL1** and **CTRL2**) are stable at **-20 °C** until the expiry date stated on the label. Before use, they have to be reconstituted each with **x µl** (x = see the enclosed product specification for the volume needed) reconstitution solution (**RECSOL**). The concentration of coenzyme Q<sub>10</sub> slightly changes from lot to lot, the exact concentration is stated on the specification data sheet.
- The internal standard (**INTSTD**) is stable at **-20 °C** until the expiry date stated on the label.
- Before using the assay for the first time, add the complete content of the internal standard (**INTSTD**) to the precipitation reagent vial (**PREC**) and mix. The mixture must be stored at **-20 °C**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 7. SPECIMEN COLLECTION AND PREPARATION

EDTA-whole blood, serum or EDTA-plasma are suited for this test system.

The samples are stable for 1 day at **2–8 °C** in the dark. For longer storage, samples should be frozen at **-20 °C**.

## 8. WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components must be treated as potentially infectious.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.

## 9. GENERAL NOTES ON THE TEST AND TEST PROCEDURE


- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Do not use reagents beyond the expiration date shown on kit label.
- Do not mix plugs and caps from different reagents.
- Do not interchange different lot numbers of any kit component within the same assay.
- Always perform the assay according to the enclosed manual.
- Follow the guidelines for medical laboratories.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not in agreement with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be made within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

## 10. SYMBOL EXPLANATION

### Symbol explanation:

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)