



Data Sheet

ADALIMUMAB

ANTIBODY, RECOMBINANT MONOCLONAL, COUPLED TO HSA

Catalog no.:	A 1234.3 / A 1234.4
Immunogen:	Adalimumab
Synonyms:	Humira®
DrugBank:	DB00051
Host:	Human cell line
Clone no.:	AF311-B4-HSA
Form:	Monovalent single chain variable fragment (scFv) coupled to HSA
Matrix:	PBS pH 7.4
Specificity:	Adalimumab variable domain
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA (capture)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at:	2-8°C
Last updated on:	31 October 2024

For research use only

Publishing research using A 1234? Please let us know so that we can cite your publication as a reference.

