



Data Sheet

HUMAN ADIPONECTIN



ANTIBODY, MONOCLONAL

Catalog no.:	A 1060.1 / A 1060.2
Immunogen:	Recombinant human Adiponectin
Synonyms:	Adipocyte, C1q and collagen domain-containing protein; 30 kDa adipocyte complement-related protein; Adipocyte complement-related 30 kDa protein (ACRP30); Adipose most abundant gene transcript 1 protein (apM-1); Gelatin-binding protein
Swiss-Prot No:	Q15848
Gene Information:	Gene Name: ADIPOQ, ACDC, ACRP30, APM1, GBP28 Gene ID: 9370
Host:	Mouse Balb/c
Clone no.:	M14-6 (A)-H11
Isotype:	IgG
Matrix:	Cell culture supernatant, Protein G purified, PBS pH 7.4
Specificity:	Human Adiponectin
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only



Antibodies



Publishing research using A 1060? Please let us know so that we can cite your publication as a reference.



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 6251 70190-0 · Fax: +49 6251 70190-363 · dept.immuochemicals@immundiagnostik.com · www.immundiagnostik.com