



Data Sheet

HUMAN α 1-ANTITRYPSIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 1021.1 / A 1021.2
Immunogen:	Human α 1-Antitrypsin
Synonyms:	Alpha-1 protease inhibitor, Alpha-1-antiproteinase, Serpin A1
Swiss-Prot No:	P01009
Gene Information:	Gene Name: SERPINA1, AAT, PI GeneID: 5265
Host:	Rabbit
Matrix:	Serum
Specificity:	Enteral/hepatic α 1-Antitrypsin
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1021? Please let us know so that we can cite your publication as a reference.

