



Data Sheet

HUMAN α 1-ANTITRYPSIN

ANTIBODY, MONOCLONAL

Catalog no.:	A 2022.1 / A 2022.2
Immunogen:	Native α 1-Antitrypsin-Trypsinogen II-complex purified from human pancreas
Synonyms:	Alpha-1 protease inhibitor, Alpha-1-antiproteinase, Serpin A1
Swiss-Prot No:	P01009
Gene Information:	Gene Name: SERPINA1, AAT, PI GenelD: 5265
Host:	Mouse Balb/c
Isotype:	IgG ₁
Matrix:	50 mM Na-citrate pH 6.0, 0.9% NaCl, 0.1% NaN ₃
Specificity:	Human hepatic α 1-Antitrypsin
Contents:	10 μ g / 100 μ g (lyophilized) Resuspend in 10 μ l / 100 μ l aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	18 March 2022

For research use only

Publishing research using A 2022? Please let us know so that we can cite your publication as a reference.

