

Data Sheet

ANGIOTENSIN II-BIOTIN SYNTHETIC

Catalog no.:	A 1035AG.3
Name:	Angiotensin II linked to biotin
Source:	Synthetic
Purity:	> 90%
Contents:	10 µg (1 mmol/l)
Store at:	- 20 °C
Known applications:	ELISA (coating antigen 10 ng/ml)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.

Last updated on: 04 May 2022

For research use only

Publishing research using A 1035AG? Please let us know so that we can cite your publication as a reference.

