



## Data Sheet

# CHICKEN IgY

## ANTIBODY, POLYCLONAL

|                            |   |
|----------------------------|---|
| <b>Catalog no.:</b>        | G 4018  |
| <b>Immunogen:</b>          | Purified IgY  |
| <b>Host:</b>               | Goat  |
| <b>Matrix:</b>             | IgG, Protein G purified from serum, then affinity purified using chicken IgY;<br>10 mM PBS, 0.15 M NaCl pH 7.2, 0.05% NaN <sub>3</sub>  |
| <b>Specificity:</b>        | The antibody reacts with heavy chain on chicken IgY and light chain on all chicken immunoglobulins.   |
| <b>Contents:</b>           | 1 mg (10.1 mg/ml)   |
| <b>Known applications:</b> | ELISA<br><br>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| <b>Store at:</b>           | 2-8° C  |
| <b>Last updated on:</b>    | 4 May 2022  |

### For research use only

**Publishing research using G 4018? Please let us know so that we can cite your publication as a reference.**

