



## Data Sheet

# CHICKEN IgY

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	G 4014
<b>Immunogen:</b>	Purified IgY
<b>Host:</b>	Donkey
<b>Matrix:</b>	Serum
<b>Specificity:</b>	The antibody reacts specifically with IgY
<b>Contents:</b>	10 ml
<b>Known applications:</b>	RIA (second antibody, 1:32)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

<b>Store at:</b>	- 20 °C
	Repeated thawing and freezing must be avoided

<b>References:</b>	1. Karmatschek M, Maier I, Seibel MJ, Woitge HW, Ziegler R, Armbruster FP (1997). Improved purification of human bone sialoprotein and development of a homologous radioimmunoassay. <i>Clin Chem</i> <b>43</b> (11): 2076-2082.
--------------------	--

<b>Last updated on:</b>	4 May 2022
-------------------------	------------

### For research use only

**Publishing research using G 4014? Please let us know so that we can cite your publication as a reference.**

