



Data Sheet

RAT ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 1-28)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | AE1022.1 / AE1022.2 |
| Immunogen: | Synthetic rat ANP (aa 1-28) poly Lysin conjugated (SLRRSSCFGGRIDRIGAQSGLGCNSFRY) |
| Synonyms: | Atrial natriuretic factor, ANF, Natriuretic peptides A (aa 1-28), Prepronatriodilatin (aa 1-28) |
| Swiss-Prot No: | P01161 |
| Gene Information: | Gene Name: Nppa GenelD: 24602 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Synthetic rat ANP (aa 1-28), rat ANP |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | RIA (1:3000), immunocytochemistry (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 20 April 2022 |

For research use only

Publishing research using AE1022? Please let us know so that we can cite your publication as a reference.

