



Data Sheet

HUMAN ATRIAL NATRIURETIC PEPTIDE

[ANP] (aa 1-28)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1020.1 / AE1020.2
Immunogen:	Synthetic human ANP (aa 1-28) poly Lysin conjugated (SLRRSSCFGGRMDRIGASGLGNSFRY)
Synonyms:	Atrial natriuretic factor, ANF
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA GenelD: 4878
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human ANP (aa 1-28)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1020? Please let us know so that we can cite your publication as a reference.

