



Data Sheet

HUMAN PRO-ATRIAL NATRIURETIC PEPTIDE [PRO-ANP] (aa 40-55)

ANTIBODY, MONOCLONAL

Catalog no.:	AE 6206.1 / AE 6206.2
Immunogen:	Synthetic human pro-ANP (aa 40-55) poly-Lysine conjugated ((C)EPNEEAGAALSPLPEV)
Synonyms:	Atrial Natriuretic Factor, ANF
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA GeneID: 4878
Host:	Mouse Balb/c
Clone no.:	A137-M6-19-B6
Isotype:	IgG _{G3 kappa}
Matrix:	Cell culture supernatant, Prot. G purified, PBS pH 7.4
Specificity:	Synthetic human pro-ANP (aa 40-55)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE6206? Please let us know so that we can cite your publication as a reference.

