



Data Sheet

HUMAN PRO-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 1-30) SYNTHETIC

Catalog no.:	AE1026AG.1 / AE1026AG.2
Name:	Human pro-ANP (aa 1-30) (NPMYNAVSNADLMDFKNLLDHLEEKMPLED)
Synonyms:	Cardiodilatin-related peptide
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GeneID: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Reference:	1. Hartter E, Khalafpour S, Missbichler A, Hawa G, Woloszczuk W (2000). Enzyme immunoassays for fragments (epitopes) of human proatrial natriuretic peptides. <i>Clin Chem Lab Med</i> 38 (1): 27-32.
Last updated on:	12 May 2022

For research use only

Publishing research using AE1026AG? Please let us know so that we can cite your publication as a reference.

