

Data Sheet

HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [proANP] (aa 79-98) SYNTHETIC

Catalog no.:	AE1025AG.1 / AE1025AG.2
Name:	Human pro-ANP (aa 79-98) (SSDRSALLKSKLRALLTAPR)
Synonyms:	Atrial natriuretic factor (ANF)
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GenelD: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Reference:	1. Missbichler A, Hawa G, Schmal N, Woloszczuk W (2001). Sandwich ELISA for proANP 1-98 facilitates investigation of left ventricular dysfunction. <i>Eur J Med Res</i> 6(3): 105-111.
Last updated on:	12 May 2022

For research use only

Publishing research using AE1025AG? Please let us know so that we can cite your publication as a reference.

