



## Data Sheet

# HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 31-67)

## SYNTHETIC

<b>Catalog no.:</b>	AE1024AG.1 / AE1024AG.2
<b>Name:</b>	Human pro-ANP (aa 31-67) (EVVPPQVLSEPNEEAGAALSPLPEVPPWTGEVSPAQR)
<b>Synonyms:</b>	Atrial natriuretic factor (ANF)
<b>Swiss-Prot No:</b>	P01160
<b>Gene Information:</b>	Gene Name: NPPA, ANP, PND GeneID: 4878
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	12 May 2022

**For research use only**

**Publishing research using AE1024AG? Please let us know so that we can cite your publication as a reference.**

