



## Data Sheet

# HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 40-55) SYNTHETIC

<b>Catalog No.:</b>	AE6202AG.1 / AE6202AG.2
<b>Name:</b>	Human pro-ANP (aa 40-55) ( <u>C</u> SDRSALLKSKLRAL) <small>Underlined amino acids differ from the original sequence; they were added or substituted for technical reasons.</small>
<b>Synonyms:</b>	Atrial natriuretic factor (ANF)
<b>Swiss-Prot No:</b>	P01160
<b>Gene Information:</b>	Gene Name: NPPA, ANP, PND GenelD: 4878
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	02 April 2025

**For research use only**

**Publishing research using AE6203AG? Please let us know so that we can cite your publication as a reference.**

