



Data Sheet

RAT ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 1-28) SYNTHETIC

Catalog No.:	AE1022AG.1 / AE1022AG.2
Name:	Rat ANP (aa 1-28) (SLRRSSCFGGRIDRIGASGLGCNSFRY)
Synonyms:	Atrial natriuretic factor (ANF)
Swiss-Prot No:	P01161
Gene Information:	Gene Name: Nppa GeneID: 24602
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	25 September 2023

For research use only

Publishing research using AE1022AG? Please let us know so that we can cite your publication as a reference.

