



Data Sheet

HUMAN β -DEFENSIN 6 (aa 1-38)

ANTIBODY, POLYCLONAL

Catalog no.:	AE 1042.1 / AE 1042.2
Immunogen:	Synthetic human β -Defensin 6 (aa 1-38) (FFDEKCNKLGKGTCKNNCGKNEELIALCQKSLKCCRTIQ)
Synonyms:	Beta-defensin 106, DEFB6
Swiss-Prot No:	Q8N104
Gene Information:	Gene Name: DEFB106A, BD6, DEFB106, DEFB6 and DEFB106B GeneID: 245909 and 503841
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human β -Defensin 6 (aa 1-38)
Contents:	20 μ / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1042? Please let us know so that we can cite your publication as a reference.

