



Data Sheet

HUMAN β -DEFENSIN 2 (aa 4-41)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1081.3 / A 1081.4
Immunogen:	Synthetic human β -Defensin 2 (aa 4-41) (DPVTCLKSGAICHVPFCPRRYKQIGTCGLPGTKCKKPK)
Synonyms:	Beta-Defensin 2, beta 2, BD-2, Skin-antimicrobial peptide, SAP1
Swiss-Prot No:	O15263
Gene Information:	Gene Name: DEFB4 GenelD: 1673
Host:	Sheep
Matrix:	IgG, Protein G purified from serum; 50 mM TRIS pH 7.4
Specificity:	Synthetic human β -Defensin 2 (aa 4-41)
Contents:	10 μ g / 100 μ g (lyophilized) Resuspend in 10 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1081? Please let us know so that we can cite your publication as a reference.

