



## Data Sheet

# HUMAN $\beta$ -DEFENSIN 5 (aa 17-51)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE 1040.1 / AE 1040.2
<b>Immunogen:</b>	Synthetic human $\beta$ -Defensin 5 (aa 17-51) (ESCKLGRGKCRKECLENEKPDGNCRLNFLCCRQRI)
<b>Synonyms:</b>	Beta-defensin 105, DEFB-5
<b>Swiss-Prot No:</b>	Q8NG35
<b>Gene Information:</b>	Gene Name: DEFB105A, BD5, DEFB105, DEFB5 and DEFB105B GeneID: 245908 and 504180
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human $\beta$ -Defensin 5 (aa 17-51)
<b>Contents:</b>	20 $\mu$ l / 100 $\mu$ l (lyophilized) Resuspend in 20 $\mu$ l / 100 $\mu$ l aqua bidest.
<b>Known applications:</b>	ELISA (1:2000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	20 April 2022

### For research use only

**Publishing research using AE1040? Please let us know so that we can cite your publication as a reference.**

