



## Data Sheet

# HUMAN $\beta$ -DEFENSIN 1 (aa 18-26)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE1049.1 / AE1049.2
<b>Immunogen:</b>	Synthetic human b-Defensin 1 (aa 18-26) (PIFTKIQGT)
<b>Synonyms:</b>	Beta-Defensin 1 (BD-1), Defensin beta 1
<b>Swiss-Prot No:</b>	P60022
<b>Gene Information:</b>	Gene Name: DEFB1 Gene ID: 1672
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human b-Defensin 1 (aa 18-26)
<b>Contents:</b>	20 $\mu$ l / 100 $\mu$ l (lyophilized) Resuspend in 20 $\mu$ l / 100 $\mu$ l aqua bidest.
<b>Known applications:</b>	Western Blot  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	20 April 2022

### For research use only

**Publishing research using AE1049? Please let us know so that we can cite your publication as a reference.**

