



Data Sheet

HUMAN β -DEFENSIN 1 (aa 28-34)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1048.1 / AE1048.2
Immunogen:	Synthetic human β -Defensin 1 (aa 28-34) (YRGKAKC(ACM)G)
Synonyms:	β -Defensin 1 (BD-1), Defensin β 1
Swiss-Prot No:	P60022
Gene Information:	Gene Name: DEFB1 Gene ID: 1672
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human β -Defensin 1 (aa 28-34)
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA, Western Blot This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1048? Please let us know so that we can cite your publication as a reference.

