



Data Sheet

HUMAN BONE SIALOPROTEIN [BSP] (aa 262-279)

ANTIBODY, POLYCLONAL

Catalog no.:	A 4218.1 / A 4218.2
Immunogen:	Synthetic human Bone Sialoprotein (aa 262-279) (YESENGEPRGDNYRAYED)
Synonyms:	Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein
Swiss-Prot No:	P21815
Gene Information:	Gene Name: IBSP Gene ID: 3381
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Bone Sialoprotein (aa 262-279) Cross reactivity with human Bone Sialoprotein. There was no cross reactivity obtained with human Osteonectin and human Osteopontin (aa 25-40).
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:1000), immunocytochemistry (1:100) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	06 April 2022

For research use only

Publishing research using A 4218? Please let us know so that we can cite your publication as a reference.

