



## Data Sheet

# HUMAN BONE SIALOPROTEIN [BSP] (aa 108-122)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 4219.1 / A 4219.2
<b>Immunogen:</b>	Synthetic human Bone Sialoprotein (aa 108-122) (YTGLAAIQLPKKAGD)
<b>Synonyms:</b>	Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein
<b>Swiss-Prot No:</b>	P21815
<b>Gene Information:</b>	Gene Name: IBSP Gene ID: 3381
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Bone Sialoprotein (aa 108-122)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:1000), Western Blot (1:250), immunocytochemistry (1:100) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 April 2022

### For research use only

**Publishing research using A 4219? Please let us know so that we can cite your publication as a reference.**

