



Data Sheet

HUMAN BONE SIALOPROTEIN [BSP] (aa 39-47)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|---|
| Catalog no.: | A 4220.1 / A 4220.2 |
| Immunogen: | Synthetic human Bone Sialoprotein (aa 39-47) (LKRFPVQGG) |
| Synonyms: | Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein |
| Swiss-Prot No: | P21815 |
| Gene Information: | Gene Name: IBSP Gene ID: 3381 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human Bone Sialoprotein (aa 39-47) |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | ELISA (1:1000), immunocytochemistry (1:100) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 06 April 2022 |

For research use only

Publishing research using A 4220? Please let us know so that we can cite your publication as a reference.

