

Data Sheet

HUMAN BONE SIALOPROTEIN (aa 108-122)

SYNTHETIC

Catalog no.:	A 4219AG.1 / A 4219AG.2
Peptide:	Human Bone Sialoprotein (aa 108-122) (YTGLAAIQLPKKAGD)
Synonyms:	Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein
Swiss-Prot No:	P21815
Gene Information:	Gene Name: IBSP GeneID: 3381
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	04 May 2022

For research use only

Publishing research using A 4219AG? Please let us know so that we can cite your publication as a reference.

