

Data Sheet

HUMAN BONE SIALOPROTEIN

RECOMBINANT

Catalog no.:	A 2001AG.1 / A 2001AG.2
Protein:	Human Bone Sialoprotein, C-terminal His-tag
Synonyms:	Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein
Swiss-Prot No:	P21815
Gene Information:	Gene name: IBSP GenelD: 3381
Source:	Produced serum-free in CHO cell line
Purification:	Affinity chromatography
Storage buffer:	PBS, pH 7.4
Purity:	> 80% (SDS-PAGE), reduced band at ca. 80 kDa
Identity:	Immunoassay with specific monoclonal Ab
Contents:	10 µg / 100 µg (2.6 mg/ml)
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	05 October 2023

For research use only

Publishing research using A 2001AG? Please let us know so that we can cite your publication as a reference.

