



## Data Sheet

# HUMAN PRO-BRAIN NATRIURETIC PEPTIDE [proBNP] (aa 32-57)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE1027.1 / AE1027.2
<b>Immunogen:</b>	Synthetic human proBNP (aa 32-57) (QVEQTSLEPLQESPRPTGVVKSREVA)
<b>Swiss-Prot No:</b>	P16860
<b>Gene Information:</b>	Gene Name: NPPB GenelD: 4879
<b>Host:</b>	Goat
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human proBNP (aa 32-57) There was no cross reactivity obtained with human pro-BNP (aa 8-29).
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:5000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	20 April 2022

### For research use only

**Publishing research using AE1027? Please let us know so that we can cite your publication as a reference.**

